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REF 10012 - 4x80+4x20 ml

10010 - 4x250 ml

H10012 - 4x80+2x40 ml

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Kit for measurement of glucose in serum, plasma and urine

Colorimetric enzymatic method GOD-POD

PRINCIPLE

Glucose is oxidized, in presence of glucose oxidase (GOD), into gluconic acid and hydrogen peroxide. This one reacts, by peroxidase (POD), with 4-aminophenazone and phenol giving a coloured compound whose colour intensity is directly proportional to the glucose concentration in the tested sample.

REAGENTS

- R1** Phosphate buffer pH 7.4 100.0 mmol/l; phenol 9.0 mmol/l;
 GOD ≥ 15000 U/l; POD ≥ 1200 U/l
R2 4-aminophenazone 80.0 mmol/l
R3 Glucose standard 100 mg/dl (5.55 mmol/l)

SAMPLE

- Serum-heparinized plasma or EDTA plasma. Diluted urine 1:10.

Note

- Do not use samples with haemolysis.
- Centrifuge the samples to separate the serum from the cells. If it is not possible, draw the samples into test tubes containing fluoride or iodoacetate.
- The glucose is stable in the samples up to 7 days at 2-8°C, after the addition of a glycolytic inhibitor as NaF, KF.

REFERENCE VALUES

Serum - plasma	70 - 105 mg/dl (3.9 - 5.8 mmol/l)
Urine	< 0.5 g/24h (<28 mmol/24h)

References values are considered indicatives since each laboratory should establish references ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

STORAGE AND STABILITY

- Store the kit at 2-8°C. Do not freeze the reagents.
- After opening, the vials **R1**, **R2** and **R3 (STD)** are stable up to the expiry date if recapped immediately and protected from contamination, evaporation, direct light, and stored at the correct temperature.
- Working solution stability (**R1+ R2**): 60 days at 2-8°C.

PREPARATION OF REAGENTS

Reagents are liquid and ready to use. About preparing working solution add to every 4 ml of **R1** reagent, 1 ml of **R2** reagent. Prepare the amount for use according to the stability of working solution. Keep out the reagents from refrigerator only for the use and recap them immediately.

NOTE

- The kit, according to this method, must be used in manual procedures. About automatic using follow specific applications.
- Evaluate carefully the results if working reagent absorbance is > 0.300 at 510 nm.
- Avoid direct light, contamination and evaporation.
- The volumes in the procedure can be changed proportionally.
- In case of complaint or quality control request, refer to the lot number on the package or the lot number on the singles vials.

AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

Suggested serum:

REF 20350 Precise Norm **REF** 20360 Precise Path

PRECAUTION IN USE

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of inactive components (preservatives, surfactants and so on) is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances.

However the reagent should be handled with caution, according to good laboratory practice.

Caution: the reagents contain Sodium Azide (0.095%) as preservative. Avoid swallowing and contacting with skin, eyes and mucous membranes.

WASTE MANAGEMENT

Please refer to the local legal requirements.

PROCEDURE

Wavelength	λ: 510 (500-550) nm
Working temperature	37°C
Optical path	1 cm
Reaction	"end point"

Bring the reagents at 15-25°C before using them.

- MONOREAGENT PROCEDURE "sample starter"

	BLANK	STD	SAMPLE
WORKING REAGENT	1000 µl	1000 µl	1000 µl
DISTILLED WATER	10 µl	--	--
SAMPLE	--	--	10 µl
STANDARD	--	10 µl	--

Mix, then incubate for 10' at 37°C. Measure the absorbance of sample (EC) and standard (ESTD) against the reagent blank. The colour obtained is stable for about 60' at 15-25°C and protected from light.

- BIREAGENT PROCEDURE "substrate starter"

	BLANK	STD	SAMPLE
REAGENT R1	800 µl	800 µl	800 µl
DISTILLED WATER	10 µl	--	--
SAMPLE	--	--	10 µl
STANDARD	--	10 µl	--

Mix, incubate at 37°C for 1' and then add:
REAGENT R2 200 µl 200 µl 200 µl
 Mix, then incubate for 10' at 37°C. Measure the absorbance of sample (EC) and standard (ESTD) against the reagent blank. The colour obtained is stable for about 60' if stand at 15-25°C and protected from light.

CALCULATION

<u>Serum</u>	Glucose [mg/dl] o [mmol/l]	=	EC/ESTD x Conc. STD
<u>Urine</u>	Glucose [mg/dl] o [mmol/l]	=	EC/ESTD x Conc. STD x 10

CONVERSION FACTOR

Glucose [mg/dl] x 0.05551 = Glucose [mmol/l]

ANALYTICAL PERFORMANCES

Interferences

Bilirubin does not interfere up to concentration of 10 mg/dl.
 Triglycerides do not interfere up to concentration of 1000 mg/dl.
 Hemoglobin does not interfere up to concentration of 500 mg/dl.
 Ascorbate acid does not interfere up to concentration of 40 mg/dl

Linearity

Reaction is linear between 2-400 mg/dl (0.11-22.2 mmol/l). Samples with values exceeding 400 mg/dl must be diluted with saline solution. Multiply, then, the result for diluting factor.

"Intra-Assay" precision (within-Run)

Determined on 30 samples for each control (L-N-H) (Low-Normal-High).

Results:

MEAN	[mg/dl]	L = 49.40	N = 107.60	H = 267.50
S.D.		0.97	1.54	6.10
C.V.%		1.97	1.44	2.28

"Inter-Assay" precision (between-Run)

Determined on 15 samples for each control (L-H-N) for 3 days.

Results:

MEAN	[mg/dl]	L = 53.44	N = 106.33	H = 266.64
S.D.		1.16	1.48	5.93
C.V.%		2.19	1.40	2.20

Analytical sensitivity

The test sensitivity in terms of detection limit is 2 mg/dl (0.11 mmol/l).

Correlation

A study based comparing this method with a similar method on 21 samples has given a correlating factor **r = 0.98**

The urine performances are available on request.

BIBLIOGRAPHY

Henry, R.J., Clin.Chem., Hoeber N.Y., 625. (1964).
 Trinder, P., Ann. Clin. Biochem., 6, 24, (1969).
 Sharp, P., Clin. Chem. Acta, 40, 115, (1972).
 Kaplan, L.A., Pesce, A.J.: "Clinical Chemistry", Mosby Ed. (1996).

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer