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BIBLIOGRAFIA

Trinder P: J. Clin. Path., 22, 158 (1969).

Allain CC, Poon LS, Chan CS, Richmond W, Fu PC: Clin. Chem., 20, 470 (1974).

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Fabbricante

Kit for measurement of cholesterol in serum or plasma

Colorimetric enzymatic Trinder method CHOD-PAP

PRINCIPLE

Free cholesterol and cholesterol released from its esters after enzymatic hydrolysis are oxidized enzymatically (CHE). The free cholesterol is oxidized by cholesterol oxidase (CHOD) to colest-4-en-3-one giving hydrogen peroxide. The indicator quinoneimine is formed from hydrogen peroxide that in presence of peroxidase (POD), reacts with 4-aminophenazone and phenol to form a coloured compound. The colour intensity is directly proportional to the cholesterol total concentration in the tested sample.

REAGENTS

R1 Goods buffer pH 6.8 > 75.0 mmol/l; cholic acid > 0.4 mmol/l; phenol > 8 mmol/l; cholesterol esterase ≥ 500 U/l; cholesterol oxidase ≥ 300 U/l; peroxidase ≥ 300 U/l; conservatives and preservatives

R2 4-aminophenazone > 4 mmol/l

STD cholesterol 200 mg/dl (5.17 mmol/l)

SAMPLE

- Serum-heparinized plasma or EDTA plasma.

Note

- After picked up, serum and plasma must be separated, as soon as possible, from red cells. Avoid anticoagulants as fluoride, citrate and oxalate.
- The cholesterol is stable in the samples for 7 days if stand at 15-25°C or 28 days at 2-8°C.

REFERENCE VALUES

Serum - plasma	Values recommended:	< 200 mg/dl (5.17 mmol/l)
	Values borderline higher:	200-239 mg/dl (5.17 - 6.18 mmol/l)
	Values high:	≥ 240 mg/dl (6.21 mmol/l)

References values are considered indicatives since each laboratory should establish references ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

STORAGE AND STABILITY

- Store the kit at 2-8°C. Do not freeze the reagents.
- After opening, the reagents **R1**, **R2** and the **standard (R3)** are stable up to the expiry date if recapped immediately and protected from contamination, evaporation, direct light, and stored at the correct temperature.
- Working solution stability (**R1**+ **R2**): 60 days at 2-8°C.

PREPARATION OF REAGENTS

Reagents are liquid and ready to use. About using as monoreagent ("sample-starter" procedure) add to every 4 ml of **R1** reagent, 1 ml of **R2** reagent. Keep out the reagents from refrigerator only for the use and recap them immediately.

NOTE

- The kit, according to this method, must be used in manual procedures. About automatic using follow specific applications.
- Evaluate carefully the results if working reagent absorbance is > 0.300 at 510 nm.
- Avoid direct light, contamination and evaporation.
- The volumes in the procedure can be changed proportionally.
- In case of complaint or quality control request, refer to the lot number on the package or the lot number on the singles vials.
- A light pink color may appear in R1 reagent without affecting the use of the kit.

AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

Suggested serum:

REF 20350 Precise Norm **REF** 20360 Precise Path

PRECAUTION IN USE

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of components is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances.
However the reagent should be handled with caution, according to good laboratory practice.

Caution: the reagents contain Sodium Azide (0.095%) as preservative. Avoid swallowing and contacting with skin, eyes and mucous membranes.

WASTE MANAGEMENT

Please refer to the local legal requirements.

PROCEDURE

Wavelength λ: 510 (500-546) nm
Working temperature 37°C
Optical path 1 cm
Reaction "end point" (increasing)

Bring the reagents at 15-25°C before using them.

- MONOREAGENT PROCEDURE "sample starter"

WORKING REAGENT	BLANK	STD	SAMPLE
DISTILLED WATER	1000 µl	1000 µl	1000 µl
SAMPLE	10 µl	--	--
STANDARD	--	--	10 µl
	--	10 µl	--

Mix, then incubate for 5' at 37°C. Measure the absorbance of sample (EC) and standard (ESTD) against the reagent blank. The colour obtained is stable for about 60' if stand at 15-25°C and protected from light.

- BIREAGENT PROCEDURE "substrate starter"

REAGENT R1	BLANK	STD	SAMPLE
DISTILLED WATER	800 µl	800 µl	800 µl
SAMPLE	10 µl	--	--
STANDARD	--	--	10 µl
	--	10 µl	--

Mix, incubate at 37°C for 1' and then add:

REAGENT R2	BLANK	STD	SAMPLE
	200 µl	200 µl	200 µl

Mix, then incubate for 5' at 37°C. Measure the absorbance of sample (EC) and standard (ESTD) against the reagent blank. The colour obtained is stable for about 60' if stand at 15-25°C and protected from light.

CALCULATION

$$\text{Cholesterol [mg/dl] o [mmol/l]} = \text{EC/ESTD} \times \text{Conc. STD}$$

CONVERSION FACTOR

$$\text{Cholesterol [mg/dl]} \times 0.02586 = \text{Cholesterol [mmol/l]}$$

ANALYTICAL PERFORMANCES

Interferences

Bilirubin does not interfere up to concentration of 15 mg/dl.
Triglycerides do not interfere up to concentration of 1000 mg/dl.
Hemoglobin does not interfere up to concentration of 500 mg/dl.

Linearity

Reaction is linear between 2-800 mg/dl (0.052-20.69 mmol/l). Samples with values exceeding 800 mg/dl must be diluted with saline solution. Multiply, then, the result for diluting factor.

"Intra-Assay" precision (within-Run)

Determined on 30 samples for each control (L-N-H) (Low-Normal-High).

Results:

MEAN	(mg/dl)	L = 87.10	N = 182.80	H = 220.90
S.D.		1.81	3.50	3.46
C.V.%		2.08	1.92	1.57

"Inter-Assay" precision (between-Run)

Determined on 15 samples for each control (L-N-H) for 3 days.

Results:

MEAN	(mg/dl)	L = 90.07	N = 178.53	H = 225.40
S.D.		2.08	2.67	3.66
C.V.%		2.31	1.49	1.62

Analytical sensitivity

The test sensitivity in terms of detection limit is 2 mg/dl (0.052 mmol/l).

Correlation

A study based comparing this method with a similar method on 21 samples has given a correlating factor **r = 0.99**

BIBLIOGRAPHY

Trinder P: J. Clin. Path., 22, 158 (1969).
Allain CC, Poon LS, Chan CS, Richmond W, Fu PC: Clin. Chem., 20, 470 (1974).
Talke H, Schubert GE: Klin. Wchens., 43, 174 (1965).
Kaplan LA, Pesce AJ: "Clinical Chemistry", Mosby Ed. (1996).

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 CE regulation)

Storing temperature limits

In vitro medical device

Producer