Per uso diagnostico in vitro
10031 - 2x250ml
V10031 - 6x25 ml
H10031 - 4x100 ml

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Kit for measurement of total protein in serum or plasma Biuret method

PRINCIPLE

In alkaline environment proteins react with bivalent ions of copper giving a coloured complex whose colour intensity is directly proportional to the total proteins concentration in the tested sample.

REAGENTS

R1 Sodium hydroxide 350.0 mmol/l; sodium potassium tartrate 20.0 mmol/l;
 iodide potassium 5.20 mmol/l; sulphate copper 4.8 mmol/l
R2 Protein standard 4 g/dl

SAMPLE

- Serum or heparinized plasma.

Note

- Do not use sample with haemolysis.
- The proteins are stable in the samples up to 7 days at 15-25°C and 30 days at 2-8°C.

REFERENCE VALUES

Serum - plasma	6.6 - 8.3 g/dl
----------------	----------------

References values are considered indicatives since each laboratory should establish references ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

STORAGE AND STABILITY

- Store the kit at 2-8°C (the reagent R1 can be stored at 15-25°C).
- After opening, the vials **R1** and **R2** are stable up to the expiry date if recapped immediately and protected from contamination, evaporation, direct light, and stored at the correct temperature.

PREPARATION OF REAGENTS

Reagents are liquid and ready to use.

NOTE

- The kit, according to this method, must be used in manual procedures. About automatic using follow specific applications.
- Evaluate carefully the results if reagent R1 absorbance is > 0.400 at 540 nm.
- Avoid direct light, contamination and evaporation.
- The volumes in the procedure can be changed proportionally.
- In case of complaint or quality control request, refer to the lot number on the package or the lot number on the singles vials.

AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

Suggested serum:

REF 20350 Precise Norm **REF** 20360 Precise Path

PRECAUTION IN USE

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of inactive components (preservatives, surfactants and so on) is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances.

However the reagent should be handled with caution, according to good laboratory practice.

WASTE MANAGEMENT

Please refer to the local legal requirements.

PROCEDURE

Wavelength λ : 540 (530-570) nm
 Working temperature 37°C (15-25°C)
 Optical path 1 cm
 Reaction "end point"

Bring the reagents at 15-25°C before using them.

- PROCEDURE

	BLANK	STD	SAMPLE
REAGENT R1	1000 μ l	1000 μ l	1000 μ l
DISTILLED WATER	10 μ l		
SAMPLE	--	--	10 μ l
STANDARD	--	10 μ l	--

Mix, then incubate 10' at 37°C (15-25°C). Measure the absorbance of sample (EC) and standard (ES) against the reagent blank.

CALCULATION

Total proteins [g/dl] = EC / ES x Conc. STD

The reagent performances are related to 37°C and 540 nm.

ANALYTICAL PERFORMANCES

Interferences

Bilirubin does not interfere up to concentration of 20 mg/dl.
 Hemoglobin does not interfere up to concentration of 2.5 g/dl.

Linearity

Reaction is between 0.3-10 g/dl. Samples with values exceeding 10 g/dl must be diluted with saline solution. Multiply, then, the result for diluting factor.

"Intra Assay" precision (within-Run)

Determined on 30 samples for each control (L-N-H) (Low-Normal-High).

Results:

MEAN	[mg/dl]	L = 3.3	N = 5.9	H = 7.2
S.D.		0.15	0.11	0.22
C.V.%		4.54	1.84	3.04

"Inter Assay" precision (between-Run)

Determined on 15 samples for each control (L-N-H) for 3 days.

Results:

MEAN	[mg/dl]	L = 2.76	N = 5.13	H = 7.18
S.D.		0.07	0.36	0.10
C.V.%		2.65	7.20	1.39

Analytical sensitivity

The test sensitivity in terms of detection limit is 0.3 g/dl.

Correlation

A study based comparing this method with a similar method on 21 samples has given a correlating factor **r = 0.97**

BIBLIOGRAPHY

Henry, R.J.: Clinical Chemistry, Hoeber, N.Y. 413, (1976)
 Tietz, N.W.: Fundamentals of Clinical Chemistry, Saunders Co., Philadelphia, PA 302 (1970)
 Kaplan LA, Pesce AJ: "Clinical Chemistry", Mosby Ed. (1996).

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer