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 Fabbricante

Kit for measurement of calcium in serum, plasma and urine Colorimetric o-cresolphthalein complexone method CPC

PRINCIPLE

O-cresolphthalein complexone (CPC) reacts with calcium ions and magnesium at alkaline pH giving a strongly coloured complex: the magnesium interference is inhibited by 8-hydroxyquinoline presence. The colour intensity is directly proportional to the calcium concentration in the tested sample.

REAGENTS

R1 AMP buffer pH 10.3 330.0 mmol/l
R2 O-cresolphthalein complexone 0.16 mmol/l; 8-hydroxyquinoline 9.0 mmol/l
Calcium standard 10 mg/dl (2.50 mmol/l)

SAMPLE

- Serum or plasma. Diluted urine 1:5.

Note

- Do not use samples with haemolysis. Do not use anticoagulants as EDTA, oxalate, citrate or fluoride.
- Avoid venous standstill. The use of tourniquet can rise the calcium level in the withdrawal in 0.5 mg/dl (0.12 mmol/l) as well.
- After picked up, serum and plasma must be separated, as soon as possible, from red cells to avoid calcium absorption by red cells.
- The calcium is stable in the serum or plasma 1 day at 2-8°C or 8 months at -20°C.
- Serum strongly jaundiced or lipoemic must be tested using a blank sample.

REFERENCE VALUES

Serum Adults	8.5 - 10.8 mg/dl
Serum Children	8.5 - 12 mg/dl
Urine	100 - 300 mg/24h

References values are considered indicatives since each laboratory should establish references ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

STORAGE AND STABILITY

- Store the kit at 15-25°C. Do not freeze the reagents.
- After opening, the vials **R1**, **R2** and **STD** are stable up to the expiry date if recapped immediately and protected from contamination, evaporation, direct light, and stored at the correct temperature.
- Working solution stability (**R1**+ **R2**): 6 hours at 15-25°C.

PREPARATION OF REAGENTS

Reagents are liquid and ready to use. About using as monoreagent ("sample-starter" procedure) mix the reagents **R1** and **R2** in equal parts.

NOTE

- The kit, according to this method, must be used in manual procedures. About automatic using follow specific applications.
- Evaluate carefully the results if working reagent absorbance is > 0.500 at 570 nm.
- Avoid direct light, contamination and evaporation.
- The volumes in the procedure can be changed proportionally.
- In case of complaint or quality control request, refer to the lot number on the package or the lot number on the singles vials.

AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

Suggested serum:

REF 20350 Precise Norm **REF** 20360 Precise Path

PRECAUTION IN USE

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of components is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances.

However the reagent should be handled with caution, according to good laboratory practice.

Caution: the reagents contain Sodium Azide (0.095%) as preservative. Avoid swallowing and contacting with skin, eyes and mucous membranes.

WASTE MANAGEMENT

Please refer to the local legal requirements.

PROCEDURE

Wavelength λ : 570 (550-590) nm
Working temperature 37°C (15-25°C)
Optical path 1 cm
Reaction "end point" (increasing)

- MONOREAGENT PROCEDURE "sample starter"

	BLANK	STD	SAMPLE
WORKING REAGENT	1000 μ l	1000 μ l	1000 μ l
DISTILLED WATER	25 μ l	--	--
SAMPLE	--	--	25 μ l
STANDARD	--	25 μ l	--

Mix, then incubate for 5' at 37°C (15-25°C). Measure the absorbance of the sample (EC) and standard (ESTD) against the reagent blank.

- BIREAGENT PROCEDURE "R2 starter"

	BLANK	STD	SAMPLE
REAGENT R1	500 μ l	500 μ l	500 μ l
DISTILLED WATER	25 μ l	--	--
SAMPLE	--	--	25 μ l
STANDARD	--	25 μ l	--

Mix, incubate at 15-25°C for 1' and then add:
REAGENT R2 500 μ l 500 μ l 500 μ l
Mix, then incubate for 5' at 37°C (15-25°C). Measure the absorbance of the sample (EC) and standard (ESTD) against the reagent blank.

CALCULATION

$$\text{Calcium [mg/dl] o [mmol/l]} = \text{EC/ESTD} \times \text{Conc. STD}$$

Diluted urines: multiply the result for diluting factor.

CONVERSION FACTOR

Calcium [mg/dl] \times 0.2495 = Calcium [mmol/l]
Calcium [mg/dl] \times 0.4990 = Calcium [mEq/l]

ANALYTICAL PERFORMANCES

The reagent performances are related to 37°C and 570 nm.

Interferences

Bilirubin does not interfere up to concentration of 40 mg/dl.
Triglycerides do not interfere up to concentration of 2000 mg/dl.
Hemoglobin does not interfere up to concentration of 500 mg/dl.
Magnesium does not interfere up to concentration of 15 mg/dl.
Ascorbate acid does not interfere up to concentration of 30 mg/dl.

Linearity

Reaction is linear between 1.5-15 mg/dl (0.37-3.74 mmol/l). Samples with values exceeding 15 mg/dl must be diluted with saline solution. Multiply, then, the result for diluting factor.

"Intra-Assay" precision (within-Run)

Determined on 30 samples for each control (L-N-H) (Low-Normal-High).

Results:

MEAN (mg/dl)	L = 4.10	N = 8.50	H = 12.10
S.D.	0.16	0.25	0.37
C.V.%	4.00	2.94	3.04

"Inter-Assay" precision (between-Run)

Determined on 15 samples for each control (L-H-N) for 3 days.

Results:

MEAN (mg/dl)	L = 3.57	N = 7.78	H = 12.22
S.D.	0.17	0.20	0.43
C.V.%	4.68	2.63	3.46

Analytical sensitivity

The test sensitivity in terms of detection limit is 1.5 mg/dl (0.37 mmol/l).

Correlation

A study based comparing this method with a similar method on 21 samples has given a correlating factor **r = 0.89**

The urine performances are available on request.

BIBLIOGRAPHY

Ray Sarkar, B.C. et al., Anal. Bioch., 20,155, (1967).
Barnett, R.N. et al., Amer. J. Clin. Oath., 59, 836, (1973).
Kaplan, L.A., Pesce, A.: "Clinical Chemistry", Mosby Ed. (1996).

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer