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Fabbricante

Kit for measurement of cholinesterase in serum or plasma Kinetic butyrylthiocholine method

PRINCIPLE

Cholinesterase catalyzes the butyrylthiocholine hydrolysis into thiocoline, which, in the presence of 5-5'-dithiobis-2-nitrobenzoic acid, forms a yellow compound which absorbs at 405 nm. The increase of absorbance is proportional to the cholinesterase activity in the sample.

REAGENTS

- R1** 5-5'-dithiobis-2-nitrobenzoic acid 0.26 mmol/l;
phosphate buffer 42.0 mmol/l
R2 Butyrylthiocholine iodide 420.0 mmol/l

In case of unitary test the components are joined in a single vial.

SAMPLE

- Serum -heparinized plasma or EDTA plasma.

Note

- Do not use samples with haemolysis.
- Separate, as soon as possible, serum and plasma from red cells after picking because the cholinesterase activity can increase till 25% if the serum is left one day into contact with the corpuscolate part.
- Do not use sodium fluoride as anticoagulant because it inhibits the enzyme activity.
- The cholinesterase is stable in the samples up to 15 days at 2-8°C.

REFERENCE VALUES

Serum - plasma (U/l) 37°C	4850 -12000
------------------------------	-------------

Children till 6 months have a cholinesterase activity higher than the adults from 40% to 50%.

Reference values are considered indicative since each laboratory should establish reference ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

PREPARATION OF REAGENTS

Dissolve the content of **R1** and **R2** vials with deionized water according to the quantity indicated on the label. Mix a portion of **R2** reagent with 20 parts of **R1** solution.

In case of **unitary test (Ref. 10067)** dissolve the content of **R1** reagent with 3 ml of deionized water.

Keep out the reagents from refrigerator only for the use and recap them immediately.

STORAGE AND STABILITY

- Store the kit at 2-8°C. Do not freeze the reagents.
- Once reconstituted, the **R1-R2** reagents are stable 28 days at 2-8°C.
- Working solution is stable 3 hours at 15-25°C.
- In case of unitary test the **R1** reagent, once reconstituted, is stable 3 hours at 15-25°C.

NOTE

- The kit, according to this method, must be used in manual procedures. About automatic using follow specific applications.
- Evaluate carefully the results if working reagent absorbance is > 0.400 at 405 nm.
- Avoid direct light, contamination and evaporation.
- The volumes in the procedure can be changed proportionally.
- In case of complaint or quality control request, refer to the lot number on the package or the lot number on the single vials.

AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

Suggested serum:

REF 20350 Precise Norm **REF** 20360 Precise Path

PRECAUTION IN USE

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of components is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances. However the reagent should be handled with caution, according to good laboratory practice.

WASTE MANAGEMENT

Please refer to the local legal requirements.

PROCEDURE

Wavelength	λ: 405 (400-440) nm
Working temperature	37°C
Optical path	1 cm
Reaction	"kinetic" (increasing)

Bring the reagents at 15-25°C before using them.

	BLANK	SAMPLE
WORKING REAGENT	1000 µl	1000 µl
DISTILLED WATER	5 µl	--
SAMPLE (*)	--	5 µl

Mix, then incubate for 1'. Measure the absorbance of sample (EC) and blank (EB) in increasing after 30", 60" and 90". Calculate the average absorbance difference per minute ($\Delta E/\text{min}$).

In case of unitary test: (Ref. 10067):

SAMPLE REAGENT	15 µl 3000 µl
-----------------------	------------------

CALCULATION

$$\text{Cholinesterase [U/l]} = \Delta E/60\text{sec} \times 15113$$

The factor and the reagent performances are related to 405 nm and 37°C.

ANALYTICAL PERFORMANCES

Interferences

Bilirubin does not interfere up to concentration of 200 mg/dl.
Triglycerides do not interfere up to concentration of 1000 mg/dl.

Linearity

Reaction is linear up to a concentration of 12000 U/l. Samples with values exceeding 12000 U/l must be diluted with saline solution. Multiply, then, the result for diluting factor.

"Intra-Assay" precision (within-Run)

Determined on 30 samples for each control (L-N-H) (Low-Normal-High).

Results:

MEAN [U/l]	L = 3721.47	N = 5464.37	H = 7548.67
S.D.	42.48	71.08	132.65
C.V.%	1.14	1.30	1.76

"Inter-Assay" precision (between-run)

Determined on 15 samples for each control (L-N-H) for 3 days.

Results:

MEAN [U/l]	L = 3728.84	N = 5549.87	H = 7649.94
S.D.	63.92	66.55	71.02
C.V.%	1.71	1.21	0.93

Correlation

A study based comparing this method with a similar method on 21 samples has given a correlating factor **r = 1.00**

BIBLIOGRAPHY

- Knedel B., Boettger R., Klin. Wschr., 45, 325, (1967).
 Jakobs D.S., Kasten Jr. B.L., Demmott W.R., Wolfson W.L.: "Laboratory Test Handbook", Lexi-Comp and Williams & Wilkins Ed. (2nd Edition - 1990).
 Deutsche Gesellschaft für Klinische Chemie: "Proposal of Standard Methods for the determination of enzyme catalytic concentrations in serum and plasma at 37°C, II. Cholinesterase (Acylcholine Acylhydrolase)". Eur. J. Clin. Chem. Clin. Biochem. 30, 163 (1992).
 Kaplan, L.A., Pesce, A..J.: "Clinical Chemistry", Mosby Ed. (1996).

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer