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Fabbricante

## Kit for measurement of cholinesterase in serum or plasma

### Colorimetric method optimized DGKC\*

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#### PRINCIPLE

Cholinesterase catalyzes the butyrylthiocholine hydrolysis at thiocholine in presence of potassium hexacyanoferrate (III) which is reduced into potassium hexacyanoferrate (II). The reduction absorbance velocity is proportional to the cholinesterase activity in the sample.

#### REAGENTS

**R1** Pirophosphate pH 7.7 65.0 mmol/l; hexacyanoferrate (III) 2.0 mmol/l  
**R2** Goods buffer pH 4.0 20.0 mmol/l;  
butyrylthiocholine iodide 65.0 mmol/l

#### SAMPLE

- Serum -heparinized plasma or EDTA plasma.

#### Note

- Do not use samples with haemolysis.
- Separate, as soon as possible, serum and plasma from red cells after picking because the cholinesterase activity can increase till 25% if the serum is left one day into contact with the corpuscolate part.
- Do not use sodium fluoride as anticoagulant because inhibit the enzyme activity.
- The cholinesterase is stable in the samples up to 15 days at 2-8°C.

#### REFERENCE VALUES

Serum - plasma	Adults (to 37°C)
Males:	5100 - 11700 U/l
Females:	4000 - 12600 U/l

Children till 6 months have a cholinesterase activity higher than the adults from 40% to 50%.

References values are considered indicatives since each laboratory should establish references ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

#### STORAGE AND STABILITY

- Store the kit at 2-8°C. Do not freeze the reagents.
- After opening, the vials **R1** and **R2** are stable up to the expiry date if recapped immediately and protected from contamination, evaporation, direct light, and stored at the correct temperature.
- Working solution stability (**R1**+ **R2**): 30 days at 2-8°C.

#### PREPARATION OF REAGENTS

Reagents are liquid and ready to use. About using as monoreagent (sample-starter procedure) add to every 4 ml of **R1** reagent, 1 ml of **R2** reagent. Keep out the reagents from refrigerator only for the use and recap them immediately.

#### NOTE

- The kit, according to this method, must be used in manual procedures. About automatic using follow specific applications.
- Evaluate carefully the results if working reagent absorbance is < 1.000 at 405 nm.
- Avoid direct light, contamination and evaporation.
- The volumes in the procedure can be changed proportionally.
- In case of complaint or quality control request, refer to the lot number on the package or the lot number on the singles vials.

#### AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

#### QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

#### Suggested serum:

REF 20350 Precise Norm REF 20360 Precise Path

#### PRECAUTION IN USE

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of components is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances. However the reagent should be handled with caution, according to good laboratory practice.

#### WASTE MANAGEMENT

Please refer to the local legal requirements.

#### PROCEDURE

Wavelength λ: 405 (400-440) nm  
Working temperature 37°C  
Optical path 1 cm  
Reaction "kinetic" (decreasing)

Bring the reagents at 15-25°C before using them.

#### - MONOREAGENT PROCEDURE "sample starter"

	BLANK	SAMPLE
WORKING REAGENT	1000 µl	1000 µl
DISTILLED WATER	15 µl	--
SAMPLE	--	15 µl

Mix, then incubate for 1'. Measure the absorbance of sample (EC) and blank (EB) in decreasing after 30", 60" and 90". Calculate the average absorbance difference per minute (ΔE/min).

#### - BIREAGENT PROCEDURE "substrate starter"

	BLANK	SAMPLE
REAGENT R1	800 µl	800 µl
DISTILLED WATER	15 µl	--
SAMPLE	--	15 µl

Mix, incubate for 1' and add:

	BLANK	SAMPLE
REAGENT R2	200 µl	200 µl

Mix, then incubate for 1'. Measure the absorbance of sample (EC) and blank (EB) in decreasing after 30", 60" and 90". Calculate the average absorbance difference per minute (ΔE/min).

#### CALCULATION

$$\text{Cholinesterase [U/l]} = (\Delta E/\text{min C} - \Delta E/\text{min B}) \times 72000$$

The factor and the reagent performances are related to 405 nm and 37°C.

#### ANALYTICAL PERFORMANCES

##### Interferences

Bilirubin does not interfere up to concentration of 200 mg/dl.  
Triglycerides do not interfere up to concentration of 1000 mg/dl.

##### Linearity

Reaction is linear up to a value of 12000 U/l (corresponding to a ΔE/min of 0.160 at 405 nm). Samples with values exceeding 12000 U/l must be diluted with saline solution. Multiply, then, the result for diluting factor.

##### "Intra-Assay" precision (within-Run)

Determined on 30 samples for each control (L-N-H) (Low-Normal-High).

Results:

MEAN	[U/l]	L = 3709.00	N = 5387.43	H = 7487.07
S.D.		65.67	81.87	65.49
C.V.%		1.77	1.52	0.87

##### "Inter-Assay" precision (between-run)

Determined on 15 samples for each control (L-N-H) for 3 days.

Results:

MEAN	[U/l]	L = 3691.91	N = 5399.84	H = 7455.38
S.D.		76.76	59.41	68.41
C.V.%		2.08	1.10	0.92

##### Analytical sensitivity

The test sensitivity in terms of detection limit is 97 U/l.

##### Correlation

A study based comparing this method with a similar method on 21 samples has given a correlating factor **r = 1.00**

##### BIBLIOGRAPHY

Knedel B., Boettger R., Klin. Wschr., 45, 325, (1967).  
Jakobs D.S., Kasten Jr. B.L., Demmott W.R., Wolfson W.L.: "Laboratory Test Handbook", Lexi-Comp and Williams & Wilkins Ed. (2nd Edition - 1990).  
Kaplan, L.A., Pesce, A..J.: "Clinical Chemistry", Mosby Ed. (1996).  
Deutsche Gesellschaft für Klinische Chemie: "Proposal of Standard Methods for the determination of enzyme catalytic concentrations in serum and plasma at 37°C, II. Cholinesterase (Acylcholine Acylhydrolase)". Eur. J. Clin. Chem. Clin. Biochem. 30, 163 (1992).

##### SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Cholinesterase LR

30192 Rev. 2 del 2012/01