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Fabbricante

## Kit for measurement of creatinine in serum, plasma and urine Colorimetric method Jaffé without deproteinization

### PRINCIPLE

Creatinine reacts in alkaline environment with picrate to give a coloured compound whose intensity is proportional to the creatinine concentration in the sample.

### REAGENTS

**R1** Lithium hydroxide 120.0 mmol/l; boric acid 80.0 mmol/l  
**R2** Picric acid 67.0 mmol/l  
**R3 Creatinine standard** 2 mg/dl (176.8 µmol/l)

### SAMPLE

- Serum or plasma.  
- Urine. (dilute the sample 1:10 prior to testing and multiply the result by the diluting factor)

### Note

- Do not use samples with haemolysis.  
- The creatinine is stable in the samples 7 days at 4°C.

### REFERENCE VALUES

Sample	Gender	Reference Range
Serum - plasma	Men	0.9 - 1.3 mg/dl (80 - 115 µmol/l)
	Women	0.6 - 1.1 mg/dl (53 - 97 µmol/l)
Urine		0.80 - 1.80 g/24h

References values are considered indicatives since each laboratory should establish references ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

### STORAGE AND STABILITY

- Store the kit at 15-25°C.  
- After opening, the vials **R1**, **R2** and the **Standard (R3)** are stable up to the expiry date if recapped immediately and protected from contamination, evaporation, direct light, and stored at the correct temperature.  
- Working solution stability (**R1+ R2**): 14 days at 2-8°C.

### PREPARATION OF REAGENTS

Reagents are liquid and ready to use. About using as monoreagent ("sample-starter" procedure) mix the reagents **R1** and **R2** in equal parts.

### NOTE

- The kit, according to this method, must be used in manual procedures. About automatic using follow specific applications.  
- Evaluate carefully the results if working reagent absorbance is > 0.400 at 510 nm.  
- Avoid direct light, contamination and evaporation.  
- The reaction is very sensitive to the temperature, so the temperatures, after set up, must be kept constant during the dosage duration.  
- The volumes in the procedure can be changed proportionally.  
- In case of complaint or quality control request, refer to the lot number on the package or the lot number on the singles vials.

### AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

### QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

#### Suggested serum:

REF 20350 Precise Norm                      REF 20360 Precise Path

### PRECAUTION IN USE

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of inactive components (preservatives, surfactants and so on) is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances.  
However the reagent should be handled with caution, according to good laboratory practice.

### WASTE MANAGEMENT

Please refer to the local legal requirements.

### PROCEDURE

Wavelength λ: 510 (500-550) nm  
Working temperature 37°C (15-25°C)  
Optical path 1 cm  
Reaction "fixed time"

#### - MONOREAGENT PROCEDURE "sample starter"

	BLANK	STD	SAMPLE
WORKING REAGENT	1000 µl	1000 µl	1000 µl
DISTILLED WATER	100 µl	--	--
SAMPLE	--	--	100 µl
STANDARD	--	100 µl	--

Mix and after 30" at 37°C (15-25°C) measure the absorbance of sample (E1C) and standard (E1STD). After other 2' make the second reading (E2C), (E2STD).

#### - BIREAGENT PROCEDURE "R2 starter"

	BLANK	STD	SAMPLE
REAGENT R1	500 µl	500 µl	500 µl
DISTILLED WATER	100 µl	--	--
SAMPLE	--	--	100 µl
STANDARD	--	100 µl	--

Mix, incubate at 37°C for 1', then add:

	BLANK	STD	SAMPLE
REAGENT R2	500 µl	500 µl	500 µl

Mix and after 30" at 37°C (15-25°C) measure the absorbance of sample (E1C) and standard (E1STD). After other 2' make the second reading (E2C), (E2STD).

### CALCULATION

$$\text{Creatinine [mg/dl] o } [\mu\text{mol/l}] = \frac{(E2C - E1C)}{(E2STD - E1STD)} \times \text{Conc. STD}$$

The reagent performances are related to 37°C and 510 nm.

### CONVERSION FACTOR

Creatinine [mg/dl] x 88.4 = Creatinine [µmol/l]

### ANALYTICAL PERFORMANCES

#### Interferences

Bilirubin does not interfere up to concentration of 10 mg/dl.  
Triglycerides do not interfere up to concentration of 2000 mg/dl.  
Hemoglobin does not interfere up to concentration of 100 mg/dl.

#### Linearity

Reaction is linear up to a concentration of 25 mg/dl (2210 µmol/l) with a range of 0.3-25 mg/dl (27-2210 µmol/l). Samples with values exceeding 25 mg/dl must be diluted with saline solution. Multiply, then, the result for diluting factor.

#### "Intra-Assay" precision (within-Run)

Determined on 30 samples for each control (L-N-H) (Low-Normal-High).

Results:

MEAN	[mg/dl]	L = 0.95	N = 1.75	H = 4.94
S.D.		0.06	0.09	0.20
C.V.%		6.32	4.97	4.00

#### "Inter Assay" precision (between-Run)

Determined on 15 samples for each control (L-N-H) for 3 days.

Results:

MEAN	[mg/dl]	L = 0.86	N = 1.32	H = 5.18
S.D.		0.05	0.08	0.21
C.V.%		6.12	6.51	4.00

#### Analytical sensitivity

The test sensitivity in terms of detection limit is 0.3 mg/dl (27 µmol/l).

#### Correlation

A study based comparing this method with a similar method on 21 samples has given a correlating factor **r = 0.97**

The urine performances are available on request.

#### BIBLIOGRAPHY

Henry, R.J.: Clinical Chemistry: Principles and Techniques, Harper & Row.N.Y., p. 287, (1964).

Young et al., Clinical Chemistry, 18, (1972).

Newman, D.J., Price, C.P.: Renal function and nitrogen metabolites. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3<sup>rd</sup> ed. Philadelphia: W.B. Saunders Company; p. 1204, (1999).

Kaplan, L.A., Pesce, A.J.: "Clinical Chemistry", Mosby Ed. (1996).

#### SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer