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Fabbricante

Kit for measurement of alkaline phosphatase in serum or plasma Colorimetric optimized method DEA DGKC*

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PRINCIPLE

Alkaline phosphatase (ALP) catalyzes the hydrolysis, in alkaline environment, of p-nitrophenylphosphate into p-nitrophenol and phosphate. The absorbance increasing is directly proportional to the alkaline phosphatase activity in the sample.

REAGENTS

R1 Diethanolamine buffer pH 9.8 2.2 mol/l; magnesium chloride 1.0 mmol/l
R2 P-nitrophenylphosphate 38.0 mmol/l

SAMPLE

- Serum or heparinized plasma.

Note

- Do not use samples with haemolysis. The withdrawal must be made, preferably, after 8 hours of fasting. The anticoagulants as EDTA, citrate and oxalate must not be used because they inhibit the enzyme activity.
 - The ALP is stable in the samples up to 7 days at 2-8°C.

REFERENCE VALUES

Serum - plasma [U/l]	25°C	30°C	37°C
Children	110-720	145-950	180-1200
Adults	60-170	80-220	100-290

References values are considered indicatives since each laboratory should establish references ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

STORAGE AND STABILITY

- Store the kit at 2-8°C. Do not freeze the reagents.
 - After opening, the vials **R1** and **R2** are stable up to the expiry date if recapped immediately and protected from contamination, evaporation, direct light, and stored at the correct temperature.
 - Working solution stability (**R1+ R2**): 30 days at 2-8°C.

PREPARATION OF REAGENTS

Reagents are liquid and ready to use. About using as monoreagent (sample-starter procedure) add to every 4 ml of **R1** reagent, 1 ml of **R2** reagent. Keep out the reagents from refrigerator only for the use and recap them immediately.

NOTE

- The kit, according to this method, must be used in manual procedures. About automatic using follow specific applications.
 - Evaluate carefully the results if working reagent absorbance is > 0.600 at 405 nm.
 - Avoid direct light, contamination and evaporation.
 - The volumes in the procedure can be changed proportionally.
 - In case of complaint or quality control request, refer to the lot number on the package or the lot number on the singles vials.

AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

Suggested calibrator:

REF 20102 Biocal H Human

QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

Suggested serum:

REF 20350 Precise Norm REF 20360 Precise Path

PRECAUTION IN USE

H318: Causes serious eye damage.

P280: Wear eye protection/ face protection.

P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310: Immediately call a POISON CENTER or doctor/physician.

The reagents contain Sodium Azide (0.095%) as preservative. Avoid swallowing and contacting with skin, eyes and mucous membranes.

WASTE MANAGEMENT

Please refer to the local legal requirements.

PROCEDURE

Wavelength λ : 405 (400-410) nm
 Working temperature 37°C (25° o 30°C)
 Optical path 1 cm
 Reaction kinetic (increasing)

Bring the reagents at 15-25°C before using them.

- MONOREAGENT PROCEDURE "sample starter"

WORKING REAGENT	BLANK	SAMPLE
DISTILLED WATER	1000 μ l	1000 μ l
SAMPLE	20 μ l	--
	--	20 μ l

Mix, then incubate for 1' a 37°C. Measure the absorbance of sample (EC) against distilled water decrease per minute during 3 minutes. Calculate the average absorbance difference per minute ($\Delta E/\text{min}$).

- BIREAGENT PROCEDURE "substrate starter"

REAGENT R1	BLANK	SAMPLE
DISTILLED WATER	800 μ l	800 μ l
SAMPLE	20 μ l	--
	--	20 μ l

Mix, incubate at 37°C for 1' and then add:
REAGENT R2 200 μ l 200 μ l
 Mix, then incubate for 1' a 37°C. Measure the absorbance of sample (EC) against distilled water decrease per minute during 3 minutes. Calculate the average absorbance difference per minute ($\Delta E/\text{min}$).

CALCULATION with calibration

$$\text{ALP [U/l]} = \Delta E/\text{min sample} / \Delta E/\text{min calib} \times \text{conc calib}$$

CALCULATION with factor (indicative)

$$\text{ALP [U/l]} = \Delta E/\text{min} \times 2757$$

The factor and the reagent performances are related to 37°C and 405 nm.

ANALYTICAL PERFORMANCES

Interferences

Bilirubin does not interfere up to concentration of 40 mg/dl.
 Triglycerides do not interfere up to concentration of 2000 mg/dl.
 Hemoglobin does not interfere up to concentration of 150 mg/dl.
 Ascorbate acid does not interfere up to concentration of 30 mg/dl.

Linearity

Reaction is linear up to a concentration of 2000 U/l (corresponding to $\Delta E/\text{min}$ of 0.725). Samples with values exceeding 2000 U/l must be diluted with saline solution. Multiply, then, the result for diluting factor.

"Intra-Assay" precision (within-Run)

Determined on 30 samples for each control (L-N-H) (Low-Normal-High).

Results:

MEAN (U/l)	L = 83.3	N = 166.7	H = 611.0
S.D.	3.76	4.15	5.62
C.V.%	4.52	2.49	0.92

"Inter-Assay" precision (between-run)

Determined on 15 samples for each control (L-N-H) for 3 days.

Results:

MEAN (U/l)	L = 86.58	N = 164.22	H = 568.98
S.D.	2.36	4.02	9.72
C.V.%	2.73	2.45	1.71

Analytical sensitivity

The test sensitivity in terms of detection limit is 19 U/l.

Correlation

A study based comparing this method with a similar method on 21 samples has given a correlating factor $r = 0.98$

BIBLIOGRAPHY

Deutsche Gesellschaft für klinische Chemie (DGKC). Recommendation of the German Society of Clinical Chemistry. Standardization of methods for measurement of enzymatic activities in biological fluids. Z. Klin. Chem. Klin. Biochem., 10, 182 (1972). Kubler W: Symp. D. Deutsch Ges. für Lab. Med., 8 (1973). Fischbach F, Zawta B: Klin. Lab., 38, 555 (1992). Szasz G. et al.: Z. Kndrheilik, 111, 233 (1971). Kaplan LA, Pesce AJ: "Clinical Chemistry", Mosby Ed. (1996).

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer