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Fabbricante

Kit for measurement of total bilirubin in serum or plasma
Colorimetric method Jendrassik-Gróf modified

Xn

R22 - Harmful if swallowed

PRINCIPLE

Total bilirubin reacts in acid environment with sulphaniilic acid and, in presence of specific reaction catalyzers, gives a coloured azocompound (azobilirubin) whose colour intensity is proportional to the total bilirubin concentration in the sample.

REAGENTS

R1 Sulphanilic acid 9.98 mmol/l; citric acid 0.95 mol/l; caffeine 0.49 mol/l;
urea 2.6 mol/l; preservatives and surface-active agents not anionic
R2 Nitrite sodium 12.4 mmol/l

SAMPLE

- Serum or EDTA-Na₂ plasma.

Note

- Do not use samples with haemolysis.
- Keep the samples far from light and heat because the bilirubin is a photosensitive pigment. Perform the samples as soon as possible.
- It's advisable to perform a blank sample for those samples which are hyperlipopemic and not clear.

REFERENCE VALUES

Adults	up to 1.2 mg/dl (17.1 µmol/l)
--------	-------------------------------

References values are considered indicatives since each laboratory should establish references ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

STORAGE AND STABILITY

- Store the kit at 15-25°C. Do not freeze the reagents.
- After opening, the vials **R1** and **R2** are stable up to the expiry date if recapped immediately and protected from contamination, evaporation, direct light, and stored at the correct temperature.
- Working solution stability (**R1** + **R2**): 21 days at 2-8°C.

PREPARATION OF REAGENTS

Reagents are liquid and ready to use. About using as monoreagent ("sample-starter" procedure) add to every 4 ml of **R1** reagent, 1 ml of **R2** reagent. Prepare the amount for use according to the stability of working solution.

NOTE

- The kit, according to this method, must be used in manual procedures. About automatic using follow specific applications.
- Evaluate carefully the results if working reagent absorbance is > 0.100 at 570 nm.
- Avoid direct light, contamination and evaporation.
- The volumes in the procedure can be changed proportionally.
- In case of complaint or quality control request, refer to the lot number on the package or the lot number on the singles vials.

AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

Suggested serum:

REF 20350 Precise Norm REF 20360 Precise Path

PRECAUTION IN USE

Look at the caution advices indicated nearby the symbol.

WASTE MANAGEMENT

Please refer to the local legal requirements.

PROCEDURE

Wavelength λ: 570 (550-580) nm
Working temperature 37°C
Optical path 1 cm
Reaction "end point"

- MONOREAGENT PROCEDURE "sample starter"

	BLANK	SAMPLE
WORKING REAGENT	1500 µl	1500 µl
DISTILLED WATER	100 µl	--
SAMPLE	--	100 µl

Mix, then incubate for 5 min. at 37°C. Measure the final absorbance of sample (EC) and blank (EBC) against water.

- BIREAGENT PROCEDURE "substrate starter"

	BLANK	SAMPLE
REAGENT R1	1000 µl	1000 µl
DISTILLED WATER	80 µl	--
SAMPLE	--	80 µl

Mix, incubate for 5 min. and then add:

	BLANK	SAMPLE
REAGENT R2	250 µl	250 µl

Mix, then incubate for 5 min. at 37°C. Measure the final absorbance of sample (EC) and blank (EBC) against water.

CALCULATION

Total bilirubin [mg/dl] = (EC - EBC) x 14.5
[µmol/l] = (EC - EBC) x 248

CONVERSION FACTOR

Bilirubin [mg/dl] x 17.0 = Bilirubin [µmol/l]

ANALYTICAL PERFORMANCES

Interferences

Glycerides do not interfere up to concentration of 500 mg/dl.
Hemoglobin does not interfere up to concentration of 1000 mg/dl.

Linearity

Reaction is linear between 0.05-25 mg/dl (0.85-425 µmol/l). Samples with values exceeding 25 mg/dl must be diluted with saline solution. Multiply, then, the result for diluting factor.

"Intra Assay" precision (within-Run)

Determined on 30 samples for each control (L-N-H) (Low-Normal-High).
Results:

	[mg/dl]	L = 0.43	N = 1.66	H = 3.55
MEAN				
S.D.		0.02	0.12	0.32
C.V.%		5.12	7.25	9.09

"Inter Assay" precision (between-Run)

Determined on 15 samples for each control (L-N-H) for 3 days.
Results:

	[mg/dl]	L = 0.84	N = 1.67	H = 4.09
MEAN				
S.D.		0.03	0.07	0.16
C.V.%		4.08	4.26	4.00

Analytical sensitivity

The test sensitivity in terms of detection limit is 0.05 mg/dl (0.85 µmol/l).

Correlation

A study based comparing this method with a similar method on 21 samples has given a correlating factor **r = 1.00**

BIBLIOGRAPHY

Jendrassik, L., Gróf, P., Biochem. Z., 297, 81 (1938).
Jakobs, D.S., Kasten, Jr. B.L., Demmott, W.R., Wolfson, W.L.: "Laboratory Test Handbook", Lexi-Comp and Williams & Wilkins Ed. 2nd Edition (1990).
Kaplan, L.A., Pesce, A.J.: "Clinical Chemistry", Mosby Ed. (1996).

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer