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Fabbricante

Kit for measurement of direct bilirubin in serum or plasma Colorimetric method Jendrassik-Grof modified

PRINCIPLE

Direct (or conjugate) bilirubin reacts, in acid environment, with diazotised sulphaniic acid giving a coloured azocompound (azobilirubin) whose colour intensity is proportional to the direct bilirubin concentration in the tested sample.

REAGENTS

R1 Sulphanilic acid 25.0 mmol/l; surface-active agents
R2 Nitrite sodium 0.32 mmol/l

SAMPLE

- Serum or EDTA-Na₂ plasma.

Note

- Do not use samples with haemolysis.
- Keep the samples far from light and heat because the bilirubin is a photosensitive pigment. Perform the samples as soon as possible.
- It's advisable to perform a blank sample for those samples which are hyperlipoeemic and not clear.

REFERENCE VALUES

Adults: up to 0.35 mg/dl (5.1 µmol/l)

References values are considered indicatives since each laboratory should establish references ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

STORAGE AND STABILITY

- Store the kit at 15-25°C. Do not freeze the reagents.
- After opening, the vials **R1** and **R2** are stable up to the expiry date if recapped immediately and protected from contamination, evaporation, direct light, and stored at the correct temperature.
- Working solution stability (**R1**+ **R2**): 21 days at 2-8°C.

PREPARATION OF REAGENTS

Reagents are liquid and ready to use. About using as monoreagent ("sample-starter" procedure) add to every 4 ml of **R1** reagent, 1 ml of **R2** reagent. Prepare the amount for use according to the stability of working solution.

NOTE

- The kit, according to this method, must be used in manual procedures. About automatic using follow specific applications.
- Evaluate carefully the results if working reagent absorbance is > 0.100 at 570 nm.
- Avoid direct light, contamination and evaporation.
- The volumes in the procedure can be changed proportionally.
- In case of complaint or quality control request, refer to the lot number on the package or the lot number on the singles vials.

AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

Suggested serum:

REF 20350 Precise Norm REF 20360 Precise Path

PRECAUTION IN USE

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of components is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances.
However the reagent should be handled with caution, according to good laboratory practice.
It's recommended to use disposable test tubes and glassware washed by HCl 1N and distilled water.

WASTE MANAGEMENT

Please refer to the local legal requirements.

PROCEDURE

Wavelength λ: 570 (550-580) nm
Working temperature 37°C
Optical path 1 cm
Reaction "end point"

- MONOREAGENT PROCEDURE "sample starter"

	BLANK	SAMPLE
WORKING REAGENT	1500 µl	1500 µl
DISTILLED WATER	100 µl	--
SAMPLE	--	100 µl

Mix, then incubate for 5' at 37°C. Measure the final absorbance of sample (EC) and blank (EBC) against water.

- BIREAGENT PROCEDURE "substrate starter"

	BLANK	SAMPLE
REAGENT R1	1000 µl	1000 µl
DISTILLED WATER	80 µl	--
SAMPLE	--	80 µl

Mix, incubate for 5' and then add:
REAGENT R2 250 µl 250 µl
Mix, then incubate for 5' at 37°C. Measure the final absorbance of sample (EC) and blank (EBC) against water.

CALCULATION

$$\begin{aligned} \text{Direct bilirubin [mg/dl]} &= (\text{EC} - \text{EBC}) \times 14.5 \\ [\mu\text{mol/l}] &= (\text{EC} - \text{EBC}) \times 248 \end{aligned}$$

CONVERSION FACTOR

$$\text{Bilirubin [mg/dl]} \times 17.0 = \text{Bilirubin } [\mu\text{mol/l}]$$

ANALYTICAL PERFORMANCES

Interferences

Triglycerides do not interfere up to concentration of 500 mg/dl.
Hemoglobin does not interfere up to concentration of 1000 mg/dl.

Linearity

Reaction is linear between 0.049-10 mg/dl (0.83-170 µmol/l). Samples with values exceeding 10 mg/dl must be diluted with saline solution. Multiply, then, the result for diluting factor.

"Intra Assay" precision (within-Run)

Determined on 30 samples for each control (N-H) (Normal-High).
Results:

MEAN	[mg/dl]	N = 0.30	H = 1.10
S.D.		0.03	0.06
C.V.%		8.20	5.72

"Inter Assay" precision (between-Run)

Determined on 15 samples for each control (N-H) for 3 days.
Results:

MEAN	[mg/dl]	N = 0.10	H = 1.25
S.D.		0.06	0.05
C.V.%		3.58	3.73

Analytical sensitivity

The test sensitivity in terms of detection limit is 0.049 mg/dl (0.83 µmol/l).

Correlation

A study based comparing this method with a similar method on 21 samples has given a correlating factor **r = 0.99**

BIBLIOGRAPHY

Jendrassik, L., Gróf, P., Biochem. Z., 297, 81 (1938).
Tiez, N.W., Textbook of Clinical Chemistry. W.B. Saunders Co (1983).
Fossati, P., Ponti, M., Principe, L., Tarengi, G., Clin. Chem., 35, 173 (1989).
Kaplan, LA, Pesce, A.J.: "Clinical Chemistry", Mosby Ed. (1996).

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer