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Fabbricante

For in vitro medical device
10118 - 4x20 ml
V10118 - 5x20+5x5 ml

Kit for measurement of lactate dehydrogenase in serum or plasma Kinetic optimized method SCE*

*Scandinavian Committee on Enzymes

PRINCIPLE

Lactate dehydrogenase (LDH) catalyzes the reaction between pyruvic acid and NADH. The NADH oxidization into NAD⁺ allows to determine the LDH in the tested sample.

REAGENTS

R1 Goods buffer pH 7.3 70.0 mmol/l; EDTA Na₂ 4.0 mmol/l;
pyruvic acid 1.2 mmol/l
R2 Goods buffer pH 7.3 70.0 mmol/l; NADH 1.20 mmol/l

SAMPLE

- Serum or heparinized plasma or EDTA plasma.

Note

- Do not use samples with haemolysis.
- The serum must be separated, as soon as possible, from corpuscular part.
- The LDH activity decreases in 8% within 3 days at 2-8°C.

REFERENCE VALUES

Adults [U/l]	25°C	30°C	37°C
	< 240	< 346	< 480

The LDH level in children (till 12 years old) is 10-15% higher than in the adults. Reference values are considered indicatives since each laboratory should establish reference ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

STORAGE AND STABILITY

- Store the kit at 2-8°C. Do not freeze the reagents.
- After opening, the vials **R1** and **R2** are stable up to the expiry date if recapped immediately and protected from contamination, evaporation, direct light, and stored at the correct temperature.
- Working solution stability (**R1**+ **R2**): 15 days at 2-8°C.

PREPARATION OF REAGENTS

Reagents are liquid and ready to use. About using as monoreagent ("sample-starter" procedure) add to every 4 ml of **R1** reagent, 1 ml of **R2** reagent. Keep out the reagents from refrigerator only for the use and recap them immediately.

NOTE

- The kit, according to this method, must be used in manual procedures. About automatic using follow specific applications.
- Evaluate carefully the results if working reagent absorbance is < 1.000 at 340 nm.
- Avoid direct light, contamination and evaporation.
- The volumes in the procedure can be changed proportionally.
- In case of complaint or quality control request, refer to the lot number on the package or the lot number on the singles vials.

QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

Suggested serum:

REF 20350 Precise Norm REF 20360 Precise Path

PRECAUTION IN USE

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of inactive components (preservatives, surfactants and so on) is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances.

However the reagent should be handled with caution, according to good laboratory practice.

Caution: the reagents contain Sodium Azide (0.095%) as preservative. Avoid swallowing and contacting with skin, eyes and mucous membranes.

WASTE MANAGEMENT

Please refer to the local legal requirements.

PROCEDURE

Wavelength λ : 340 nm (334 o 365 nm)
Working temperature 37°C (25 o 30°C)
Optical path 1 cm
Reaction kinetic (decreasing)

Bring the reagents at 15-25°C before using them.

- MONOREAGENT PROCEDURE "sample starter"

WORKING REAGENT	BLANK	SAMPLE
DISTILLED WATER	1000 μ l	1000 μ l
SAMPLE	20 μ l	--
	--	20 μ l

Mix, then incubate for 1' a 37°C. Measure the absorbance of sample (EC) against distilled water decrease per minute during 3 minutes. Calculate the absorbance variation $\Delta E/\text{min}$ from readings executed.

- BIREAGENT PROCEDURE "substrate starter"

REAGENT R1	BLANK	SAMPLE
DISTILLED WATER	800 μ l	800 μ l
SAMPLE	20 μ l	--
	--	20 μ l

Mix, incubate at 37°C for 1' and then add:

REAGENT R2	BLANK	SAMPLE
	200 μ l	200 μ l

Mix, then incubate for 1' a 37°C. Measure the absorbance of sample (EC) against distilled water decrease per minute during 3 minutes. Calculate the absorbance variation $\Delta E/\text{min}$ from readings executed.

CALCULATION

$$\text{LDH [U/l]} = \Delta E/\text{min} \times 8095$$

The factor and the reagent performances are related to 37°C and 340 nm.

ANALYTICAL PERFORMANCES

Interferences

Bilirubin does not interfere up to concentration of 40 mg/dl.
Triglycerides do not interfere up to concentration of 2000 mg/dl.
Ascorbate acid does not interfere up to concentration of 30 mg/dl.

Linearity

Reaction is linear up to an activity of 800 U/l (corresponding to $\Delta E/\text{min}$ of 0.099 at 340 nm). Samples with values exceeding 800 U/l must be diluted with saline solution. Multiply, then, the result for diluting factor.

"Intra-Assay" precision (within-Run)

Determined on 30 samples for each control (L-N-H) (Low-Normal-High).

Results:

MEAN (U/l)	L = 91.00	N = 184.90	H = 345.80
S.D.	3.36	5.35	8.10
C.V.%	3.70	2.89	2.34

"Inter-Assay" precision (between-run)

Determined on 15 samples for each control (L-N-H) for 3 days.

Results:

MEAN (U/l)	L = 109.47	N = 209.20	H = 385.07
S.D.	3.73	4.01	6.73
C.V.%	3.39	1.93	1.75

Analytical sensitivity

The test sensitivity in terms of detection limit is 17 U/l.

Correlation

A study based comparing this method with a similar method on 21 samples has given a correlating factor $r = 0.97$

BIBLIOGRAPHY

Kreutzer, H. H. et al., Clin. Chim. Acta, 9, 64, (1964).
Young, D.S. et al., Clin. Chem., 21 ID 432D, (1975).
Fischbach, F., Zawta, B., Klin. Lab., 38:555-61, (1992).
Kaplan, L.A., Pesce, A.: "Clinical Chemistry", Mosby Ed. (1996)

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer