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## Kit for measurement of alanine-aminotransferase in serum or plasma – Kinetic UV optimized method IFCC\*

\* International Federation of Clinical Chemistry and Laboratory Medicine

### PRINCIPLE

The enzyme alanine-aminotransferase (ALT) (or glutamic-pyruvic transaminase/GPT) catalyzes reaction between alpha-ketoglutarate and L-alanine giving glutamate and pyruvic acid. In presence of lactate dehydrogenase (LDH), pyruvic acid reacts with NADH giving lactic acid and NAD<sup>+</sup>. The absorbance variation is proportional to the ALT activity. In the reagent is contained LDH to convert the endogenous pyruvate into lactate during the preincubation.

### REAGENTS

- R1** Goods buffer pH 7.5 80.0 mmol/l; L-alanine 500.0 mmol/l;  
 LDH ≥ 1500 U/l;  
**R2** Goods buffer pH 7.5 80.0 mmol/l; alpha-ketoglutarate 50.0 mmol/l;  
 NADH ≥ 1.18 mmol/l

### SAMPLE

- Serum-heparinized plasma or EDTA plasma.

### Note

- Do not use samples with haemolysis because this one could cause wrongly positive results.
- The ALT activity tends to decrease (< 10%) after 3 days at 2-8 °C.

### REFERENCE VALUES

Serum - plasma [U/l]	37°C
Women	31
Men	41

References values are considered indicatives since each laboratory should establish references ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

### STORAGE AND STABILITY

- Store the kit at 2-8 °C. Do not freeze the reagents.
- After opening, the vials **R1** and **R2** are stable up to the expiry date if recapped immediately and protected from contamination, evaporation, direct light, and stored at the correct temperature.
- Working solution stability (**R1**+ **R2**): 15 days at 2-8 °C.

### PREPARATION OF REAGENTS

Reagents are liquid and ready to use. About using as monoreagent ("sample-starter" procedure) add to every 4 ml of **R1** reagent, 1 ml of **R2** reagent. Keep out the reagents from refrigerator only for the use and recap them immediately.

### NOTE

- The kit, according to this method, must be used in manual procedures. About automatic using follow specific applications.
- Evaluate carefully the results if working reagent absorbance is < 1.000 at 340 nm.
- Avoid direct light, contamination and evaporation.
- The volumes in the procedure can be changed proportionally.
- In case of complaint or quality control request, refer to the lot number on the package or the lot number on the singles vials.

### AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

### QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

### Suggested serum:

**REF** 20350 Precise Norm **REF** 20360 Precise Path

### PRECAUTION IN USE

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of inactive components (preservatives, surfactants and so on) is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances.

However the reagent should be handled with caution, according to good laboratory practice.

**Caution:** the reagents contain Sodium Azide (0.095%) as preservative. Avoid swallowing and contacting with skin, eyes and mucous membranes.

### WASTE MANAGEMENT

Please refer to the local legal requirements.

### PROCEDURES

Wavelength	λ: 340 nm
Working temperature	37 °C
Optical path	1 cm
Reaction	kinetic (decreasing)

Bring the reagents at 15-25 °C before using them.

### - MONOREAGENT PROCEDURE "sample starter"

	BLANK	SAMPLE
WORKING REAGENT	1000 µl	1000 µl
DISTILLED WATER	100 µl	--
SAMPLE	--	100 µl

Mix, then incubate for 1' a 37 °C. Measure the initial absorbance and at 1 minute intervals thereafter for 3 minutes. Calculate the difference between consecutive absorbances, and the average absorbance difference per minute ( $\Delta E/min$ ).

### - BIREAGENT PROCEDURE "substrate starter"

	BLANK	SAMPLE
REAGENT R1	800 µl	800 µl
DISTILLED WATER	100 µl	--
SAMPLE	--	100 µl

Mix, incubate at 37 °C for 1' and then add:

	BLANK	SAMPLE
REAGENT R2	200 µl	200 µl

Mix, then incubate for 1' a 37 °C. Measure the initial absorbance and at 1 minute intervals thereafter for 3 minutes. Calculate the difference between consecutive absorbances, and the average absorbance difference per minute ( $\Delta E/min$ ).

### CALCULATION

$$ALT [U/l] = \Delta E/min \times 1746$$

### ANALYTICAL PERFORMANCES

#### Interferences

Bilirubin does not interfere up to concentration of 40 mg/dl.  
 Triglycerides do not interfere up to concentration of 2000 mg/dl.  
 Hemoglobin does not interfere up to concentration of 400 mg/dl.  
 Ascorbate acid does not interfere up to concentration of 30 mg/dl.

#### Linearity

Reaction is linear up to 440 U/l. Samples with values exceeding this range must be diluted with saline solution. Multiply, then, the result for diluting factor.

#### "Intra-Assay" precision (within-Run)

Determined on 30 samples for each control (L-N-H) (Low-Normal-High).

Results:

	MEAN (U/l)	L = 25.4	N = 53.2	H = 111.0
S.D.	1.07	1.14	3.00	
C.V.%	4.22	2.13	2.70	

#### "Inter-Assay" precision (between-run)

Determined on 15 samples for each control (L-N-H) for 3 days.

Results:

	MEAN (U/l)	L = 26.93	N = 54.11	H = 108.62
S.D.	0.74	1.14	1.63	
C.V.%	2.75	2.10	1.51	

#### Analytical sensitivity

The test sensitivity in terms of detection limit is 4 U/l.

#### Correlation

A study based comparing this method with a similar method on 21 samples has given a correlating factor  $r = 0.99$

#### BIBLIOGRAPHY

- Karmen A. et al.: J. Clin. Invest. 34, 126 (1955).
- Young DS, Pestaner LC and Gibberman V: Clin. Chem., 21(5), 1D-432D (1975).
- Lorentz K, Röhle G, Siekmann L.: DG Klinische Chemie Mitteilungen, 26, 190 (1995).
- Bergmeyer HU, Horder M, Rej R.: International Federation of Clinical Chemistry (IFCC) Scientific Committee. J. Clin. Chem. Clin. Biochem., 24, 497 (1986).
- Kaplan LA, Pesce AJ: "Clinical Chemistry", Mosby Ed. (1996).

### SYMBOLS

-  Read instruction for use
-  CE mark (requirement of 98/79 regulation)
-  Storing temperature limits
-  In vitro medical device