IVD Per uso diagnostico in vitro
REF 10139 - 5x50 ml
 V10139 - 5x20+5x5 ml
 H10139 - 4x80+2x40 ml

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Fabbricante

IVD For in vitro medical device
REF 10139 - 5x50 ml
 V10139 - 5x20+5x5 ml
 H10139 - 4x80+2x40 ml

Kit for measurement of aspartate-aminotransferase in serum or plasma - Kinetic UV optimized method IFCC*

* International Federation of Clinical Chemistry and Laboratory Medicine

PRINCIPLE

The enzyme aspartate-aminotransferase (AST) (or glutamic-oxaloacetic transaminase/GOT) catalyzes reaction between alpha-ketoglutarate and L-aspartate giving glutamate and oxaloacetate. In presence of malate dehydrogenase (MDH), oxaloacetate reacts with NADH giving malate and NAD⁺. The change of absorbance measured at 340 nm is proportional to the AST activity of the sample.

REAGENTS

R1 Goods buffer pH 7.8 80.0 mmol/l; L-aspartate 240.0 mmol/l; LDH ≥ 1800 U/l; MDH ≥ 800 U/l

R2 Goods buffer pH 7.8 80.0 mmol/l; alpha-ketoglutarate 65.0 mmol/l; NADH ≥ 1.18 mmol/l

SAMPLE

- Serum-heparinized plasma or EDTA plasma.

Note

- Do not use samples with haemolysis because this one could cause results wrongly positive. The anticoagulants containing ammonium salt (es. ammonium heparinate) mustn't be used.
- The AST activity tends to decrease (< 8%) after 3 days at 2-8°C.

REFERENCE VALUES

Serum - plasma [U/l]	37°C
Women	≤ 31
Men	≤ 38

References values are considered indicatives since each laboratory should establish references ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

STORAGE AND STABILITY

- Store the kit at 2-8°C. Do not freeze the reagents.
- After opening, the vials **R1** and **R2** are stable up to the expiry date if recapped immediately and protected from contamination, evaporation, direct light, and stored at the correct temperature.
- Working solution stability (**R1+ R2**): 30 days at 2-8°C.

PREPARATION OF REAGENTS

Reagents are liquid and ready to use. About using as monoreagent ("sample-starter" procedure) add to every 4 ml of **R1** reagent, 1 ml of **R2** reagent. Keep out the reagents from refrigerator only for the use and recap them immediately.

NOTE

- The kit, according to this method, must be used in manual procedures. About automatic using follow specific applications.
- Evaluate carefully the results if working reagent absorbance is < 1.000 at 340 nm.
- Avoid direct light, contamination and evaporation.
- The volumes in the procedure can be changed proportionally.
- In case of complaint or quality control request, refer to the lot number on the package or the lot number on the singles vials.

AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

Suggested serum:

REF 20350 Precise Norm **REF** 20360 Precise Path

PRECAUTION IN USE

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of inactive components (preservatives, surfactants and so on) is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances.

However the reagent should be handled with caution, according to good laboratory practice.

Caution: the reagents contain Sodium Azide (0.095%) as preservative. Avoid swallowing and contacting with skin, eyes and mucous membranes.

WASTE MANAGEMENT

Please refer to the local legal requirements.

PROCEDURES

Wavelength λ: 340 nm
 Working temperature 37°C
 Optical path 1 cm
 Reaction kinetic (decreasing)

Bring the reagents at 15-25°C before using them.

- MONOREAGENT PROCEDURE "sample starter"

WORKING REAGENT	BLANK	SAMPLE
DISTILLED WATER	1000 µl	1000 µl
SAMPLE	100 µl	--
	--	100 µl

Mix, then incubate for 1' a 37°C. Measure the absorbance of sample (EC) against distilled water decrease per minute during 3 minutes. Calculate the average absorbance difference per minute (ΔE/min).

- BIREAGENT PROCEDURE "substrate starter"

REAGENT R1	BLANK	SAMPLE
DISTILLED WATER	800 µl	800 µl
SAMPLE	100 µl	--
	--	100 µl

Mix, incubate at 37°C for 1' and then add:

REAGENT R2	BLANK	SAMPLE
	200 µl	200 µl

Mix, then incubate for 1' a 37°C. Measure the absorbance of sample (EC) against distilled water decrease per minute during 3 minutes. Calculate the average absorbance difference per minute (ΔE/min).

CALCULATION

$$AST [U/l] = \Delta E/min \times 1746$$

ANALYTICAL PERFORMANCES

Interferences

Bilirubin does not interfere up to concentration of 40 mg/dl.
 Triglycerides do not interfere up to concentration of 2000 mg/dl.
 Hemoglobin interferes also at minimum concentrations.
 Ascorbate acid does not interfere up to concentration of 30 mg/dl.

Linearity

Reaction is linear up to a concentration of 440 U/l. Samples with values exceeding this range must be diluted with saline solution. Multiply, then, the result for diluting factor.

"Intra-Assay" precision (within-Run)

Determined on 30 samples for each control (L-N-H) (Low-Normal-High).

Results:

MEAN (U/l)	L = 16.20	N = 35.00	H = 179.50
S.D.	0.66	1.03	4.93
C.V.%	4.10	2.96	2.74

"Inter-Assay" precision (between-run)

Determined on 15 samples for each control (L-N-H) for 3 days.

Results:

MEAN (U/l)	L = 17.58	N = 35.51	H = 174.22
S.D.	0.59	0.78	2.86
C.V.%	3.37	2.20	1.64

Analytical sensitivity

The test sensitivity in terms of detection limit is 0.3 U/l.

Correlation

A study based comparing this method with a similar method on 21 samples has given a correlating factor **r = 0.99**

BIBLIOGRAPHY

- Karmen A. et al.: J. Clin. Invest. 34, 126 (1955).
- Young DS, Pestaner LC and Gibberman V: Clin. Chem., 21(5),1D-432D (1975).
- Lorentz K, Röhle G, Siekmann L.: DG Klinische Chemie Mitteilungen, 26, 190 (1995).
- Bergmeyer HU, Horder M, Rej R.: International Federation of Clinical Chemistry (IFCC) Scientific Committee. J. Clin. Chem. Clin. Biochem., 24, 497 (1986).
- Kaplan LA, Pesce AJ: "Clinical Chemistry", Mosby Ed. (1996).

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer