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Fabbricante

Kit for measurement of chlorides in serum, plasma and urine Colorimetric method mercurious thiocyanate

PRINCIPLE

Chlorides react with mercurious thiocyanate releasing thiocyanate ions which form, together with iron trivalent, a coloured compound whose colour intensity is proportional to chlorine concentration.

REAGENTS

R1 Mercurious thiocyanate 1.0 mmol/l; iron nitrate 10.0 mmol/l;
nitric acid 62.0 mmol/l; nitrate mercurious 0.15 mmol/l
Standard chlorides 100 mEq/l (100 mmol/l)

SAMPLE

- Serum or plasma. Diluted urine 1:2.

Note

- Do not use samples with haemolysis.
- Samples strongly jaundiced, lipoemic or with haemolysis require the use of a blank sample against distilled water.
- The chlorides are stable in the samples up to 7 days at 2-8°C.

REFERENCE VALUES

Serum - plasma	98 - 110 mEq/l
Urine	170 - 250 mEq/24h

References values are considered indicatives since each laboratory should establish references ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

STORAGE AND STABILITY

- Store the kit at 15-25°C.
- After opening, the reagent and the standard are stable up to the expiry date if recapped immediately and protected from contamination, evaporation, direct light, and stored at the correct temperature.

PREPARATION OF REAGENTS

Reagents are liquid and ready to use.

NOTE

- The kit, according to this method, must be used in manual procedures. About automatic using follow specific applications.
- Evaluate carefully the results if reagent R1 absorbance is > 0.300 at 480 nm.
- The volumes in the procedure can be changed proportionally.
- In case of complaint or quality control request, refer to the lot number on the package or the lot number on the singles vials.

AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

Suggested serum:

REF 20350 Precise Norm REF 20360 Precise Path

PRECAUTION IN USE

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of components is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances. However the reagent should be handled with caution, according to good laboratory practice.

Caution: the reagents contain Sodium Azide (0.095%) as preservative. Avoid swallowing and contacting with skin, eyes and mucous membranes.

WASTE MANAGEMENT

Please refer to the local legal requirements.

PROCEDURE

Wavelength λ : 480 (460-500) nm
Working temperature 37°C (15-25°C)
Optical path 1 cm
Reaction "end point"

- PROCEDURE

	BLANK	STD	SAMPLE
REAGENT R1	1000 μ l	1000 μ l	1000 μ l
DISTILLED WATER	10 μ l	--	--
SAMPLE	--	--	10 μ l
STANDARD	--	10 μ l	--

Mix, then incubate for 5' at 37°C (15-25°C). Measure the absorbance of sample (EC) and standard (ESTD) against the reagent blank.

CALCULATION

Chlorides [mEq/l] = EC/ESTD x Conc. STD

Diluted urines: multiply the result for diluting factor.

ANALYTICAL PERFORMANCES

The reagent performances are related to 37°C and 480 nm.

Interferences

Bilirubin interferes also at minimum concentrations.
Triglycerides interfere also at minimum concentrations.
Hemoglobin interferes also at minimum concentrations.

Linearity

Reaction is linear between 28-130 mEq/l. Samples with values exceeding 130 mEq/l must be diluted with distilled water. Multiply, then, the result for diluting factor.

"Intra-Assay" precision (within-Run)

Determined on 30 samples for each control (L-N-H) (Low-Normal-High).

Results:

MEAN	[mEq/l]	L = 63.3	N = 95.6	H = 103.0
S.D.		0.80	0.68	0.98
C.V.%		1.27	0.71	0.95

"Inter-Assay" precision (between-Run)

Determined on 15 samples for each control (L-H-N) for 3 days.

Results:

MEAN	[mEq/l]	L = 59.16	N = 90.67	H = 103.76
S.D.		0.70	0.91	0.89
C.V.%		1.18	1.00	0.85

Analytical sensitivity

The test sensitivity in terms of detection limit is 28 mEq/l.

Correlation

A study based comparing this method with a similar method on 21 samples has given a correlating factor **r = 0.94**

The urine performances are available on request.

BIBLIOGRAPHY

Schoenfeld R.G., loewellw C.S.: Clinical Chemistry (1964), 10, 533.
Kaplan LA, Pesce AJ: "Clinical Chemistry", Mosby Ed. (1996).

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer