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Kit for measurement of uric acid in serum, plasma and urine Colorimetric enzymatic method Uricase-POD-PAP

PRINCIPLE

Uric acid is converted by uricase and hydrogen peroxide which, under the catalytic influence of peroxidase (POD), oxidizes compound, reacts with 4-aminophenazone and 3,5-dichlorophenol-sulphonate giving a red coloured compound, whose colour intensity is directly proportional to the uric acid concentration in the tested sample.

REAGENTS

- R1** TRIS buffer 60.0 mmol/l; 2,4-dichlorophenol-sulphonate 2.0 mmol/l
R2 Uricase ≥ 90 U/l; peroxidase ≥ 1000 U/l;
4-aminophenazone 0.8 mmol/l; no reactive stabilizers
R3 Uric acid standard 6 mg/dl (297.5 μ mol/l)

SAMPLE

- Serum or heparinized plasma. Diluted urine 1:10.

Note

- Do not use samples with haemolysis. Do not use oxalate as anticoagulant. It's not advisable using EDTA and fluoride because they could cause positive interferences.
- The acid uric in the serum and in the plasma is stable for 3 days if stand at 2-8°C or 6 months at -20°C.

REFERENCE VALUES

Serum - plasma	Men 3.4 - 7.0 mg/dl (202 - 416 μ mol/l)
	Women 2.4 - 5.7 mg/dl (143 - 339 μ mol/l)
Urine	250 - 750 mg/24h (1.49 - 4.46 mmol/24h)

References values are considered indicatives since each laboratory should establish references ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

PREPARATION OF REAGENTS

Dissolve the content of **R2** vial with **R1** solution. Keep out the reagents from refrigerator only for the use and recap them immediately.

STORAGE AND STABILITY

- Store the kit at 2-8°C. Do not freeze the reagents.
- Reconstituted reagent stability: 60 days at 2-8°C.

NOTE

- The kit, according to this method, must be used in manual procedures. About automatic using follow specific applications.
- Evaluate carefully the results if working reagent absorbance is > 0.300 at 510 nm.
- Avoid direct light, contamination and evaporation.
- The volumes in the procedure can be changed proportionally.
- In case of complaint or quality control request, refer to the lot number on the package or the lot number on the singles vials.

AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

Suggested serum:

REF 20350 Precise Norm REF 20360 Precise Path

PRECAUTION IN USE

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of components is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances.

However the reagent should be handled with caution, according to good laboratory practice.

Caution: the reagents contain Sodium Azide (0.095%) as preservative. Avoid swallowing and contacting with mucous membranes.

WASTE MANAGEMENT

Please refer to the local legal requirements.

PROCEDURE

Wavelength	λ : 510 (500-550) nm
Working temperature	37°C
Optical path	1 cm
Reaction	"end point" (increasing)

Bring the reagents at 15-25°C before using them.

	BLANK	STD	SAMPLE
SAMPLE	--	--	25 μ l
DISTILLED WATER	25 μ l	--	--
STANDARD	--	25 μ l	--
REAGENT	1000 μ l	1000 μ l	1000 μ l

Mix, then incubate for 5' at 37°C. Measure the absorbance of sample (EC) and standard (ESTD) against the reagent blank. The colour obtained is stable for about 60' if stand at 15-25°C and protected from light.

CALCULATION

$$\text{Uric acid [mg/dl] o } [\mu\text{mol/l}] = \text{EC/ESTD} \times \text{Conc. STD}$$

Diluted urines: multiply the result for diluting factor.

The factor and the reagent performances are related to 510 nm and 37°C.

CONVERSION FACTOR

$$\text{Uric acid [mg/dl]} \times 59.48 = \text{Uric acid } [\mu\text{mol/l}]$$

ANALYTICAL PERFORMANCES

Interferences

- Bilirubin does not interfere up to concentration of 20 mg/dl.
- Triglycerides do not interfere up to concentration of 2000 mg/dl.
- Hemoglobin does not interfere up to concentration of 50 mg/dl.
- Ascorbate acid does not interfere up to concentration of 30 mg/dl.

Linearity

Reaction is linear between 0.28-20 mg/dl (16.6-1189.6 μ mol/l). Samples with values exceeding 20 mg/dl must be diluted with saline solution. Multiply, then, the result for diluting factor.

"Intra-Assay" precision (within-Run)

Determined on 30 samples for each control (L-N-H) (Low-Normal-High).

Results:

MEAN	(mg/dl)	L = 2.46	N = 4.90	H = 8.16
S.D.		0.06	0.12	0.20
C.V.%		2.53	2.40	2.40

"Inter-Assay" precision (between-Run)

Determined on 15 samples for each control (L-H-N) for 3 days.

Results:

MEAN	(mg/dl)	L = 2.48	N = 5.07	H = 8.17
S.D.		0.07	0.09	0.15
C.V.%		2.68	1.73	1.85

Analytical sensitivity

The test sensitivity in terms of detection limit is 0.28 mg/dl (16.6 μ mol/l).

Correlation

A study based comparing this method with a similar method on 21 samples has given a correlating factor $r = 0.99$

BIBLIOGRAPHY


- Trivedi R.C., Rebar L., Berka E., Strong L., Clin Chem (1978), 24, 1908.
 Barham D, Trinder P: Analyst, 97 142 (1972).
 Fossati P, Prencipe L, Berti G: Clin. Chem., 26(2) 227 (1980).
 Kaplan LA, Pesce AJ: "Clinical Chemistry", Mosby Ed. (1996).

SYMBOLS

 Read instruction for use

 CE mark (requirement of 98/79 CE regulation)

 Storing temperature limits

 In vitro medical device

 Producer