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Fabbricante

Kit for measurement of triglycerides in serum or plasma Colorimetric enzymatic method GPO-PAP

PRINCIPLE

Triglycerides are hidrolized, in presence of lypoproteinlipase (LPL), into fat acid and glycerol which is transformed, by glycerolchinase (GK), ATP and glicerol-3-P-oxidase (GPO), into diidroxiaceton-phosphate and H₂O₂. The hydrogen peroxid catalyzed from peroxidase (POD) reacts with 4- aminophenazone and 4-phenol-chloride giving a coloured compound whose intensity is proportional to the concentration of triglycerides in the sample.

REAGENTS

- R1** PIPES buffer 40.0 mmol/l; 4-chlorophenol 4.0 mmol/l.
R2 Lipoproteinlipase ≥ 800 U/l; glicerolkinase ≥ 1200 U/l;
 glicerol-3-phosphate-oxidase ≥ 2500 U/l; peroxidase ≥ 400 U/l;
 4-aminophenazone 0.7 mmol/l; ATP 0.18 mmol/l; no reactive stabilizers.
R3 Glycerol standard equivalent to 200 mg/dl (2.26 mmol/l) of triglycerides.

SAMPLE

- Serum-heparinized plasma or EDTA plasma.

Note

- In the samples stored at 2-8°C and at -20°C the triglycerides are stable, respectively, up to 3 days and 12 months.
- It's advisable, in presence of strongly lipoemic, jaundiced or turbid serum, to prepare a blank of the sample using saline solution.

REFERENCE VALUES

Men	60-165 mg/dl (0.68 - 1.86 mmol/l)
Women	40-140 mg/dl (0.45 - 1.58 mmol/l)

It has been noticed that the values obtained using the plasma as sample are, from 2% to 4%, lower than values obtained using the serum.

References values are considered indicatives since each laboratory should establish references ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

PREPARATION OF REAGENTS

Dissolve the content of **R2** vial with **R1** solution.
 Keep out the reagents from refrigerator only for the use and recap them immediately.

STORAGE AND STABILITY

- Store the kit at 2-8°C. Do not freeze the reagents.
- Reconstituted reagent stability: 40 days at 2-8°C.

NOTE

- The kit, according to this method, must be used in manual procedures. About automatic using follow specific applications.
- Evaluate carefully the results if working reagent absorbance is >0.300 at 510nm.
- Avoid direct light, contamination and evaporation.
- The volumes in the procedure can be changed proportionally.
- In case of complaint or quality control request, refer to the lot number on the package or the lot number on the singles vials.

AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

Suggested serum:

REF 20350 Precise Norm **REF** 20360 Precise Path

PRECAUTION IN USE

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of components is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances.

However the reagent should be handled with caution, according to good laboratory practice.

Caution: the reagents contain Sodium Azide (0.095%) as preservative. Avoid swallowing and contacting with skin, eyes and mucous membranes.

WASTE MANAGEMENT

Please refer to the local legal requirements.

PROCEDURE

Wavelength	λ: 510 (500-550) nm
Working temperature	37°C
Optical path	1 cm
Reaction	"end point"

Bring the reagents at 15-25°C before using them.

	BLANK	STD	SAMPLE
WORKING REAGENT	1000 µl	1000 µl	1000 µl
DISTILLED WATER	10 µl	--	--
SAMPLE	--	--	10 µl
STANDARD	--	10 µl	--

Mix, then incubate for 5' at 37°C. Measure the absorbance of sample (EC) and standard (ESTD) against the reagent blank. The colour obtained is stable for about 90' if stand at 15-25°C and protected from light.

CALCULATION

$$\text{Triglycerides [mg/dl]} \text{ o [mmol/l]} = \text{EC/ES} \times \text{Conc. STD}$$

CONVERSION FACTOR

$$\text{Triglycerides [mg/dl]} \times 0.01126 = \text{Triglycerides [mmol/l]}$$

ANALYTICAL PERFORMANCES

Interferences

Bilirubin does not interfere up to concentration of 40 mg/dl.
 Ascorbate acid does not interfere up to concentration of 6 mg/dl.
 Hemoglobin does not interfere up to concentration of 2.5 mg/dl.

Linearity

Reaction is linear between 7.3-1000 mg/dl (0.082-11.26 mmol/l). Samples with values exceeding 1000 mg/dl must be diluted with saline solution. Multiply, then, the result for diluting factor.

"Intra Assay" precision (within-Run)

Determined on 30 samples for each control (L-N-H) (Low-Normal-High).
 Results:

MEAN	(mg/dl)	L = 50.40	N = 95.40	P = 212.00
S.D.		1.16	2.09	5.17
C.V.%		2.31	2.20	2.44

"Inter Assay" precision (between-Run)

Determined on 15 samples for each control (L-N-H) for 3 days.
 Results:

MEAN	(mg/dl)	L = 51.07	N = 106.08	P = 230.76
S.D.		1.35	1.85	5.06
C.V.%		2.64	1.75	2.19

Analytical sensitivity

The test sensitivity in terms of detection limit is 7.3 mg/dl (0.082 mmol/l).

Correlation

A study based comparing this method with a similar method on 21 samples has given a correlating factor **r = 1.00**

BIBLIOGRAPHY

- Trinder P., Ann Clin Biochem (1969), 6, 24.
 Bowie L., Gochman N. Clin.Chem (1973), 19, 56.
 Bucolo G, David M: Clin. Chem., 19, 476 (1973).
 McGowan MW, Artiss JD, Standbergh DR, Zak B: Clin. Chem., 29, 538 (1983).
 Kaplan LA, Pesce AJ: "Clinical Chemistry", Mosby Ed. (1996).

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer