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BIBLIOGRAFIA

- NCCLS Document: "Procedures for the Collection of Arterial Blood Specimens; Approved Standard - Third Edition (1999)".
- Kaplan, L.A., Pesce, A.J.: "Clinical Chemistry", Mosby Ed. (1996).
- EU-Dir 1999/11 Commission Directive of 8 March 1999 adapting to technical progress the principles of Good Laboratory Practice as specified in Council Directive 87/18/EEC.
- Jakobs, D.S., Kasten, Jr., B.L., Demmott, W.R., Wolfson, W.L.: "Laboratory Test Handbook", Lexi-Comp and Williams & Wilkins Ed. (2nd Edition - 1990).
- M. Okada et al.: Low-density lipoprotein cholesterol can be chemically measured: A new superior method, J. Lab. Clin. Med., in press.
- Bachorik, P.S., Ross, J.W., "national Cholesterol Education Program Recommendations for measurement of low-density lipoprotein cholesterol: executive summary". Clin. Chem. 41: 1414-1420 (1995).

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Fabbricante

Kit for direct enzymatic measurement of LDL cholesterol in serum Elimination method

PRINCIPLE

This assay can directly measure LDL cholesterol level in serum thanks to an elimination method that consists of two specific steps. In the first step, chylomicron, VLDL and HDL are eliminated under specific conditions so cholesterol is derived only from LDL. In fact these fractions are oxidized (cholesterol esterase, cholesterol oxidase) to cholestenone and H₂O₂, that is subsequently degraded by catalase. In the second step, after various enzymatic reactions and in the presence of specific surfactants, remaining LDL cholesterol can be measured specifically as a colour formation (quinone pigment), which intensity is proportional to the concentration of LDL cholesterol contained in the sample.

REAGENTS

- R1** Mes buffer (pH 6.3); cholesterol esterase ≥ 1500 U/l;
 colesterolo oxidase ≥ 1500 U/l; peroxidase ≥ 1300 ppg U/l;
 4-aminoantipyrine ≥ 0.1%; ascorbate oxidase ≥ 3000 U/l; preservatives
- R2** Mes buffer (pH 6.3); N,N-bis(4-sulfobutyl)-toluidine ≥ 1.0 mM; preservatives

SAMPLE

- Serum.

Note

- The LDL cholesterol is stable in the samples up to 3 days at 2-8°C.

REFERENCE VALUES

	Recommended levels	Moderate risk	High risk
LDL cholesterol:	< 130 mg/dl (< 3.37 mmol/l)	130 - 159 mg/dl (3.37 - 4.12 mmol/l)	≥ 160 mg/dl (≥ 4.14 mmol/l)

References values are considered indicatives since each laboratory should establish references ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

STORAGE AND STABILITY

- Store the kit at 2-8°C. Do not freeze the reagents.
- After opening, the vials **R1** and **R2** are stable up to the expiry date if recapped immediately and protected from contamination, evaporation, direct light, and stored at the correct temperature.
- Working solution stability (**R1**+ **R2**): 30 days at 2-8°C.

PREPARATION OF REAGENTS

Reagents are liquid and ready to use.

NOTE

- The kit, according to this method, must be used in manual procedures. About automatic using follow specific applications.
- The reagents must be limpid; do not use if turbid.
- A possible green coloration of reagent R1 does not affect the reagent performances.
- Do not use calibration and control materials that contain sodium azide as it is an inhibitor for the catalase enzyme contained in reagent R1.
- Avoid direct light, contamination and evaporation.
- The volumes in the procedure can be changed proportionally.
- In case of complaint or quality control request, refer to the lot number on the package or the lot number on the singles vials.

AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

Suggested serum:

REF 20120 HDL/LDL calibrator HDL /LDL control lipidic

PRECAUTION IN USE

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of components is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances.

However the reagent should be handled with caution, according to good laboratory practice.

Caution: the reagents contain Sodium Azide (0.095%) as preservative. Avoid swallowing and contacting with skin, eyes and mucous membranes.

WASTE MANAGEMENT

Please refer to the local legal requirements.

PROCEDURE

Wavelength	λ: 600 (570-600) nm
Working temperature	37°C
Optical path	1 cm
Reaction	"end point" (increasing)

Bring the reagents at 37°C before using them.

- BIREAGENT PROCEDURE "substrate starter"

	BLANK	SAMPLE	CALIBRATOR
DISTILLED WATER	3 µl	--	--
SAMPLE	--	3 µl	--
CALIBRATOR	--	--	3 µl
REAGENT R1	300 µl	300 µl	300 µl
Mix, incubate at 37°C for 5' and then add:			
	BLANK	SAMPLE	CALIBRATOR
REAGENT R2	100 µl	100 µl	100 µl
Mix and incubate at 37°C for 5'. Read the absorbance of sample (AC) and calibrator (Acal) against the reagent blank.			

CALCULATION

$$\text{LDL cholesterol [mg/dl] o [mmol/l]} = \text{AC} / \text{Acal} \times \text{CAL}^*$$

*CAL = concentration of LDL cholesterol of the calibrator

CONVERSION FACTOR

LDL cholesterol [mg/dl] x 0.0259 = LDL cholesterol [mmol/l]

ANALYTICAL PERFORMANCES

Interferences

Bilirubin does not interfere up to concentration of 20 mg/dl.
 Triglycerides do not interfere up to concentration of 1293 mg/dl.
 Hemoglobin does not interfere up to concentration of 500 mg/dl.
 Ascorbate acid does not interfere up to concentration of 50 mg/dl.
 Gamma-globulins do not interfere up to concentration of 5000 mg/dl.

Linearity

Reaction is linear up to a concentration of 992 mg/dl with a range of 0.278-992 mg/dl. Samples with values exceeding 992 mg/dl must be diluted with saline solution. Multiply, then, the result for diluting factor.

"Intra-Assay" precision (within-Run)

Determined on 20 samples for each control (L-N-H) (Low-Normal-High).

Results:

MEAN	[mg/dl]	L = 98.1	N = 146.5	H = 209.8
S.D.		0.72	0.96	1.31
C.V.%		0.73	0.66	0.62

"Inter-Assay" precision (between-Run)

Determined on 40 samples for each control (L-N-H) for 10 days.

Results:

MEAN	[mg/dl]	L = 98.1	N = 142.7	H = 207.3
S.D.		2.2	2.8	3.6
C.V.%		2.27	1.95	1.73

Analytical sensitivity

The test sensitivity in terms of detection limit is 0.278 mg/dl.

Correlation

A study based comparing this method with a similar method on 92 samples has given a correlating factor **r = 0.97**

BIBLIOGRAFIA


NCCLS Document: "Procedures for the Collection of Arterial Blood Specimens; Approved Standard – Third Edition (1999)".
 Kaplan, L.A., Pesce, A.J.: "Clinical Chemistry", Mosby Ed. (1996).
 EU-Dir 1999/11 Commission Directive of 8 March 1999 adapting to technical progress the principles of Good Laboratory Practice as specified in Council Directive 87/18/EEC.
 Jakobs, D.S., Kasten, Jr., B.L., Demmott, W.R., Wolfons, W.L.: "Laboratory Test Handbook", Lexi-Comp and Williams & Wilkins Ed. (2nd Edition - 1990).
 M. Okada et al.: Low-density lipoprotein cholesterol can be chemically measured: A new superior method, J. Lab. Clin. Med., in press.
 Bachorik, P.S., Ross, J.W., "national Cholesterol Education Program Recommendations for measurement of low-density lipoprotein cholesterol: executive summary". Clin. Chem. 41: 1414-1420 (1995).

SYMBOLS

 Read instruction for use

 CE mark (requirement of 98/79 regulation)

 Storing temperature limits

 In vitro medical device

 Producer