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BIBLIOGRAFIA

Grove TH. Effect of reagent pH on determination of high-density lipoprotein cholesterol by precipitation with sodium phosphotungstate-magnesium. *Clin. Chem.* 1979; 25: 560-564.

Burstein M., Scholnick HR and Morfin R. rapid method for the isolation of lipoproteins from human serum by precipitation with polyanions. *Scand J Clin Lab Invest* 1980 ; 40 : 583-595.

National Cholesterol Education Program Expert Panel. Third report of the national Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (ATP III). NIH Publication. Bethesda: national Heart, Lung, and Blood Institute; 2001.

Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Press, 1995.

Tietz NW. Clinical guide to laboratory tests, 2nd ed. Saunders Co, 1991.

Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997.

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Fabbricante

Precipitating reagent for measurement of HDL cholesterol in serum or plasma - PEG 6000 precipitating method

The following procedure and performances refer to the use of precipitating joined to reagent "CHOLESTEROL" ref. 10027

PRINCIPLE

Low density lipoproteins (LDL and VLDL) and chylomicron fractions precipitate quantitatively by adding buffered PEG 6000. After centrifugation, is determined the concentration of HDL (high density lipoprotein) cholesterol fraction which remains in the supernatant.

REAGENTS

- R1 Buffered PEG 6000 20.0 g/l
R2 HDL cholesterol standard 50 mg/dl (1.29 mmol/l).

SAMPLE

- Serum or EDTA plasma.

Note

- Anticoagulants as heparin, EDTA, oxalate or fluoride do not interfere.
- The HDL cholesterol is stable in the samples for 7 days if stand at 2-8°C.

REFERENCE VALUES

	Normal	Moderate risk	High risk
Man	> 55 mg/dl	35 - 55 mg/dl	< 35 mg/dl
Woman	> 65 mg/dl	45 - 65 mg/dl	< 45 mg/dl

References values are considered indicatives since each laboratory should establish references ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

PREPARATION OF REAGENTS

Reagents are liquid and ready to use. Keep out the reagents from refrigerator only for the use and recap them immediately.

STORAGE AND STABILITY

- Store the kit at room temperature. Do not freeze the reagents.

NOTE

- Evaluate carefully the results if working reagent absorbance is > 0.600 at 510 nm.
- Avoid direct light, contamination and evaporation.
- The volumes in the procedure can be changed proportionally.
- In case of complaint or quality control request, refer to the lot number on the package or the lot number on the singles vials.

AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

PRECAUTION IN USE

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of components is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances.

However the reagent should be handled with caution, according to good laboratory practice.

Caution: the reagents contain Sodium Azide (0.095%) as preservative. Avoid swallowing and contacting with skin, eyes and mucous membranes.

WASTE MANAGEMENT

Please refer to the local legal requirements.

PROCEDURE

Pipette into centrifuge tubes:

SAMPLE	500 µl
REAGENT R1	500 µl
Mix and after 5' centrifuge at 300 rpm for 10'. Check up the supernatant is clear. Determine HDL cholesterol on supernatant with cholesterol reagent.	

Wavelength	λ: 510 nm
Working temperature	37°C
Optical path	1 cm
Reaction	"end point"

Bring the reagents at 15-25°C before using them.

- PROCEDURE

	BLANK	STD	SAMPLE
SUPERNATANT	--	--	50 µl
STANDARD	--	50 µl	--
REAGENT CHOLESTEROL	1000 µl	1000 µl	1000 µl
Mix then incubate for 5' at 37°C or 10' at 15-25°C. Measure the absorbance of the sample (EC) and standard (ES) against the reagent blank.			

CALCULATION

$EC/ES \times \text{Conc STD} \times 2 = \text{mg (mmol) of cholesterol HDL/dl (L) of sample}$

CONVERSION FACTOR

$\text{Cholesterol [mg/dl]} \times 0.02586 = \text{Cholesterol [mmol/l]}$

ANALYTICAL PERFORMANCES

Interferences

Bilirubin does not interfere up to concentration of 15 mg/dl.
Hemoglobin does not interfere up to concentration of 500 mg/dl.
Triglycerides do not interfere up to concentration of 1000 mg/dl.

Linearity

Reaction is linear up to a concentration of 150 mg/dl with a range of 5.03-150 mg/dl. Samples with values exceeding 150 mg/dl must be diluted with saline solution. Multiply, then, the result for diluting factor.

"Intra-Assay" precision (within-Run)

Determined on 30 samples for each control (N-H) (Normal-High).

Results:

MEAN (mg/dl)	N = 26.20	H = 88.80
S.D.	0.84	1.30
C.V.%	3.19	1.47

"Inter-Assay" precision (between-Run)

Determined on 15 samples for each control (N-H) for 3 days.

Results:

MEAN (mg/dl)	N = 26.07	H = 86.73
S.D.	0.95	5.63
C.V.%	3.63	6.73

Analytical sensitivity

The test sensitivity in terms of detection limit is 5.03 mg/dl.

Correlation

A study based comparing this method with a similar method on 21 samples has given a correlating factor $r = 0.95$

BIBLIOGRAPHY

Grove TH. Effect of reagent pH on determination of high-density lipoprotein cholesterol by precipitation with sodium phosphotungstate-magnesium. *Clin. Chem.* 1979; 25: 560-564.
Burstein M., Scholnick HR and Morfin R. rapid method for the isolation of lipoproteins from human serum by precipitation with polyanions. *Scand J Clin Lab Invest* 1980 ; 40 : 583-595.
National Cholesterol Education Program Expert Panel. Third report of the national Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (ATP III). NIH Publication. Bethesda: national Heart, Lung, and Blood Institute; 2001.
Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Press, 1995.
Tietz NW. Clinical guide to laboratory tests, 2nd ed. Saunders Co, 1991.
Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997.

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 CE regulation)

Storing temperature limits

In vitro medical device

Producer