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BIBLIOGRAFIA

Henderson, A.R., Donald W.M., Enzymes, Tietz Fundamentals of Clinical Chemistry, 5th Ed. Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA) (2001), 352
 Sanhai, W.R., Christenson, R.H., Cardiac and muscle disease, Clinical Chemistry; Theory, Analysis, Correlation, 4th Ed. Kaplan, L.A. Pesce, A.J. Kazmierczak, S.C.(Mosby Inc. eds St. Louis USA) (2003), 566 and Appendix.
 Schumann, G. et al. Clin. Chem. Lab. Med. (2002), 40, 635

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 **Fabbricante**

Kit for measurement of CK MB in serum UV method with immunoinhibition of M fraction - IFCC

PRINCIPLE

An antibody is incorporated in the CK reagent. This antibody will bind to and inhibit the activity of the M subunit of CK-MB. This means that only the activity of the B subunit in serum is measured. If this activity is multiplied by a factor of 2, it will give the activity of CK-MB in serum.

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- R1** Imidazole buffer pH 7.20 100 mmol/l; magnesium acetate 10 mmol/l; antibody CK-M;
- R2** N- Acetyl-L-Cysteine 20.0 mmol/l; ADP 2.0 mmol/l; AMP 5.0 mmol/l; NADP 2.0 mmol/l; D-Glucose 20.0 mmol/l; diadenosine 5 phosphate 10.0 µmol/l; EDTA 2.0 mmol/l; creatine phosphate 29.0 mmol/l; HK ≤ 3500 U/l; G6P-DH 2000 U/l;
- R3** Control serum. Add water as reported in the label.

SAMPLE

- Serum not hemolysed.

NOTE

- Sample must be analyzed immediately or stored, protected from air and light, for 2 days at 2-8°C or 1 months at -20°C.
- Do not use hemolysed samples because hemolysis may increase CK-MB concentration due to adenylate kinase release.
- The method also measure CK-BB isoenzyme present in serum alone or complexed with immunoglobulins (macro-CK). The activity of the isoenzymes is negligible. However, if a significant amount of CK-BB activity is present the CK-MB activity will be overestimated.

REFERENCE VALUES

Serum (U/l)	37°C
CK man	< 190 U/l
CK woman	< 170 U/l
CK MB	< 25 U/l

Reference values are considered indicative since each laboratory should establish reference ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

PREPARATION OF REAGENTS

Dissolve the content of a **R2** vial with 3 ml of **R1** buffer.
 Wait about 15 minutes before use.
 Keep out the reagents from refrigerator only for the use and recap them immediately.

STORAGE AND STABILITY

- Store the kit at 2-8°C.
- Reagents, if stored at 2-8°C and protected from light are stable up to expiry date on the label.
- Reconstituted reagent stability: 3 days at 2-8°C – 24 hours at 15-25°C.

NOTE

- The kit, according to this method, must be used in manual procedures. About automatic using follow specific applications.
- Avoid direct light, contamination and evaporation.
- The volumes in the procedure can be changed proportionally.
- In case of complaint or quality control request, refer to the lot number on the package or the lot number on the singles vials.

AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

PRECAUTION IN USE

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of components is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances.

However the reagent should be handled with caution, according to good laboratory practice.

Caution: reagent R1 contains Sodium Azide (0.095%) as preservative. Sodium azide can react with copper and lead plumbing to form explosive compounds. Avoid swallowing and contacting with skin, eyes and mucous membranes.

WASTE MANAGEMENT

Please refer to the local legal requirements.

PROCEDURE

Wavelength	λ: 340 nm
Working temperature	37°C
Optical path	1 cm
Reaction	"kinetic" (increasing)

Bring the reagents at 15-25°C before using them.

WORKING REAGENT	STD	SAMPLE
200 µl	200 µl	200 µl
SAMPLE	--	10 µl

Mix, then incubate for 5'. Read absorbance values every minute during 3 minutes. Calculate ΔE/min.

CALCULATION

$$\begin{aligned} \text{CK-B (U/l)} &= \Delta/\text{min} \times 3.376 \\ \text{CK-MB} &= \text{CK-B} \times 2 \text{ or } \text{CK-MB} = \Delta/\text{min} \times 6752 \\ \% \text{ CK-MB} &= \frac{\text{CK-MB}}{\text{Total - CK}} \times 100 \end{aligned}$$

ANALYTICAL PERFORMANCES

The following data were obtained using COBAS MIRA analyser (37°C).

Interferences

Unconjugated bilirubin does not interfere up to concentration of 206 µmol/l.
 Conjugated bilirubin does not interfere up to concentration of 102 µmol/l.
 Turbidity do not interfere up to concentration of 6 g/l.

Linearity

Reaction is linear up to a concentration of 600 U/l. Samples with values exceeding 600 U/l must be diluted with saline solution. Multiply, then, the result for diluting factor.

"Intra-Assay" precision (within-Run)

Determined on 20 samples for each control (L-N-H) (Low-Normal-High).

MEAN	L = 27 U/l	N = 77 U/l	H = 202 U/l
C.V.%	3.1 %	2.2 %	1.2 %

"Inter-Assay" precision (between-run)

Determined on 80 samples for each control (L-N-H) for 3 days.

MEAN	L = 24 U/l	N = 72 U/l	H = 199 U/l
C.V.%	6.9 %	2.8 %	2.1 %

Analytical sensitivity

The test sensitivity in terms of detection limit is 4 U/l.

Correlation

A study based comparing this method with a similar method on 21 samples has given a correlating factor **r = 0.99**

BIBLIOGRAPHY

Henderson, A.R., Donald W.M., Enzymes, Tietz Fundamentals of Clinical Chemistry, 5th Ed. Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA) (2001), 352
 Sanhai, W.R., Christenson, R.H., Cardiac and muscle disease, Clinical Chemistry: Theory, Analysis, Correlation, 4th Ed. Kaplan, L.A. Pesce, A.J. Kazmierczak, S.C. (Mosby Inc. eds St. Louis USA) (2003), 566 and Appendix.
 Schumann, G. et al. Clin. Chem. Lab. Med. (2002), 40, 635

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer