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Fabbricante

**Kit for measurement of creatin kinase in serum or plasma
 UV kinetic optimized method DGKC-IFCC***

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PRINCIPLE

Creatin kinase (CK) catalyzes reaction between creatine phosphate and ADP, giving creatine and ATP. This one, in presence of glucose and hexokinase (HK) is transformed into ADP and glucose-6-phosphate which, with glucose-6-phosphate-dehydrogenase (G6P-DH) intervention, forms glucose-6-phosphogluconate while NADP⁺ is reduced to NADPH. The CK activity in the sample is calculated measuring the absorbance variation, due to the transformation of NADP⁺ into NADPH.

REAGENTS

R1 Imidazole buffer pH 7.10 100.0 mmol/l; magnesium acetate 10.0 mmol/l
R2 ADP 2.0 mmol/l; G6P-DH ≥ 2000 U/l; creatine phosphate 30.0 mmol/l;

EDTA 2.0 mmol/l; excipients

SAMPLE

- Serum or heparinized plasma.

Note

- Store the samples protected from light.
 - The CK activity is stable for 7 days at 2-8°C or 30 days at -20°C.

REFERENCE VALUES

Serum - plasma [U/l]	25°C	30°C	37°C
Women	10-70	15-110	24-170
Men	10-80	15-130	24-195

References values are considered indicatives since each laboratory should establish references ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

PREPARATION OF REAGENTS

Dissolve the content of **R2** vial with **R1** solution. In case of unitary test dissolve the content of **R2** vial with 3 ml of **R1** solution.

Keep out the reagents from refrigerator only for the use and recap them immediately.

STORAGE AND STABILITY

- Store the kit at 2-8°C. Do not freeze the reagents.
 - Reconstituted reagent stability: 7 days at 2-8°C.

NOTE

- The kit, according to this method, must be used in manual procedures. About automatic using follow specific applications.
 - Evaluate carefully the results if working reagent absorbance is > 0.200 at 340 nm.
 - Avoid direct light, contamination and evaporation.
 - The volumes in the procedure can be changed proportionally.
 - In case of complaint or quality control request, refer to the lot number on the package or the lot number on the singles vials.

AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

Suggested serum:

REF 20350 Precise Norm

REF 20360 Precise Path

PRECAUTION IN USE

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of components is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances.

However the reagent should be handled with caution, according to good laboratory practice.

Caution: the reagents contain Sodium Azide (0.095%) as preservative. Avoid swallowing and contacting with skin, eyes and mucous membranes.

WASTE MANAGEMENT

Please refer to the local legal requirements.

PROCEDURE

Wavelength λ: 340 (334 o 365) nm
 Working temperature 37°C (25° o 30°C)
 Optical path 1 cm
 Reaction kinetic (increasing)

Bring the reagents at 15-25°C before using them.

	BLANK	SAMPLE
WORKING REAGENT	1500 µl	1500 µl
DISTILLED WATER	50 µl	--
SAMPLE	--	50 µl

Mix and after 2' measure the absorbance. Make at least two readings at a distance of 60". Calculate the absorbance variation ΔE/min.

CALCULATION

340 nm	CK [U/l] = ΔE/min x 4921
365 nm	CK [U/l] = ΔE/min x 9118

The factor and the reagent performances are related to 37°C and 340 nm.

ANALYTICAL PERFORMANCES
Interferences

Bilirubin does not interfere up to concentration of 50 µmol/l.
 Triglycerides do not interfere up to concentration of 1000 mg/dl.
 Hemoglobin does not interfere up to concentration of 1.0 g/l.

Linearity

Reaction is linear up to a concentration of 1200 U/l. Samples with values exceeding 1200 U/l must be diluted with saline solution. Multiply, then, the result for diluting factor.

"Intra-Assay" precision (within-Run)

Determined on 30 samples for each control (N-H) (Normal-High).

Results:

MEAN [U/l]	N = 142.40	H = 451.90
S.D.	2.17	2.23
C.V.%	1.52	0.49

"Inter-Assay" precision (between-run)

Determined on 15 samples for each control (N-H) for 3 days.

Results:

MEAN [U/l]	N = 142.13	H = 452.93
S.D.	2.67	3.17
C.V.%	1.88	0.70

Analytical sensitivity

The test sensitivity in terms of detection limit is 9.5 U/l.

Correlation

A study based comparing this method with a similar method on 21 samples has given a correlating factor **r = 0.99**

BIBLIOGRAPHY

Ann. Biol. Clin., (1982), 40, 99.

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer