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BIBLIOGRAFIA

Stein, W. "Creatine kinase (total activity), creatine kinase isoenzymes and variants". In: Thomas L, ed. Clinical laboratory diagnostics. Frankfurt: TH-Books Verlagsgesellschaft; p. 71-80, (1998).
 Recommendations of the German Society for Clinical Chemistry. "Standardization of methods for the estimation of enzyme activities in biological fluids: Standard method for the determination of creatine kinase activity". J. Clin. Chem. Clin. Biochem., 15, 255-60, (1977).
 Moren, L.G., Clin. Chem., 23, 1569 (1977).
 Kaplan, L.A., Pesce, A..J.: "Clinical Chemistry", Mosby Ed. (1996).

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Fabbricante

Kit for measurement of creatin kinase in serum or plasma Kinetic optimized method DGKC-IFCC*

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PRINCIPLE

Creatin kinase (CK) catalyzes reaction between creatine phosphate and ADP, giving creatine and ATP. This one, in presence of glucose and hexokinase (HK) is transformed into ADP and glucose-6-phosphate which, with glucose-6-phosphate-dehydrogenase (G6P-DH) intervention, forms glucose-6-phosphogluconate while NADP⁺ is reduced to NADPH. The CK activity in the sample is calculated measuring the absorbance variation, due to the reduction of NADP⁺ into NADPH.

REAGENTS

R1 Imidazole buffer pH 6.7 125.0 mmol/l; hexokinase (HK) ≥ 6800 U/l; magnesium acetate 12.5 mmol/l; D-glucose 25.0 mmol/l; EDTA 2.02 mmol/l; NADP 2.52 mmol/l; diadenosine phosphate 10.0 mmol/l; N-acetyl-cysteine (NAC) 25.0 mmol/l
R2 G6P-DH ≥ 8800 U/l; creatine phosphate 250 mmol/l; ADP 15.2 mmol/l; AMP 25.0 mmol/l; di-Adenosine-5-pentaphosphate 103 mmol/l.

CAMPIONI

- Serum free of hemolysis or heparin plasma.

Note

- Store the samples protected from light at 2-8°C for 7 days (maximum).
 - The CK activity is stable for 7 days at 2-8°C or 30 days at -20°C.

REFERENCE VALUES

Serum - plasma [U/l]	25°C	30°C	37°C
Adults: Women	10 - 70	115 - 110	24 - 170
Men	10 - 80	15 - 125	24 - 195
Children: Babies	192 - 494	300 - 770	468 - 1200
≤ 5 days	80 - 288	125 - 450	195 - 700
< 6 months	17 - 136	27 - 212	41 - 330
> 6 months	10 - 94	15 - 147	24 - 229

References values are considered indicatives since each laboratory should establish references ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

STORAGE AND STABILITY

- Store the kit at 2-8°C. Do not freeze the reagents.
 - After opening, the vials **R1** and **R2** are stable up to the expiry date if recapped immediately and protected from contamination, evaporation, direct light, and stored at the correct temperature.
 - Working solution stability (**R1+ R2**): 21 days at 2-8°C or 5 days at room temperature (15-25°C).

PREPARATION OF REAGENTS

Reagents are liquid and ready to use. About using as monoreagent (sample-starter procedure) add to every 4 ml of **R1** reagent, 1 ml of **R2** reagent. In case of 3 ml package use 2.5 ml of **R1** and 0.5 ml of **R2**.
 Keep out the reagents from refrigerator only for the use and recap them immediately.

NOTE

- The kit, according to this method, must be used in manual procedures. About automatic using follow specific applications.
 - Evaluate carefully the results if working reagent absorbance is ≥ 1.60 at 340 nm or in presence of particles and turbidity.
 - Avoid direct light, contamination and evaporation.
 - The volumes in the procedure can be changed proportionally.
 - In case of complaint or quality control request, refer to the lot number on the package or the lot number on the singles vials.

AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

Suggested serum:

REF 20350 Precise Norm **REF** 20360 Precise Path

PRECAUTION IN USE

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of components is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances.

However the reagent should be handled with caution, according to good laboratory practice.

Caution: the reagents contain Sodium Azide (0.095%) as preservative. Avoid swallowing and contacting with skin, eyes and mucous membranes.

WASTE MANAGEMENT

Please refer to the local legal requirements.

PROCEDURE

Wavelength λ: 340 (334 o 365) nm
 Working temperature 37°C (25° o 30°C)
 Optical path 1 cm
 Reaction kinetic (increasing)

Bring the reagents at 15-25°C before using them.

- MONOREAGENT PROCEDURE "sample starter"

WORKING REAGENT	BLANK	SAMPLE
DISTILLED WATER	1000 µl	1000 µl
SAMPLE	20 µl	--
	--	20 µl

Mix and after 2 min. measure the absorbance. Make at least three readings at a distance of 60 sec. Calculate the average absorbance difference per minute (ΔE/min).

- BIREAGENT PROCEDURE "substrate starter"

REAGENT R1	BLANK	SAMPLE
DISTILLED WATER	800 µl	800 µl
SAMPLE	20 µl	--
	--	20 µl
REAGENT R2	200 µl	200 µl

Mix, incubate at 37°C for about 2 min. and then add:
 Mix and after 1 min. measure the absorbance. Make at least three readings at a distance of 60 sec. Calculate the average absorbance difference per minute (ΔE/min).

CALCULATION

$$CK [U/l] = \Delta E/min \times 8095$$

The factor and the reagent performances are related to 37°C and 340 nm.

ANALYTICAL PERFORMANCES

Interferences

Triglycerides do not interfere up to concentration of 7 mmol/l.
 Hemoglobin does not interfere up to concentration of 50 mg/dl.
 Glucose does not interfere up to concentration of 70 mg/dl.

Linearity

Reaction is linear from a concentration of 2.2U/l to a concentration of 2000 U/l. Samples with values exceeding 2000 U/l must be diluted with saline solution. Multiply, then, the result for diluting factor.

Sensitivity: 1U/l = 0.0012 ΔA/min.

"Intra-Assay" precision (within-Run)

Determined on several controls (N-H, Normal-High).

Results:

MEDIA [U/l]	N = 147	H = 494
D.S.	1.23	3.60
C.V.%	0.84	0.73

"Inter-Assay" precision (between-run)

Determined on several controls (N-H, Normal-High).

Results:

MEDIA [U/l]	N = 145	H = 485
D.S.	2.91	8.97
C.V.%	2.01	1.85

Correlation

A study based comparing this method with a similar method has given a correlating factor **r = 0.9995**.

BIBLIOGRAPHY

Stein, W. "Creatine kinase (total activity), creatine kinase isoenzymes and variants". In: Thomas L, ed. Clinical laboratory diagnostics. Frankfurt: TH-Books Verlagsgesellschaft; p. 71-80, (1998).

Recommendations of the German Society for Clinical Chemistry. "Standardization of methods for the estimation of enzyme activities in biological fluids: Standard method for the determination of creatine kinase activity". J. Clin. Chem. Clin. Biochem., 15, 255-60, (1977).

Moren, L.G., Clin. Chem., 23, 1569 (1977).

Kaplan, L.A., Pesce, A.: "Clinical Chemistry", Mosby Ed. (1996).

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer