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Fabbricante

## Kit for measurement of creatin kinase in serum or plasma Kinetic optimized method DGKC-IFCC\*

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### PRINCIPLE

Creatin kinase (CK) catalyzes reaction between creatine phosphate and ADP, giving creatine and ATP. This one, in presence of glucose and hexokinase (HK) is transformed into ADP and glucose-6-phosphate which, with glucose-6-phosphate-dehydrogenase (G6P-DH) intervention, forms glucose-6-phosphogluconate while NADP<sup>+</sup> is reduced to NADPH. The CK activity in the sample is calculated measuring the absorbance variation, due to the reduction of NADP<sup>+</sup> into NADPH.

### REAGENTS

- R1** Imidazole buffer pH 6.3 100.0 mmol/l; hexokinase (HK) 2500 U/l; magnesium acetate 10.0 mmol/l; D-glucose 20.0 mmol/l; EDTA 2.0 mmol/l; AMP 20.0 mmol/l; NADP 2.0 mmol/l; diadenosine phosphate 10.0 mmol/l; N-acetyl-cysteine (NAC) 20.0 mmol/l
- R2** Imidazole buffer pH 9.13 100.0 mmol/l; G6P-DH 8000 U/l; creatine phosphate 30.0 mmol/l; ADP 10.50 mmol/l; EDTA Na2 2.0 mmol/l

### CAMPIONI

- Serum or heparinized plasma.

### Note

- Store the samples protected from light.
- The CK activity is stable for 7 days at 2-8°C or 30 days at -20°C.

### REFERENCE VALUES

Serum - plasma [U/l]	25°C	30°C	37°C
Adults: Women	10 - 70	115 - 110	24 - 170
Men	10 - 80	15 - 125	24 - 195
Children: Babies	192 - 494	300 - 770	468 - 1200
≤ 5 days	80 - 288	125 - 450	195 - 700
< 6 months	17 - 136	27 - 212	41 - 330
> 6 months	10 - 94	15 - 147	24 - 229

References values are considered indicative since each laboratory should establish reference ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

### STORAGE AND STABILITY

- Store the kit at 2-8°C. Do not freeze the reagents.
- After opening, the vials **R1** and **R2** are stable up to the expiry date if recapped immediately and protected from contamination, evaporation, direct light, and stored at the correct temperature.
- Working solution stability (**R1+ R2**): 21 days at 2-8°C.

### PREPARATION OF REAGENTS

Reagents are liquid and ready to use. About using as monoreagent (sample-starter procedure) add to every 4 ml of **R1** reagent, 1 ml of **R2** reagent. In case of 3 ml package use 2.5 ml of **R1** and 0.5 ml of **R2**. Keep out the reagents from refrigerator only for the use and recap them immediately.

### NOTE

- The kit, according to this method, must be used in manual procedures. About automatic using follow specific applications.
- Evaluate carefully the results if working reagent absorbance is > 0.200 at 340 nm.
- Avoid direct light, contamination and evaporation.
- The volumes in the procedure can be changed proportionally.
- In case of complaint or quality control request, refer to the lot number on the package or the lot number on the singles vials.

### AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

### QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

### Suggested serum:

REF 20350 Precise Norm                      REF 20360 Precise Path

### PRECAUTION IN USE

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of components is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances. However the reagent should be handled with caution, according to good laboratory practice.

**Caution:** the reagents contain Sodium Azide (0.095%) as preservative. Avoid swallowing and contacting with skin, eyes and mucous membranes.

### WASTE MANAGEMENT

Please refer to the local legal requirements.

### PROCEDURE

Wavelength	λ: 340 (334 o 365) nm
Working temperature	37°C (25° o 30°C)
Optical path	1 cm
Reaction	kinetic (increasing)

Bring the reagents at 15-25°C before using them.

### - MONOREAGENT PROCEDURE "sample starter"

WORKING REAGENT	BLANK	SAMPLE
DISTILLED WATER	1000 µl	1000 µl
SAMPLE	25 µl	--
	--	25 µl

Mix and after 1' measure the absorbance. Make at least two readings at a distance of 60". Calculate the average absorbance difference per minute (ΔE/min).

### - BIREAGENT PROCEDURE "substrate starter"

REAGENT R1	BLANK	SAMPLE
DISTILLED WATER	800 µl	800 µl
SAMPLE	25 µl	--
	--	25 µl

Mix, incubate at 37°C for about 5' and then add:

REAGENT R2	BLANK	SAMPLE
	200 µl	200 µl

Mix and after 1' measure the absorbance. Make at least two readings at a distance of 60". Calculate the average absorbance difference per minute (ΔE/min).

### CALCULATION

$$CK [U/l] = \Delta E/min \times 6507$$

The factor and the reagent performances are related to 37°C and 340 nm.

### ANALYTICAL PERFORMANCES

#### Interferences

Bilirubin does not interfere up to concentration of 40 mg/dl.  
Triglycerides do not interfere up to concentration of 2000 mg/dl.  
Hemoglobin does not interfere up to concentration of 20 mg/dl.  
Ascorbate acid does not interfere up to concentration of 30 mg/dl.

#### Linearity

Reaction is linear up to a concentration of 1600 U/l (corresponding to ΔE/min of 0.246 a 340 nm). Samples with values exceeding 1600 U/l must be diluted with saline solution. Multiply, then, the result for diluting factor.

#### "Intra-Assay" precision (within-Run)

Determined on 30 samples for each control (L-N-H) (Low-Normal-High).

Results:

MEAN [U/l]	L = 82.7	N = 168.9	H = 413.9
S.D.	2.00	4.05	15.92
C.V.%	2.42	2.40	3.85

#### "Inter-Assay" precision (between-run)

Determined on 15 samples for each control (L-N-H) for 3 days.

Results:

MEAN [U/l]	L = 87.80	N = 166.82	H = 398.93
S.D.	2.87	4.26	9.33
C.V.%	3.27	2.55	2.34

#### Correlation

A study based comparing this method with a similar method on 21 samples has given a correlating factor **r = 0.97**

#### BIBLIOGRAPHY

Stein, W. "Creatine kinase (total activity), creatine kinase isoenzymes and variants". In: Thomas L, ed. Clinical laboratory diagnostics. Frankfurt: TH-Books Verlagsgesellschaft; p. 71-80, (1998).

Recommendations of the German Society for Clinical Chemistry. "Standardization of methods for the estimation of enzyme activities in biological fluids: Standard method for the determination of creatine kinase activity". J. Clin. Chem. Clin. Biochem., 15, 255-60, (1977).

Moren, L.G., Clin. Chem., 23, 1569 (1977).

Kaplan, L.A., Pesce, A..J.: "Clinical Chemistry", Mosby Ed. (1996).

### SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer