 Per uso diagnostico in vitro
10228 - 5x50 ml
V10228 - 5x20+5x5 ml
H10228 - 4x40+2x20 ml

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BIBLIOGRAFIA

Committee of the Scand. Soc. for Clin. Chem. - Scand. J. Clin. Lab. Invest., 36, 119 (1976).

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Fabbricante

 For in vitro medical device
10228 - 5x50 ml
V10228 - 5x20+5x5 ml
H10228 - 4x40+2x20 ml

Kit for measurement of gamma-GT in serum or plasma SZASZ method

PRINCIPLE

Gamma-glutamyl group is transferred from gamma-glutamyl-carboxy-nitroanilide to glycylglycine by gamma-glutamyltransferase (GGT). The reaction releases p-nitroaniline whose absorbance is directly proportional to the GGT activity in the sample.

REAGENTS

R1 Goods buffer pH 8.25 300.0 mmol/l; glycylglycine 120.0 mmol/l
R2 L-gamma-glutamyl-3-carboxy-4-nitroanilide 16.0 mmol/l

SAMPLE

- Serum or EDTA plasma.

Note

- Do not use samples with haemolysis.
 - The GGT activity is stable in the samples 7 days at 2-8°C.

REFERENCE VALUES

Serum - plasma [U/l]	30°C	37°C
Women	6-24	8-31
Men	7-38	10-50

References values are considered indicatives since each laboratory should establish references ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

STORAGE AND STABILITY

- Store the kit at 2-8°C. Do not freeze the reagents.
 - After opening, the vials **R1** and **R2** are stable up to the expiry date if recapped immediately and protected from contamination, evaporation, direct light, and stored at the correct temperature.
 - Working solution stability (**R1+ R2**): 21 days at 2-8°C.

PREPARATION OF REAGENTS

Reagents are liquid and ready to use. About using as monoreagent ("sample-starter" procedure) add to every 4 ml of **R1** reagent, 1 ml of **R2** reagent. Keep out the reagents from refrigerator only for the use and recap them immediately.

NOTE

- The kit, according to this method, must be used in manual procedures. About automatic using follow specific applications.
 - Evaluate carefully the results if working reagent absorbance is > 0.800 at 405 nm.
 - Avoid direct light, contamination and evaporation.
 - The volumes in the procedure can be changed proportionally.
 - In case of complaint or quality control request, refer to the lot number on the package or the lot number on the singles vials.

AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

Suggested serum:

REF 20350 Precise Norm REF 20360 Precise Path

PRECAUTION IN USE

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of inactive components (preservatives, surfactants and so on) is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances.

However the reagent should be handled with caution, according to good laboratory practice.

Caution: the reagents contain Sodium Azide (0.095%) as preservative. Avoid swallowing and contacting with skin, eyes and mucous membranes.

WASTE MANAGEMENT

Please refer to the local legal requirements.

PROCEDURE

Wavelength λ : 405 (400-420) nm
 Working temperature 37°C
 Optical path 1 cm
 Reaction Kinetic (increasing)

Bring the reagents at 15-25°C before using them.

- MONOREAGENT PROCEDURE "sample starter"

	BLANK	SAMPLE
WORKING REAGENT	1000 μ l	1000 μ l
DISTILLED WATER	100 μ l	--
SAMPLE	--	100 μ l

Mix, then incubate for 1' a 37°C. Measure the absorbance of sample (EC) against distilled water decrease per minute during 3 minutes. Calculate the average absorbance difference per minute ($\Delta E/min$).

- BIREAGENT PROCEDURE "substrate starter"

	BLANK	SAMPLE
REAGENT R1	800 μ l	800 μ l
DISTILLED WATER	100 μ l	--
SAMPLE	--	100 μ l

Mix, incubate at 37°C for 1' and then add:

	BLANK	SAMPLE
REAGENT R2	200 μ l	200 μ l

Mix, then incubate for 1' a 37°C. Measure the absorbance of sample (EC) against distilled water decrease per minute during 3 minutes. Calculate the average absorbance difference per minute ($\Delta E/min$).

CALCULATION

$$GGT [U/l] = \Delta E/min \times 1159$$

The factor and the reagent performances are related to 37°C and 405 nm.

ANALYTICAL PERFORMANCES

Interferences

Bilirubin does not interfere up to concentration of 40 mg/dl.
 Triglycerides do not interfere up to concentration of 2000 mg/dl.
 Hemoglobin does not interfere up to concentration of 400 mg/dl.

Linearity

Reaction is linear up to a concentration of 1200 U/l. Samples with values exceeding 1200 U/l must be diluted with saline solution. Multiply, then, the result for diluting factor.

"Intra-Assay" precision (within-Run)

Determined on 30 samples for each control (N-H) (Normal-High).

Results:

MEAN (U/l)	N = 28.50	H = 117.40
S.D.	0.78	2.78
V.C. %	2.73	2.36

"Inter-Assay" precision (between-run)

Determined on 15 samples for each control (N-H) for 3 days.

Results:

MEAN (U/l)	N = 29.73	H = 112.91
S.D.	1.08	2.26
V.C. %	3.65	2.00

Analytical sensitivity

The test sensitivity in terms of detection limit is 6.2 U/l.

Correlation

A study based comparing this method with a similar method on 21 samples has given a correlating factor **r = 1.00**

BIBLIOGRAPHY

Committee of the Scand. Soc. for Clin. Chem. - Scand. J. Clin. Lab. Invest., 36, 119 (1976).
 Szasz G: Clin. Chem., 22, 2051 (1976).
 Persijn JP, van der Silk W: J. Clin. Chem. Clin. Biochem., 14, 421 (1976).
 Kaplan LA, Pesce AJ: "Clinical Chemistry", Mosby Ed. (1996).

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer