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
Fabbricante


Kit for measurement of phosphorous in serum, plasma and urine Colorimetric Molybdenum Blue method

PRINCIPLE

The inorganic phosphorus reacts with ammonium molybdate in an acid medium to form a phosphomolybdate complex which absorbs light at 340 nm.

The absorbance at this wavelength is directly proportional to the amount of inorganic phosphorus present in the sample.

REAGENTS

R1 Hydrochloric acid 1.1 mmol/l; ammonium molybdate 0.68 mmol/l;
surface-active agents polyanionic
Standard phosphorous 5 mg/dl (1.61 mmol/l)

SAMPLE

- Serum or plasma. Diluted urine 1:10.

Note

- Do not use samples with haemolysis. Do not use sodium fluoride as anticoagulant.
- The phosphorous is stable in the samples up to 7 days at 2-8°C.

REFERENCE VALUES

| | | |
|----------------|----------|-----------------|
| Serum - plasma | Men | 2.1 - 5.6 mg/dl |
| | Women | 1.5 - 6.8 mg/dl |
| | Children | 4.0 - 7.0 mg/dl |
| Urine | | 0.3 - 1.0 g/24h |

References values are considered indicatives since each laboratory should establish references ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

STORAGE AND STABILITY

- Store the kit at 15-25°C.
- After opening, the vials **R1** and **R2** are stable up to the expiry date if recapped immediately and protected from contamination, evaporation, direct light, and stored at the correct temperature.

PREPARATION OF REAGENTS

Reagents are liquid and ready to use.

NOTE

- The kit, according to this method, must be used in manual procedures. About automatic using follow specific applications.
- Evaluate carefully the results if reagent R1 absorbance is > 0.600 at 340 nm.
- Avoid direct light, contamination and evaporation.
- The volumes in the procedure can be changed proportionally.
- In case of complaint or quality control request, refer to the lot number on the package or the lot number on the singles vials.

AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

Suggested serum:

REF 20350 Precise Norm REF 20360 Precise Path

PRECAUTION IN USE

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of inactive components (preservatives, surfactants and so on) is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances.

However the reagent should be handled with caution, according to good laboratory practice.

Caution: the reagents contain Sodium Azide (0.095%) as preservative. Avoid swallowing and contacting with skin, eyes and mucous membranes.

WASTE MANAGEMENT

Please refer to the local legal requirements.

PROCEDURE

Wavelength λ : 340 (334 o 365) nm
Working temperature 37°C (15 -25°C)
Optical path 1 cm
Reaction "end point" (increasing)

- MONOREAGENT

| | BLANK | STD | SAMPLE |
|-----------------|--------------|--------------|--------------|
| REAGENT R1 | 1000 μ l | 1000 μ l | 1000 μ l |
| DISTILLED WATER | 10 μ l | -- | -- |
| SAMPLE | -- | -- | 10 μ l |
| STANDARD | -- | 10 μ l | -- |

Mix, then incubate for 5' at 37°C (15-25°C). Measure the absorbance of sample (EC) and standard (ESTD) against the reagent blank.

CALCULATION

$$\text{Phosphorous [mg/dl] o [mmol/l]} = \text{EC / ESTD} \times \text{Conc. STD}$$

Diluted urines: multiply the result for diluting factor.

CONVERSION FACTOR

$$\text{Phosphorous [mg/dl]} \times 0.32 = \text{Phosphorous [mmol/l]}$$

ANALYTICAL PERFORMANCES

The reagent performances are related to 340 nm and 37°C.

Interferences

Bilirubin does not interfere up to concentration of 20 mg/dl.
Triglycerides do not interfere up to concentration of 1000 mg/dl.
Hemoglobin does not interfere up to concentration of 50 mg/dl.

Linearity

Reaction is linear between 0.15-20 mg/dl (0.048-6.46 mmol/l). Samples with values exceeding 300 mg/dl must be diluted with saline solution. Multiply, then, the result for diluting factor.

"Intra-Assay" precision (within-Run)

Determined on 30 samples for each control (L-N-H) (Low-Normal-High).

Results:

| | | | | |
|-------|---------|----------|----------|----------|
| MEAN | [mg/dl] | L = 2.00 | N = 4.20 | H = 6.70 |
| S.D. | | 0.05 | 0.21 | 0.14 |
| C.V.% | | 2.36 | 4.99 | 2.10 |

"Inter Assay" precision (between-Run)

Determined on 15 samples for each control (L-N-H) for 3 days.

Results:

| | | | | |
|-------|---------|----------|----------|----------|
| MEAN | [mg/dl] | L = 2.00 | N = 4.01 | H = 6.37 |
| S.D. | | 0.04 | 0.11 | 0.14 |
| C.V.% | | 2.05 | 2.68 | 2.24 |

Analytical sensitivity

The test sensitivity in terms of detection limit is 0.15 mg/dl (0.048 mmol/l).

Correlation

A study based comparing this method with a similar method on 21 samples has given a correlating factor $r = 0.97$

The urine performances are available on request.

BIBLIOGRAPHY

Talke H, Schubert GE: Klin Wchens., (1965), 43, 174.
Kaplan LA, Pesce AJ: "Clinical Chemistry", Mosby Ed. (1996).

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer