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## Ione exchange resin for glycohemoglobin detection in blood

### PRINCIPLE

A small quantity of whole blood is mixed with a lysing reagent to eliminate, with haemolysis red cells, all the hemoglobin present in them. A part of this hemolyzed preparation is in direct contact to a weak binding cation exchange resin which bind to the hemoglobin not glycosylated making free from the supernatant liquid the glycated part. This one is separated by filtration. The percentage of glycohemoglobin is determined by spectrophotometer.

### REAGENTS

- R1** Lysing reagent (KCN 10 nM + surfactant)  
**R2** Ione exchange resin (cation resin 8 mg/ml) buffered at pH 6.9  
**R3** Glycohemoglobin calibrator  
 Separator filters

### Note

**R3 reagent (calibrator)** is obtained using only blood of donors tested by an FDA method and found non-reactive for HbsAg and negative for antibodies to HIV-1/2 and HCV. However, no known test method can offer complete assurance that products derived from human blood will not transmit infectious diseases, this product should be handled as potentially infectious biological material.

- The ione exchange resin may cause irritation. Avoid contact with eyes, skin and clothing. In case of contact wash with large amounts of water.
- The lysing reagent contains cyanide. Do not mix acids. Wash hands after handling, discard by flushing with large amounts of water.
- The supernatant liquid above the resin and the lysing reagent should be clear and colourless. Turbidity or discoloration would indicate deterioration and the reagent should not be used.

### SAMPLE

- Whole blood collected with EDTA.

### Note

- Handle with same precautions used for all human blood samples.
- No special additives or preservatives other than anticoagulants are required.
- Gross lipemia may cause falsely high results. For grossly lipemic samples centrifugate the red cells and remove the lipemic plasma replacing it with an approximately equal amount of saline and proceed with performances of test. Glycosylated HbS and HbC bind lightly to the resin which produces falsely low results. Fetal hemoglobin (HbF) does not interfere significantly in the assay. The unstable fraction (aldimine) is eliminated during resin mixing and does not contribute to the glycohemoglobin value.
- Glycohemoglobin in the sample is stable for 7 days at 2-8°C.

### STORAGE AND STABILITY

- Store the kit at 2-8°C.
- Reagent **R1** and **R2** are stable until expiration date showed on the package.
- Reagent **R3**, reconstituted with 1 ml of H<sub>2</sub>O, is stable 14 days at 2-8°C sealed tightly or frozen for 8 weeks. Removed from refrigerator must be used immediately.

### AUXILIARY EQUIPMENT

Materials not included in the kit: spectrophotometer or colorimeter calibrated at 415 nm, micropipettes (from 0.02 ml to 0.10 ml), pipettes (from 0.50 ml, 3.0 ml and 5 ml), rocker or rotator and test tubes.

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The reliability of test results should be monitored routinely using suitable quality control materials analyzed in the same manner of samples.

### WASTE MANAGEMENT

Please refer to the local legal requirements.

## PROCEDURE

### Note

Resin and lysing reagent must be at room temperature before starting the test procedure.

### Hemolysate preparation

	CALIBRATOR	CONTROL	SAMPLE
LYSING REAGENT	500 µl	500 µl	500 µl
CALIBRATOR	100 µl	--	--
CONTROL	--	100 µl	--
SAMPLE	--	--	100 µl

Mix carefully until complete lysis is evident, then allow to stand for 5 min.

### Glycohemoglobin seaparation

### Note

Before use, mix the resin well by inverting the vial at least six times. Swirl the vial between each dispensing.

- Place into the resin tube 100 µl of the hemolysate relates to calibrator, control and samples.
- Put the filter separator in the tube so that the rubber sleeve is approximately 2 cm above the liquid level.
- Mix for 5'.
- Push the filter separator into the tube until the resin is firmly packed.
- Adjust the instrument to zero using deionized water, than read at 415 nm absorbance values for calibrator, control, samples. **These readings are for glycohemoglobin.**

### Total Hemoglobin fraction

The hemolysates obtained from calibrator, control and samples have to be analyzed as follows:

	CALIBRATOR	CONTROL	SAMPLE
H <sub>2</sub> O	5000 µl	5000 µl	5000 µl
CALIBRATOR	20 µl	--	--
CONTROL	--	20 µl	--
SAMPLE	--	--	20µl

Adjust the instrument to zero using deionized water, than read at 415 nm absorbance values for calibrator, control, samples.

### Note

The test samples should be read within 1 hour.

### CALIBRATION

Glycohemoglobin calibrator has an assigned value established by an accepted standard reference method.

### CALCULATION

For each sample and for the calibrator, calculate the ratio (**R**) of the glycohemoglobin absorbance to the total hemoglobin absorbance

$$R (\text{sample.}) / R(\text{cal}) \times \text{conc. cal} = \text{Glycohemoglobin HbA1 (\%)}$$

The ratio "R" (calibrator) is valid for all test performed with the kit.

### LIMITATIONS

- Blood samples with total hemoglobin greater than 18 g/dl should be diluted 1:2 with saline before assay.
- Samples from patients with hemoglobinopathies of decreased red cell survival times may show incorrect results (see "SAMPLES" section).

## EXPECTED VALUES

### Total hemoglobin (Hb tot)

NORMAL VALUES	6.0 - 8.3 %
NORMAL CONTROL	7.5 - 9.0 %
HIGH	9.0 - 10.0 %
LOW	> 10.0 %

### Hemoglobin HbA1C

NORMAL VALUES	4.2 - 6.2 %
NORMAL CONTROL	5.5 - 6.8 %
HIGH	6.8 - 7.6 %
LOW	> 7.6 %

## BLOOD GLUCOSE ESTIMATE

In a test study performed by Dr. D.M. Nathan the glycohemoglobin test was shown to have a linear correlation with mean blood glucose results from patients performing frequent self-monitoring of their blood glucose levels.

## BIBLIOGRAPHY

Trivelli, L.A., Ranney, H.M. and Lai, H.T. New Eng. J. Med. 284, 353 (1971).  
 Goen, B., and Rubenstein, A.H., Diabetologia 15, (1978).  
 Gabby, K.H., Hasty, K., Breslow, J.L., Ellison, R.C., Bunn, H.F., and Gallop, P.M., J. Clin. Endocrinol. Metab. 44, 859 (1977).  
 Bates, H.M., Lab Manag., Vol 16 (Jan 1978).  
 Nathan, D.M., et al., The Clinical Information Value of the Glycosylated Hemoglobin Assay, The New England Journal of Medicine 310, 341-346 (1984).

## SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer

## ANALYTICAL PERFORMANCES

### Interferences

Bilirubin does not interfere with the test.  
 Grossly lipemic or turbidity serum may interfere.  
 The hemoglobin interference is not valuable.

### Linearity

Reaction is linear up to a concentration of 20 %. Samples with values exceeding 20 % must be diluted with saline solution. Multiply, then, the result for diluting factor.

### “Intra-Assay” precision (within-Run)

Determined on 5 samples for each control (N-H) (Normal-High).  
 Results:

MEAN (%)	N = 7.56	H = 11.90
S.D.	0.09	0.10
C.V.%	1.18	0.84

### “Inter Assay” precision (between-Run)

Determined on 5 samples for each control (N-H) for 3 days  
 Results:

MEAN (%)	N = 7.58	H = 11.91
S.D.	0.08	0.09
C.V.%	1.10	0.79

### Analytical sensitivity

The test sensitivity in terms of detection limit is 0.0151%.

### Correlation

A study based comparing this method with a similar method has given a correlation of 98%.