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Fabbricante

Kit for measurement of lipase in serum or plasma UV turbidimetric method

PRINCIPLE

In the presence of lipase the triolein split in monoglyceride and oleic acid. The decrease in turbidity is measured at 340 nm.

REAGENTS

- R1** TRIS buffer pH 8,9 26.0 mmol/l
R2 Sodium deoxycholate 16,7 mmol/l; calcium chloride 0,04 mmol/l; triolein 0,3 mmol/l; colipase 4 mg/l
R3 Standard lipase

SAMPLE

- Serum or heparinized plasma.

Note

- Lipase is stable in the sample for 5 days at 2-8°C or 24 hours at 15-25°C.

REFERENCE VALUES

SERUM	up to 190 U/l (25/30/37°C)
--------------	----------------------------

References values are considered indicatives since each laboratory should establish references ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

PREPARATION OF REAGENTS

- Reagent **R1** is ready for use.
 - **R2**: reconstitute the contents of one vial of **R2** with 2.5 ml of **R1**.
 - **R3**: dissolve the contents of one vial of **R3** in 3.0 ml of redistilled water, swirling gently for 30 minutes before use.
- Keep out the reagents from refrigerator only for the use and recap them immediately.

STORAGE AND STABILITY

- Reagent **R1** is stable until expiry date when stored at 2-8°C.
- Reagent **R2** is stable for 14 days at 2-8°C or for 5 days at 15-25°C.
- Reagent **R3** is stable for 5 days at 2-8°C.

NOTE

- The kit, according to this method, must be used in manual procedures. About automatic using follow specific applications.
- Avoid direct light, contamination and evaporation.
- The volumes in the procedure can be changed proportionally.
- In case of complaint or quality control request, refer to the lot number on the package or the lot number on the singles vials.

AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

Suggested serum:

REF 20350 Precise Norm **REF** 20360 Precise Path

PRECAUTION IN USE

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of components is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances. However the reagent should be handled with caution, according to good laboratory practice.

Caution: the reagents contain Sodium Azide (0.095%) as preservative. Avoid swallowing and contacting with skin, eyes and mucous membranes.

WASTE MANAGEMENT

Please refer to the local legal requirements.

PROCEDURE

Wavelength	λ: 340 (365-334) nm
Working temperature	25/30/37°C
Optical path	1 cm
Reaction	turbidimetric

Bring the reagents at 15-25°C before using them.

- PROCEDURE

Pipette into cuvette:	STANDARD	SAMPLE
REAGENT	2.5 ml	2.5 ml
SAMPLE	--	0.1 ml
STANDARD	0.1 ml	--

Mix. Avoid the formation of foam. Read absorbance A1 of standard and sample after 4'. After another 5' read absorbance A2 of standard and sample. (ΔA of sample or standard = $A1 - A2$).

NOTE

In rare cases, a patient's serum may give an increase in absorbance rather than a decrease. The lipase activity of these samples usually falls within the normal range. Extremely high lipase activities can lead to considerable substrate consumptions, with A1 being less than 0.500. In such cases, dilute the samples 1+9 with 0.9% NaCl solution and repeat the assay. Multiply the result by 10. It is preferable to use disposable cuvettes. Glass cuvettes should be cleaned thoroughly especially after being used for triglyceride or cholesterol assays.

CALCULATION

$$\text{Factor} = \text{activity standard} / \Delta A \text{ standard}$$

$$\text{Sample lipase activity} = \text{factor} \times \Delta A \text{ sample}$$

ANALYTICAL PERFORMANCES

Interferences

Haemolysis interferes with the assay.

Linearity

Reaction is linear up to a concentration of 700 U/l. Samples with values exceeding 700 U/l must be diluted with saline solution. Multiply, then, the result for diluting factor.

"Intra-Assay" precision (within-Run)

Determined on 30 samples for each control (L-N-H) (Low-Normal-High). Results:

$$C.V.\% \quad (U/l) \quad L = 104.1 \quad N = 165.0 \quad H = 593.1$$

Analytical sensitivity

The test sensitivity in terms of detection limit is 22.82 U/l.

Correlation

A study based comparing this method with a similar method on 21 samples has given a correlating factor $r = 0.99$

BIBLIOGRAPHY

Lott, J.A., et al. Clin. Chem. 1986; 32: 1920.
 Ziegenhorn J., et al. Clin. Chem. 1979; 25: 1067
 Weisshaar, H.D., et al.. (1981). Dtsch. Med. Wschr. 106: 239.

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 CE regulation)

Storing temperature limits

In vitro medical device

Producer