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Fabbricante

Kit for measurement of total protein in urine and CSF Pyrogallol red colorimetric method

REACTION PRINCIPLE

Protein react in acid solution with pyrogallol red and molybdate to form a coloured complex. The intensity of the colour formed is proportional to the protein concentration in the sample.

REAGENTS

R1 Pyrogallol red 0.06 mmol/l; sodium molybdate 0.04 mmol/l Sodium dodecyl sulphate 0.07 mmol/l; Succinate Buffer pH 2.5 0.05 mmol/l.
R2 Albumin standard 100 mg/dl

SAMPLE

- Urine
- Cerebrospinal fluid (CSF)
- Stability: 3 days at 2-8 °C

REFERENCE VALUES

| | |
|-------------------|----------------------|
| Proteins in urine | < 120 mg/24 h |
| Proteins in CSF | < 40 mg/dl (0.4 g/l) |

References values are considered indicatives since each laboratory should establish references ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

PREPARATION AND STABILITY OF WORKING SOLUTION

Reagents are liquid and ready to use.
Bring the calibrator at room temperature before use.

STORAGE AND STABILITY OF WORKING SOLUTION

Store the reagent at room temperature and in darkness, calibrator at 2-8°C.
All the components are stable until the stated expiration date if stored tightly closed to their own temperature.

AUXILIARY EQUIPMENT

Materials not included in the kit: diluents solutions, laboratory glassware, disposable tips, photometers and calibrators

QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

PRECAUTIONS AND WARNINGS

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of components is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances.
However the reagent should be handled with caution, according to good laboratory practice.

WASTE MANAGEMENT

Please refer to the local legal requirements.

PROCEDURE

| | |
|---------------------|-------------|
| Wavelength | λ: 600 nm |
| Working temperature | 15-25°C |
| Optical path | 1 cm |
| Reaction | "end point" |

- PROCEDURE

- Adjust the instrument to zero with distilled water

- Pipette into a cuvette:

| | BLANK | STANDARD | SAMPLE |
|---------------------|-------|----------|--------|
| REAGENT (ml) | 1.0 | 1.0 | 1.0 |
| STANDARD (μl) | -- | 20 | -- |
| SAMPLE (μl) | -- | -- | 20 |
| DISTILLED WATER(μl) | 20 | -- | -- |

Mix and incubate for 5' at Room temperature. Read the absorbance of the sample (EC) and standard (ESTD) against the blank. The colour is stable for at least 30'.

CALCULATION

Urine

$EC/ESTD \times 100$ (standard value) \times vol. (dl) urine 24 h = mg protein/24 h

CSF

$EC/ESTD \times 100$ (standard value) = mg/dl protein in the sample

PERFORMANCE CHARACTERISTICS

Linearity

Reaction is linear up to a concentration of 250 mg/dl. If the results obtained were greater than linearity limit, dilute the sample 1:2 with distilled water and multiply the result by 2.

Sensitivity

Values less than 2.5 mg/dl give non - reproducible results

Precision

The tests have confirmed a CV between 2 and 7 %

Interferences

With samples very cloudy or haemolysed make a sample blank also: add 10 μl of sample to 1 ml of saline solution and read the absorbance at wavelength of 600 nm against saline solution.

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