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 **Fabbricante**

Kit for measurement of lipase in serum or plasma

Colorimetric determination

PRINCIPLE

At alkaline pH the lipase dye substrate 1,2-O-dilauryl-rac-glycero-3-glutaric acid-(6'-methylresorufin) ester is cleaved by the catalytic action of pancreatic lipase to form 1,2-O-dilauryl-rac-glycerol an unstable compound glutaric acid (6-metil-resorufin) ester. In that alkaline solution, this compound degrades to gutaric acid and methylresorufin. The colour intensity of the red dye formed is directly proportional to the lipase activity. This activity can be measured at 578 nm.

REAGENTS

R1 Bicine buffer pH 8.0 20 mmol; Colipase ≥ 1 mg/l; sodium desoxycolate 1.8 mM; calcium chloride 12 mM.

R2 Tartrate Buffer pH 4; 1,2-O-dilauryl-rac-glycero-3-glutaric acid-(6'-methylresorufin)estere 0.27mM; Taurodeoxycholate (9 mM); Surfactants and preservatives.

R3 Standard lipase

SAMPLE

Serum or heparinised plasma free of hemolysis. Do not use plasma with EDTA.

The lipase activity in serum is stable for 4 days at 2-8°C and 24 H at room temperature (20-25°C).

REFERENCE VALUES

SERUM and PLASMA	< 60 U/l
-------------------------	----------

References values are considered indicatives since each laboratory should establish references ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

PREPARATION OF REAGENTS

- Reagent **R1** is ready to use.
- Reagent **R2** is ready to use.
- Calibrator **R3** is lyophilized.

Dissolve the contents of one vial of R3 in 1.0 ml of redistilled water, swirling gently until complete dissolution. The concentration is reported on the label of the vial.

STORAGE AND STABILITY

Store the Kit at 2-8°C

All the components are stable until the stated expiration day if stored tightly closed and refrigerated. Once opened the standard is stable for 15 days at 2-8°C.

QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

Suggested serum:

REF 20350 Precise Norm REF 20360 Precise Path

PRECAUTION IN USE

The reagents contain sodium azide at 0.95%. Handle with care. The disposal of the residues has to be made according to legal local regulation.

PROCEDURE

Wavelength λ : 578 nm
 Working temperature 37°C
 Optical path 1 cm
 Reaction kinetic colorimetric

Pipette into cuvettes	Blank	Sample	Calibrator
Distilled water	10 μ l	-	-
Sample	-	10 μ l	-
Calibrator	-	-	10 μ l
R1	1 ml	1ml	1ml
Mix accurately and let stand for 5 min at 37°C then add reagent R2:			
R2	600 μ l	600 μ l	600 μ l

Mix well. After 60 seconds read (E1). After 90 seconds more, read again (E2).

CALCULATION

$\Delta A = A_2 - A_1$

$$\text{Lipase (U/l)} = \frac{(\Delta A \text{ sample}) - (\Delta A \text{ blank})}{(\Delta A \text{ standard}) - (\Delta A \text{ blank})} \times \text{standard conc.}$$

To express the result in μ K/L, multiply the U/L by the factor 0.0167

To express the result in colorimetric unit, multiply the result by the value 3,2.

PERFORMANCES CHARACTERISTICS

Linearity

Up to 300 U/l. For higher concentrations, dilute the sample 1/5 with saline (NaCl 0.9%) and assay once again. Multiply the final result by 5. The analytical performance characteristics of the product depend both on the reagent and the reading system used, manual or automatic.

"Intra-Assay" precision (within-Run): CV coefficient 2.28%

"Inter-Assay" precision (between run): CV coefficient 2.64%

Interferences

Hemolysis interferes with the assay. Concentrations of Hemoglobin >400 mg/dl and Bilirubin >60 mg/dl can interfere with the assay.

BIBLIOGRAFIA

Rick, W.,(1969), Zeitschrift Clin. Chem. Clin.Biochem.,7, 530-536.

Ziegenhorn, J.et al., (1979), Clin. Chem.,25, 1067

Neumann,U., Kaspar, P.,Ziegenhorn, J., (1984) Meth. Enz. Anal (3rd edition), 26-34.

Lott, J.A. et al., (1986). Clin Chem., 32, 1290-1302.

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 CE regulation)

Storing temperature limits

In vitro medical device