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Fabbricante

## Kit for measurement of iron in serum and plasma

### Colorimetric cromazurole method (CAB)

#### PRINCIPLE

Iron, in ionic trivalent form, reacts with cromazurole (CAB) and cetyl-trimethyl-ammonium-bromide (CTAB) to give a strongly coloured compound. The colour intensity is proportional to the iron concentration in the sample.

#### REAGENTS

**R1** Sodium acetate 450.0 mmol/l; sodium chloride 820.0 mmol/l;  
 mordant blue 0.15 mmol/l; acetic acid 225.0 mmol/l  
**R2** Iron standard 100 µg/dl (17.9 µmol/l)

#### SAMPLE

- Serum or heparinized plasma.

#### Note

- Do not use samples with haemolysis and lipoemic. Use only heparine salts as anticoagulants.  
 - The iron is stable in the samples up to 7 days at 2-8°C.

#### REFERENCE VALUES

Serum or plasma	Men	59 -158 µg/dl	(10.6 - 28.2 µmol/l)
	Women	37 -145 µg/dl	(6.5 - 26 µmol/l)

References values are considered indicatives since each laboratory should establish references ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

#### STORAGE AND STABILITY

- Store the kit at 15-25°C.  
 - After opening, the vials **R1** and **STD** are stable up to the expiry date if recapped immediately and protected from contamination, evaporation, direct light, and stored at the correct temperature.  
 - Working solution stability (**R1**): up to the expiry date at 15-25°C.

#### PREPARATION OF REAGENTS

Reagents are liquid and ready to use.

#### NOTE

- The kit, according to this method, must be used in manual procedures. About automatic using follow specific applications.  
 - Evaluate carefully the results if reagent R1 absorbance is > 0.250 at 620 nm.  
 - Avoid direct light, contamination and evaporation.  
 - The volumes in the procedure can be changed proportionally.  
 - It's advisable to use glassware washed by HCl 1N solution and distilled water and disposable test tubes to eliminate any iron contamination.  
 - In case of complaint or quality control request, refer to the lot number on the package or the lot number on the singles vials.

#### AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

#### QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

#### Suggested serum:

**REF 20350** Precise Norm                      **REF 20360** Precise Path

#### PRECAUTION IN USE

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of components is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances.  
 However the reagent should be handled with caution, according to good laboratory practice.

#### WASTE MANAGEMENT

Please refer to the local legal requirements.

#### PROCEDURE

Wavelength λ: 620 nm  
 Working temperature 37°C  
 Optical path 1 cm  
 Reaction "end point"

#### - PROCEDURE

	BLANK	STD	SAMPLE
REAGENT R1	1000 µl	1000 µl	1000 µl
DEIONIZED WATER	50 µl	--	--
SAMPLE	--	--	50 µl
STANDARD	--	50 µl	--

Mix, then incubate 5' at 37°C. Measure the absorbance of sample (EC) and standard (ES) against the reagent blank .

#### CALCULATION

$$\text{Iron } [\mu\text{g/dl}] = \text{EC/ES} \times \text{Conc. STD}$$

#### CONVERSION FACTOR

$$\text{Iron } [\mu\text{g/dl}] \times 0.1791 = \text{Iron } [\mu\text{mol/l}]$$

#### ANALYTICAL PERFORMANCES

##### Interferences

Bilirubin does not interfere up to concentration of 20 mg/dl.  
 Triglycerides always interfere also at minimum concentrations.

##### Linearity

Reaction is linear up to a concentration of 600 µg/dl (107.46 µmol/l) with a range of 4.5-600 µg/dl (0.80-107.46 µmol/l). Samples with values exceeding 600 µg/dl must be diluted with saline solution. Multiply, then, the result for diluting factor.

##### "Intra-Assay" precision (within-Run)

Determined on 30 samples for each control (L-N-H) (Low-Normal-High).

Results:

MEAN	[µg/dl]	L = 73.23	N = 143.87	H = 201.50
S.D.		2.94	4.70	17.92
C.V.%		4.02	3.26	8.90

##### "Inter-Assay" precision (between-Run)

Determined on 15 samples for each control (L-H-N) for 3 days.

Results:

MEAN	[µg/dl]	L = 73.71	N = 145.76	H = 203.69
S.D.		2.67	3.93	3.66
C.V.%		3.63	2.69	1.80

##### Analytical sensitivity

The test sensitivity in terms of detection limit is 4.5 µg/dl (0.80 µmol/l).

##### Correlation

A study based comparing this method with a similar method on 21 samples has given a correlating factor **r = 0.99**

#### BIBLIOGRAPHY

Kaplan LA, Pesce AJ: "Clinical Chemistry", Mosby Ed. (1996).

#### SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer