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### BIBLIOGRAFIA

Trinder P: J. Clin. Path., 22, 158 (1969).

Allain CC, Poon LS, Chan CS, Richmond W, Fu PC: Clin. Chem., 20, 470 (1974).

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Fabbricante

## Kit for measurement of iron in serum or plasma Colorimetric enzymatic ferene method

### PRINCIPLE

The iron released from transferrin complex at pH acid, reacts, once reduced to ferrous state, with Ferene S giving a stable complex which absorbs at 600 nm. The intensity of colour is directly proportional to iron concentration in the sample.

### REAGENTS

- R1** Acetate buffer 1.3 mol/l; thiourea 65.0 mmol/l;  
hydroxylamine sulphate 60.0 mmol/l; surfactants  
**R2** Ferene S 0.65 mmol/l;  
**R3** Iron standard 100 µg/dl (17.9 µmol/l)

### SAMPLE

- Serum or heparinized plasma.

### Note

- Do not use samples with haemolysis and lipoemic. Use only heparine salts as anticoagulants.
- The iron is stable in the samples up to 7 days at 2-8°C.

### REFERENCE VALUES

Serum or plasma	Men	59 -158 µg/dl	(10.6 - 28.2 µmol/l)
	Women	37 -145 µg/dl	(6.5 - 26 µmol/l)

References values are considered indicatives since each laboratory should establish references ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

### STORAGE AND STABILITY

- Store the kit at 2-8°C. Do not freeze the reagents.
- After opening, the vials **R1**, **R2** and **STD** are stable up to the expiry date if recapped immediately and protected from contamination, evaporation, direct light, and stored at the correct temperature.
- Working solution stability (**R1+ R2**): 60 days at 2-8°C.

### PREPARATION OF REAGENTS

Reagents are liquid and ready to use. About using as monoreagent ("sample-starter" procedure) add to every 4 ml of **R1** reagent, 1 ml of **R2** reagent. Keep out the reagents from refrigerator only for the use and recap them immediately.

### NOTE

- The kit, according to this method, must be used in manual procedures. About automatic using follow specific applications.
- Evaluate carefully the results if working reagent absorbance is > 0.300 at 620 nm.
- Avoid direct light, contamination and evaporation.
- The volumes in the procedure can be changed proportionally.
- In case of complaint or quality control request, refer to the lot number on the package or the lot number on the singles vials.

### AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

### QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

### Suggested serum:

REF 20350 Precise Norm REF 20360 Precise Path

### PRECAUTION IN USE

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of components is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances. However the reagent should be handled with caution, according to good laboratory practice.

### WASTE MANAGEMENT

Please refer to the local legal requirements.

### PROCEDURE

Wavelength	λ: 600 nm
Working temperature	37°C
Optical path	1 cm
Reaction	"end point"

Bring the reagents at 15-25°C before using them.

### - MONOREAGENT PROCEDURE "sample starter"

	BLANK	STD	SAMPLE
WORKING REAGENT	1000 µl	1000 µl	1000 µl
DISTILLED WATER	200 µl	--	--
SAMPLE	--	--	200 µl
STANDARD	--	200 µl	--

Mix, then incubate for 5' at 37°C. Measure the absorbance of sample (EC) and standard (ESTD) against the reagent blank.

### - BIREAGENT PROCEDURE "substrate starter"

	BLANK	STD	SAMPLE
REAGENT R1	800 µl	800 µl	800 µl
DISTILLED WATER	200 µl	--	--
SAMPLE	--	--	200 µl
STANDARD	--	200 µl	--

Mix, incubate at 37°C for 1' and then add:

	BLANK	STD	SAMPLE
REAGENT R2	200 µl	200 µl	200 µl

Mix, then incubate for 5' at 37°C. Measure the absorbance of sample (EC) and standard (ESTD) against the reagent blank.

### CALCULATION

$$EC/ES \times \text{conc. STD} = \mu\text{g } (\mu\text{mol/l}) \text{ of iron /dl (L) of sample}$$

### CONVERSION FACTOR

$$\text{Iron } [\mu\text{g/dl}] \times 0.1791 = \text{Iron } [\mu\text{mol/l}]$$

### ANALYTICAL PERFORMANCES

#### Interferences

Bilirubin interferes up to concentration of 20 mg/dl.  
Triglycerides always interfere also at minimum concentrations.

#### Linearity

Reaction is linear up to a concentration of 600 µg/dl (107.46 µmol/l) with a range of 5-600 µg/dl (0.89-107.46 µmol/l). Samples with values exceeding 600 µg/dl must be diluted with saline solution. Multiply, then, the result for diluting factor.

#### "Intra-Assay" precision (within-Run)

Determined on 30 samples for each control (L-N-H) (Low-Normal-High).

Results:

MEAN	[µg/dl]	L = 70.67	N = 152.43	H = 185.20
S.D.		2.51	6.20	3.52
C.V.%		3.55	4.07	1.90

#### "Inter-Assay" precision (between-Run)

Determined on 15 samples for each control (L-H-N) for 3 days.

Results:

MEAN	[µg/dl]	L = 72.00	N = 132.56	H = 163.43
S.D.		1.79	3.23	3.58
C.V.%		2.48	2.43	2.22

#### Analytical sensitivity

The test sensitivity in terms of detection limit is 5 µg/dl (0.89 µmol/l).

#### Correlation

A study based comparing this method with a similar method on 21 samples has given a correlating factor **r = 0.99**

#### BIBLIOGRAPHY

- Trinder P: J. Clin. Path., 22, 158 (1969).  
Allain CC, Poon LS, Chan CS, Richmond W, Fu PC: Clin. Chem., 20, 470 (1974).  
Talke H, Schubert GE: Klin. Wchens., 43, 174 (1965).  
Kaplan LA, Pesce AJ: "Clinical Chemistry", Mosby Ed. (1996).

#### SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 CE regulation)

Storing temperature limits

In vitro medical device

Producer