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### BIBLIOGRAFIA

- Baadenhuinjsen H et al. Modification in Ramsay's method for correct measurement of total iron-binding capacity. Clin. Chim 1988; (175): 9-16.  
Perrotta G. Iron and iron-binding capacity. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1063-1065.  
Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.  
Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.  
Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.  
Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

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Fabbricante

## Total iron-binding capacity (TIBC) – Saturating and precipitating reagent

### PRINCIPLE

Serum transferrin is saturated with an excess of Fe<sup>3+</sup> and the unbound portion is precipitated with magnesium carbonate. The total amount of iron is then determined. The difference between the total iron-binding capacity (TIBC) and initial serum iron (SI) yields the unsaturated iron-binding capacity (UIBC).

### REAGENTS

- R1** Saturating solution: iron solution 500 µg/dl  
**R2** Precipitating reagent: magnesium carbonate

### SAMPLE

- Serum or heparinized plasma.

### Note

- Use samples free of haemolysis.
- Separate from cells as rapidly as possible.
- The iron is stable in the samples up to 7 days at 2-8°C.

### REFERENCE VALUES

<b>SERUM or PLASMA:</b>	200 - 400 µg/dl	36 - 72 µmol/l
-------------------------	-----------------	----------------

Reference values are considered indicative since each laboratory should establish reference ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

### STORAGE AND STABILITY

- All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8 °C, protected from light and contaminations prevented during their use.
- Do not use reagents over the expiration date.

### Note

Presence of particles and turbidity indicates signs of reagent deterioration.

### PREPARATION OF REAGENTS

The reagents are ready to use.

### NOTE

- The kit, according to this method, must be used in manual procedures. About automatic use follow specific applications.
- Avoid direct light, contamination and evaporation.
- In case of complaint or quality control request, refer to the lot number on the package or the lot number on the single vials.

### AUXILIARY EQUIPMENT

Materials not included in the kit: iron determination reagent, diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

### Note

It is recommended to use disposable material. If glassware is used the material should be soaking for 6 hours in diluted HCl (20% v/v) and then thoroughly rinsed with distilled water and dried before use.

### QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

### Suggested serum:

**REF** 20350 Precise Norm                      **REF** 20360 Precise Path

### PRECAUTION IN USE

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of components is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances.

However the reagent should be handled with caution, according to good laboratory practice.

**Caution:** the reagents contain Sodium Azide (0.095%) as preservative. Avoid swallowing and contacting with skin, eyes and mucous membranes.

### WASTE MANAGEMENT

Please refer to the local legal requirements.

### PROCEDURE

Pipette into the tubes:

<b>SAMPLE</b>	0.5 ml
<b>SATURATING SOLUTION</b>	1.0 ml
Mix well and incubate for 10' at 15-25°C, then add to each tube:	
<b>PRECIPITATING REAGENT</b>	(*) 3 spoonful
Mix well and incubate for 10' at 15-25°C. Centrifuge 15' at 3000/rpm. Collect the supernatant carefully and measure the iron concentration (**) with iron reagent kit.	

(\*) Powder: dispense using the enclosed spoon (dosage: aprox. 70 mg).

(\*\*) The supernatant is stable up to 1 hour at 15-25°C. If appear turbid, centrifuge again.

### CALCULATION

The calculation are indicated in the Iron insert determination.

$$\text{TIBC} = \text{iron concentration in the supernatant} \times 3 \text{ (dilution factor)}$$

### ANALYTICAL PERFORMANCES

#### Interferences

Haemolyzed samples are rejected, since erythrocytes contain iron and therefore falsely elevate the serum results.

A list of drugs and other interfering substances with iron determination has been reported by Young et. al. (see "Bibliography").

#### Linearity

Reaction is linear up to a concentration of 1000 µg/dl with a range of 3 µg/dl - 1000 µg/dl. Samples with values exceeding 1000 µg/dl must be diluted with saline solution. Multiply, then, the result for diluting factor.

#### "Intra-Assay" precision (within-Run)

Determined on 20 samples for each control (N-H) (Normal-High).

Results:

MEAN (µg/dl)	N = 367	H = 567
S.D.	2.82	9.43
C.V.%	0.77	1.66

#### "Inter Assay" precision (between-Run)

Determined on 20 samples for each control (N-H).

Results:

MEAN (µg/dl)	N = 359	H = 565
S.D.	7.16	9.46
C.V.%	1.99	1.67

#### Analytical sensitivity

The test sensitivity in terms of detection limit is 1 µg/dl.

### BIBLIOGRAPHY

- Baadenhuinjen H et al. Modification in Ramsay's method for correct measurement of total iron-binding capacity. Clin. Chim 1988; (175): 9-16.  
Perrotta G. Iron and iron-binding capacity. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984; 1063-1065.  
Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.  
Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.  
Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.  
Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.

### SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer