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Fabbricante

Kit for measurement of acid phosphatase in serum Kinetic naphthylphosphate method

PRINCIPLE

Acid phosphatase (ACP), in acid environment, catalyzes the alpha-naphthylphosphate (alpha-NP) hydrolysis in alpha-naphtol and phosphate. Alpha-naphtol reacts with diazo-2-chloro-5-toluene (Fast Red TR Salt) and forms an azo dye compound with consequent absorbance increasing at 405 nm which is proportional to the total acid phosphatase activity in the tested sample. Acid phosphatase of prostatic origin is inhibited by tartrate and is measured by indirect method for the difference between total (ACP) and non-prostatic acid phosphatase.

REAGENTS

- R1** Citrate buffer 150.0 mmol/l;
R2 1-Naphthylphosphate 10.0 mmol/l; fast red TR salt 2.2 mmol/l
R3 Sodium tartrate 530.0 mmol/l

SAMPLE

- Serum.

Note

- Do not use plasma.
- Do not use samples with haemolysis. The withdrawal must be made, preferably, after 8 hours of fasting.
- In the samples previously acidified the acid phosphatase is stable up to 6 days at 2-8°C.

REFERENCE VALUES

Serum [U/l]	MAN	WOMAN	PROSTATIC
37°C up to	5.4	4.2	1.7

References values are considered indicatives since each laboratory should establish references ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

PREPARATION OF REAGENTS

- *Total acid phosphatase*: dissolve the content of **R2** vial with 3 ml of **R1** solution and mark **A**.
 - *Non-prostatic acid phosphatase*: as preceding point, then add 0.03 ml of **R3** sodium tartrate and mark **B**.
- Keep out the reagents from refrigerator only for the use and recap them immediately.

STORAGE AND STABILITY

- Store the kit at 2-8°C. Do not freeze the reagents.
- Reconstituted reagent stability: 2 days at 2-8°C – 6 hours at 15-25°C.
- Stability of reagent with R3: 2 days at 2-8°C – 6 hours at 15-25°C.

NOTE

- The kit, according to this method, must be used in manual procedures. About automatic using follow specific applications.
- Evaluate carefully the results if working reagent absorbance is > 0.300 at 405 nm.
- Avoid direct light, contamination and evaporation.
- The volumes in the procedure can be changed proportionally.
- In case of complaint or quality control request, refer to the lot number on the package or the lot number on the singles vials.

AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

Suggested serum:

REF 20350 Precise Norm **REF** 20360 Precise Path

PRECAUTION IN USE

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of components is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances. However the reagent should be handled with caution, according to good laboratory practice.

WASTE MANAGEMENT

Please refer to the local legal requirements.

PROCEDURE

Wavelength λ : 405 nm
 Working temperature 37°C
 Optical path 1 cm
 Reaction kinetic (increasing)

Bring the reagents at 15-25°C before using them.

- "Total acid phosphatase" PROCEDURE

SAMPLE	100 μ l
REAGENT A	1000 μ l

Mix, then after 5' at 37°C measure the absorbance of sample (EC) against distilled water every minute during 3 minutes. Calculate the average absorbance difference per minute ($\Delta E/min$).

- "Non-prostatic acid phosphatase" PROCEDURE

SAMPLE	100 μ l
REAGENT B	1000 μ l

Mix, then after 5' at 37°C measure the absorbance of sample (EC) against distilled water every minute during 3 minutes. Calculate the average absorbance difference per minute ($\Delta E/min$).

CALCULATION

$$\begin{aligned} \text{Total acid phosphatase [U/l]} &= \Delta E/min \times 750 \\ \text{Non-prostatic acid phosphatase [U/l]} &= \Delta E/min. \times 750 \\ \text{Prostatic phosphatase [U/l]} &= \text{Total acid phosphatase} - \text{Non prostatic} \end{aligned}$$

ANALYTICAL PERFORMANCES

Interferences (ACP)

Triglycerides do not interfere up to concentration of 5 g/l.
 Bilirubin always interferes, also at minimum concentrations.
 Hemoglobin always interferes, also at minimum concentrations.

Interferences (ACP NON-PROSTATIC)

Triglycerides do not interfere up to concentration of 400 mg/dl.
 Bilirubin does not interfere up to concentration of 20 mg/dl.
 Hemoglobin does not interfere up to concentration of 1.5 g/dl.

Linearity (ACP)

Reaction is linear up to a concentration of 150 U/l. Samples with values exceeding 150 U/l must be diluted with saline solution. Multiply, then, the result for diluting factor.

Linearity (ACP NON-PROSTATIC)

Reaction is linear up to a concentration of 70 U/l. Samples with values exceeding 70 U/l must be diluted with saline solution. Multiply, then, the result for diluting factor.

"Intra-Assay" precision (within-Run)

Determined on 30 samples for each control (L-N-H) (Low-Normal-High).

Results:

ACP

MEAN L = 5.93 N = 12.08 H = 15.97
 S.D. 0.16 0.32 0.47
 C.V.% 2.66 2.64 2.92

"Inter-Assay" precision (between-run)

Determined on 15 samples for each control (L-N-H) for 3 days.

Results:

ACP

MEAN L = 5.93 N = 12.12 H = 16.04
 S.D. 0.15 0.30 0.38
 C.V.% 2.48 2.50 2.35

Correlation

A study based comparing this method with a similar method on 21 samples has given a correlating factor $r = 0.99$

BIBLIOGRAPHY

Z. Klin. Chem. & Klin. Biochem., (1972), 10, 182.
 Szasz G. et al., Z. Kinderheilk., (1971), 111, 233.

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer