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 Fabbricante

Kit for measurement of copper in serum

Colorimetric method without deproteinization

3-5-di-Br-PAESA method

PRINCIPLE

Copper, removed from proteins through serum acidification and reduced to copper(I) copper, reacts with 3-5-di-Br-PAESA giving a coloured complex which intensity is proportional to the copper concentration in the sample.

REAGENTS

R1 : Sodium acetate 160 mol/l; 3-5-Di-Br-PAESA 0.2 mmol/l ;
Surfactant 1%.

STANDARD: copper 100 µg/dl (15.7 µmol/l)

SAMPLE

- Fresh serum

Note

- Do not use samples with haemolysis.
- The copper is stable in the samples up to 8 days at 2-8°C.

REFERENCE VALUES

Man	80 ÷ 140 µg/dl
Woman	80 ÷ 155 µg/dl

Reference values are considered indicative since each laboratory should establish reference ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

STORAGE AND STABILITY

- Store the kit at 15-25°C.
- Reagent R1 will precipitate at a temperature less than 18°C. To keep it at preexistent physical condition warm by water bath (26°C) for 5-6 minutes.
- After opening the reagents are stable up to the expiry date if recapped immediately and protected from contamination, evaporation, direct light, and stored at the correct temperature.

PREPARATION OF REAGENTS

Reagents are liquid and ready to use. It is necessary that all reagents are at room temperature when are used.

NOTE

- The kit, according to this method, must be used in manual procedures. About automatic using follow specific applications.
- Avoid direct light, contamination and evaporation.
- The volumes in the procedure can be changed proportionally.
- In case of complaint or quality control request, refer to the lot number on the package or the lot number on the singles vials.

AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

Suggested serum:

REF 20350 Precise Norm REF 20360 Precise Path

PRECAUTION IN USE

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of components is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances. However the reagent should be handled with caution, according to good laboratory practice.

WASTE MANAGEMENT

Please refer to the local legal requirements.

ASSAY CONDITIONS

Wavelength λ: 580 (570-600) nm
Working temperature 15-25°C
Optical path 1 cm
Reaction "end point" (increasing)

- PROCEDURE

	BLANK	STD	SAMPLE
REAGENT R1	1000 µl	1000 µl	1000 µl
DISTILLED WATER	50 µl	--	--
SAMPLE	--	--	50 µl
STANDARD	--	50 µl	--

Mix carefully and incubate 5 min. at 37°C. Read the absorbance at 580nm. (A) against Blank.

CALCULATION

$$\text{Copper } [\mu\text{g/dl}] = \frac{A_{\text{Sample}}}{A_{\text{STD}}} \times 100$$

$$\text{Copper Units SI } [\mu\text{mol/l}] = [\mu\text{g/dl}] \times 0.157$$

ANALYTICAL PERFORMANCES

Interferences

Hemoglobin does not interfere up to concentration of 50 mg/dl.
Bilirubin does not interfere up to concentration of 10 mg/dl.
Triglycerides does not interfere up to concentration of 500 mg/dl.
Iron does not interfere up to concentration of 300 µg/dl.
Zinc does not interfere up to concentration of 300 µg/dl.

Linearity

Reaction is linear up to a concentration of 500 µg/dl (78.5 µmol/l). Samples with values exceeding 500 µg/dl must be diluted 1:2 with saline (NaCl 0.9%). Run the test again and calculate the result multiplying it for 2.

"Intra-Assay" precision (within-Run)

Determined on 20 samples for each control (N-H) (Normal-High).
Results:

MEAN	[µg/dl]	N = 120.9	H = 269.1
S.D.		2.548	2.844
C.V.%		2.11	1.06

"Inter-Assay" precision (between-Run)

Determined on 20 samples for each control (H-N) for 3 days.
Results:

MEAN	[µg/dl]	N = 121.1	H = 268.5
S.D.		2.58	3.27
C.V.%		2.13	1.22

Analytical sensitivity

The test sensitivity in terms of detection limit is 15 µg/dl.

Correlation

A study based comparing this method with a similar method on 20 samples has given a correlating factor **r = 0.98**

BIBLIOGRAPHY

Abe A., Yamashita S., Noma a., (1989) Clin. Chem. 35, 552-554.

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer