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Fabbricante

Kit for sodium determination in serum or plasma Enzymatic-colorimetric method


Reagent R1/A contains Tris

R36/37/38: Irritating to eyes, respiratory system and skin.
S26/36: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable protective clothing.

PRINCIPLE

This is a kinetic method for measuring serum Na⁺ concentration based on determination of Na⁺-dependent β-galactosidase activity. The serum Na⁺ concentration is measured by monitoring the rate of formation of o-nitrophenol from o-nitrophenyl-β-D-galactopyranoside (oNPG), with use of Na⁺-dependent β-galactosidase, as in the following scheme:

REAGENTS

R1/A (liquid)	Tris buffer pH 8.90
R1/B (lyophilized)	β-galactosidase
R2/A (liquid)	Tris buffer pH 6.20
R2/B (lyophilized)	O-nitrophenyl-β-D-galactopyranoside (o-NPG)
R3	Sodium Standard 140 mEq/l

SAMPLE

- Serum, plasma treated with lithium-heparinate.

Note

- Don't use sodium-EDTA as anticoagulant.

REFERENCE VALUES

Serum	135 - 150 mmol/l (311 - 345 mg/dl)
-------	------------------------------------

References values are considered indicatives since each laboratory should establish references ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

PREPARATION OF REAGENTS

- Let reagents reach working temperature before use.
- **R1:** Reconstitute each vial of R1/B with exactly 15 ml R1/A. Mix by gently swirling till to complete solubilization, avoiding the formation of foam. Wait 5 minutes before use.
- Once reconstituted, reagent R1 can appear cloudy or flocculated. This is not affecting the good performance of the kit: Filter if necessary.
- **R2:** Reconstitute R2/B with exactly 10 ml R2/A. Mix by gently swirling till to complete solubilization, avoiding the formation of foam. Wait 5 minutes before use.

STORAGE AND STABILITY

- Reagents are stable at 2-8°C up to the expiration date on the label.
- Reagent R1 (R1/A+R1/B) reconstituted is stable 2 weeks at 2-8°C.
- Reagent R2 (R2/A+R2/B) reconstituted is stable 4 weeks at 2-8°C.

AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

Suggested serum:

REF 20350 Precise Norm REF 20360 Precise Path

PRECAUTION IN USE

Look at the caution advices indicated nearby the symbol.

WASTE MANAGEMENT

Please refer to the local legal requirements.

PROCEDURE

Wavelength	λ:405 nm
Working temperature	37°C
Optical path	1 cm
Reaction	kinetic

Let reagents reach working temperature before use.

Pipette into test tubes or cuvettes labelled as follows:

	BLANK	SAMPLE	STD
DISTILLED WATER	0.025 ml	--	--
SAMPLE	--	0.025 ml	--
STANDARD	--	--	0.025 ml
REAGENT R1	0.75 ml	0.75 ml	0.75 ml
REAGENT R2	0.25 ml	0.25 ml	0.25 ml

Mix and incubate at the working temperature for 2 minutes. Read initial absorbance at 405 nm and start timer simultaneously. Read again after 1 and 2 minutes. Calculate average of absorbance/minute readings (ΔA/min) for reagent (B), sample (EC) and standard (ST).

CALCULATION

$$\text{mmol sodium/L sample} = \frac{\Delta A/\text{min EC} - \Delta A/\text{min B}}{\Delta A/\text{min ST} - \Delta A/\text{min B}} \times [\text{STD}]^*$$

* Sodium concentration in the standard in mmol/L.

ANALYTICAL PERFORMANCES

Interferences

Triglycerides do not interfere up to concentration of 2500 mg/dl.
 Bilirubin does not interfere up to concentration of 27 mg/dl.

Linearity

Reaction is linear up to a concentration of 180 mmol/l (414 mg/dl) with a range of 80-180 mmol/l (184-414 mg/dl). Samples with values exceeding 180 mmol/l dilute the sample with same distilled water volume and multiply the obtained result for two.

"Intra-Assay" precision (within-Run)

	Level 1	Level 2
MEAN [mM]	120	160
S.D.	2.48	4.59
C.V.%	2.06	2.86

"Inter Assay" precision (between-Run)

	Level 1	Level 2
MEAN [mM]	123	155
S.D.	4.8	7.31
C.V.%	3.90	4.72

Analytical sensitivity

The test sensitivity in terms of detection limit is 80 mmol/l (184 mg/dl).

Correlation with flame photometry

The kit presents a correlation coefficient of 0.985 than flame photometry.

NOTE

- The kit, according to this method, must be used in manual procedures. About automatic using follow specific applications.
- The volumes in the procedure can be changed proportionally.
- In case of complaint or quality control request, refer to the lot number on the package or the lot number on the singles vials.
- Use only water without sodium, potassium and calcium ions.
- Use very clean laboratorial equipment.
- If sodium and potassium are requested together, sodium is determined immediately prior to potassium (2-channel method).
- Carry out standard measure for each sample set.

BIBLIOGRAPHY

Berry M.N. et al., Clin. Chem., 34,2295 (1988)
 Kumar A. et al., Clin. Chem. 34,1709 (1988)

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer