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
Fabbricante

## Kit for sodium determination in serum or plasma

### Enzymatic-colorimetric method

#### PRINCIPLE

This is a kinetic method for measuring serum Na<sup>+</sup> concentration based on determination of Na<sup>+</sup>-dependent β-galactosidase activity. The serum Na<sup>+</sup> concentration is measured by monitoring the rate of formation of o-nitrophenol from o-nitrophenyl-β-D-galactopyranoside (oNPG), with use of Na<sup>+</sup>-dependent β-galactosidase, as in the following scheme:

#### REAGENTS

**R1** Good buffer pH 8.5, β-D-galactosidase 8U/ml, Criptando >0.4mM, ProClin 300 0.02%.

**R2** Good buffer pH 6.5, o-nitrophenyl-β-D-glycoside (o-NPG) >0.5 mM, ProClin 300 0.02%.

**R3** Sodium standard whose concentration is stated on the label.

#### SAMPLE

- Serum.

#### REFERENCE VALUES

Serum	136 - 146 mmol/l (313 - 336 mg/dl)
-------	------------------------------------

References values are considered indicatives since each laboratory should establish references ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

#### PREPARATION AND STABILITY OF REAGENTS

- Reagents are ready to use and stable until the date indicated on the label.  
 - Let reagents reach working temperature before use.

#### AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

#### Recommended calibrator:

REF 20105 Biocal - H human

#### QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

#### Recommended control serum:

REF 20350 Precise Norm human      REF 20360 Precise Path human

#### PRECAUTION IN USE

Look at the caution advices indicated nearby the symbol.

#### WASTE MANAGEMENT

Please refer to the local legal requirements.

#### PROCEDURE

Wavelength	λ:405 nm
Working temperature	37°C
Optical path	1 cm
Reaction	kinetic

Let reagents reach working temperature before use.

Pipette into test tubes or cuvettes labelled as follows:

	BLANK	SAMPLE	STD
REAGENT R1	200 μl	200 μl	200 μl
SAMPLE	--	8 μl	--
DISTILLED WATER	8 μl	--	--
STANDARD	--	--	8 μl
REAGENT R2	100 μl	100 μl	100 μl

Mix and incubate at the working temperature for 7 minutes then read initial absorbance at 405 nm (A1). Read again after 3 minutes (A2). Calculate ΔA/min for reagent (B), sample (EC) and standard (ST).

#### CALCULATION

$$\text{mmol sodium/L sample} = \frac{\Delta A/\text{min EC} - \Delta A/\text{min B}}{\Delta A/\text{min ST} - \Delta A/\text{min B}} \times [\text{STD}]^*$$

\* Sodium concentration in the standard in mmol/L.

#### ANALYTICAL PERFORMANCES

##### Interferences

Triglycerides do not interfere up to concentration of 1000 mg/dl.  
 Bilirubin does not interfere up to concentration of 40 mg/dl.

##### Linearity

Reaction is linear up to a concentration of 180 mmol/l (414 mg/dl) with a range of 80-180 mmol/l (184-414 mg/dl). Samples with values exceeding 180 mmol/l dilute the sample with same distilled water volume and multiply the obtained result for two.

##### "Intra-Assay" precision (within-Run)

		Level 1	Level 2
MEAN	(mM)	128.94	155.84
S.D.	(mM)	1.57	1.72
C.V.%		1.2	1.1

##### "Inter Assay" precision (between-Run)

		Level 1	Level 2
MEAN	(mM)	128.94	155.84
S.D.	(mM)	2.01	2.56
C.V.%		1.56	1.65

##### Analytical sensitivity

The test sensitivity in terms of detection limit is 80 mmol/l (184 mg/dl).

##### Correlation with flame photometry

The kit presents a correlation coefficient of 0.985 than flame photometry.

#### NOTE

- The kit, according to this method, must be used in manual procedures. About automatic using follow specific applications.
- The volumes in the procedure can be changed proportionally.
- In case of complaint or quality control request, refer to the lot number on the package or the lot number on the singles vials.
- Use only water without sodium, potassium and calcium ions.
- Use very clean laboratorial equipment.
- If sodium and potassium are requested together, sodium is determined immediately prior to potassium (2-channel method).
- Carry out standard measure for each sample set.

#### BIBLIOGRAPHY

Berry M.N. et al., Clin. Chem., 34,2295 (1988)  
 Kumar A. et al., Clin. Chem. 34,1709 (1988)

#### SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer