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Fabbricante

Kit for potassium determination in serum and plasma

Turbidimetric method

PRINCIPLE

Potassium ions precipitate with sodium tetraphenylborate giving a stable turbidity, proportional to the sample concentration.

REAGENTS

- R1** Tris buffer 50 mmol/l; tetraphenylborate sodium 102 mmol/l;
Boric acid 250 mmol/l
R2 Potassium Standard 5 mmol/l

SAMPLE

Serum or Heparinized plasma not hemolyzed.

REFERENCE VALUES

Serum	3.6 - 5.5 mmol/l
Plasma	4.0 - 4.8 mmol/l

Reference values are considered indicative since each laboratory should establish reference ranges for its own patient's population. The analytical results should be evaluated with other information coming from patient's clinical story.

PREPARATION OF REAGENTS

Reagent is liquid and ready to use.

STORAGE AND STABILITY

- Store the kit at 2-8°C.
- Reagent, if stored at 2-8°C and far from light, is stable up to expiry date indicated on the label.

NOTE

- The kit, according to this method, must be used in manual procedures. About automatic use follow specific applications.
- Avoid direct light, contamination and evaporation.
- In case of complaint or quality control request, refer to the lot number on the package or the lot number on the single vials.

AUXILIARY EQUIPMENT

Materials not included in the kit: diluent solutions, laboratory glassware, disposable tips, photometers and calibrators.

QUALITY CONTROLS

It's necessary, every time the kit is used, to make the quality controls and to check that values obtained are within the acceptance range provided in the insert.

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PRECAUTION IN USE

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of components is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances.

However the reagent should be handled with caution, according to good laboratory practice.

WASTE MANAGEMENT

Please refer to the local legal requirements.

PROCEDURE

Wavelength	λ: 578 (580) nm
Working temperature	37°C
Optical path	1 cm
Reaction	End point

- PROCEDURE

	BLANK	SAMPLE	CALIBRATOR (Std)
WORKING REAGENT	1.0 ml	1.0 ml	1.0 ml
SAMPLE	---	0.025 ml	---
CALIBRATOR (Std)	---	---	0.025 ml

Reset the instrument with distilled water.
Mix to obtain a turbidity homogeneous. Read after 5-10 minutes.
Shake again before the first reading.

CALCULATION

$$\text{Potassium (mmol/l)} = \text{OD Sample} / \text{OD Standard} \times \text{Conc. STD}$$

ANALYTICAL PERFORMANCES

Interferences

Bilirubin does not interfere up to concentration of 10 mg/dl.

Linearity

Reaction is linear between 0.6 - 10 mmol/l. Samples with values exceeding 10 mmol/l must be diluted with saline solution. Multiply, then, the result for diluting factor.

Precisione "intra-Assay" (NELLA SERIE)

Determined on 30 samples for each control (L-H) (Low-High).
Results:

MEAN (mmol/l)	L = 3.99	H = 6.10
S.D.	0.11	0.08
C.V.%	2.75	1.24

Precisione "inter-Assay" (FRA LE SERIE)

Determined on 15 samples for each control (L-H) for 3 days.
Results:

MEAN (mmol/l)	L = 3.96	H = 6.08
S.D.	0.10	0.08
C.V.%	2.56	1.33

Analytical sensitivity

The test sensitivity in terms of detection limit is 0.6 mmol/l.

Correlation

A study based comparing this method with a similar method on 21 samples has given a correlating factor $r = 0.95$

BIBLIOGRAPHY

- Tietz, N.W., Textbook of Clinical Chemistry, W.B. Saunders Co., Philadelphia, 1986, p. 1172-1175.
Sunderman, F.W. Jr., Sunderman, F.W., Am. Journ. of Clin. Pathology, vol. 29, p. 95, 1958.
Young, D.S., Pestaner, L.C. and Gibberman, V., Clin. chem. 21, 371d (1975).

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer