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Kit for quantitative determination of homocysteine in serum and plasma - Enzymatic method

PRINCIPLE

Homocysteine (HCY) is a thiol-containing amino acid produced by the intracellular demethylation of methionine. Total homocysteine represents the sum of all forms including oxidized protein bound free. In blood it is present more in oxidized form joined to proteins than in free reduced form.




Elevated level of Homocysteine has emerged as an important risk factor in the assessment of cardiovascular disease⁽¹⁻³⁾. Excess HCY in the bloodstream may cause injuries to arterial vessels due to its irritant nature, and result in inflammation and plaque formation, which may eventually cause blockage of blood flow to the heart.

Elevated Homocysteine levels are caused by four major factors, including:

- Genetic deficiencies in enzymes involved in HCY metabolisms such as cystathionine beta-synthase (CBS), methionine synthase (MS), and methylenetetrahydrofolate reductase (MTHFR).
- Nutritional deficiency in B vitamins (B6, B12 and folate)
- Renal failure for effective amino acid clearance
- Drug interactions such as nitric oxide, methotrexate and phenytoin that interfere with HCY metabolisms.

Elevated levels of T(HCY) are also linked with Alzheimer's disease⁽⁴⁾ and osteoporosis⁽⁵⁾.

The enzymatic test for the quantitative Homocysteine determination (HCY) is based on a series of enzymatic reactions causing a decrease in absorbance value due to NADH oxidation to NAD⁺. HCY concentration in the sample is directly proportional to the quantity of NADH converted to NAD⁺ (ΔA_{340nm}). The enzymatic reactions are the following:

* SAM = S-Adenosyl-methionine
**SAH = S-Adenosyl-homocysteine

COMPOSITION

R1: S-adenosylmethionine SAM 0.1 mM; NADH 0.2 mM; TCEP 0.5 mM; 2-oxoglutarate 5.0 mM.

R2: Glutamate dehydrogenase 10 KU/l; SAH hydrolase 3.0 KU/l; Adenosine deaminase 5.0 KU/l; HCY methyltransferase 5.0 KU/l.

CALIBRATORS: 2 Levels 2x1 ml

PREPARATION

The reagents are liquid ready to use.

The reagent used should be clear. It should be discarded and can not be used if it becomes turbid or the initial absorbance is less than 0.5 at 340 nm.

STORAGE AND STABILITY

Store at 2-8°C. Do not freeze the reagents. The reagents are stable up to the expiry date stated on the label if contamination and evaporation are avoided, protected from light. The vials must be open just only for the time to take the reagent, closed immediately.

SAMPLES

Use fresh serum or plasma. It is important to centrifuge blood samples immediately after collection to separate the plasma from the blood cells. If immediate centrifugation is not possible, collected blood specimens should be kept on ice and centrifuged within an hour. Hemolysed or turbid specimens or severely lipemic specimens are not recommended for HCY assay.

After separation of serum/plasma from cells, HCY is stable in the sample for 4 days if stored at room temperature, 4 weeks at 2-8°C, and 12 months at -20°C.

REFERENCE VALUES

Men - Women	≤15 µmol/L
-------------	------------

Each laboratory is recommended to establish a range of normal values for the population in their region.

QUALITY CONTROL

It is recommended to use commercial Quality Control sera (SGM-REF.20370) with known Homocysteine concentration. Check that the values obtained are within the reference range provided.

ANALYTICAL PROCEDURE

Allow the reagents to room temperature before use.

	Blank	CAL 1	CAL 2	Sample
Reagent 1	240 µl	240 µl	240 µl	240 µl
Saline	13 µl	-	-	-
Calibrator 1	-	13 µl	-	-
Calibrator 2	-	-	13 µl	-
Sample	-	-	-	13 µl
Mix and incubate for 5 minutes at 37°C.				
Then add:				
Reagent 2	65 µl	65 µl	65 µl	65 µl
Mix and incubate at 37°C for 2.5 minutes. Read against blank absorbance (A ₁) at 340nm and then read again after 2.5 minuti (A ₂).				

NOTE:

- Reaction volumes can be proportionally changed.
- Samples with values greater than 50 µmol/l should be diluted 1:2 and tested again. Multiply results by 2.

CALCULATIONS

Plot (A₂-A₁) obtained against the Hcys concentrations of both calibrators. Hcys concentration in the sample is calculated by interpolation of its (A₂-A₁) in the linear calibration curve.

AUTOMATIC ANALYZERS

The present kit can be used with any type of Clinical Chemistry automatic analyzer. Applications parameters for the most commons commercially available instruments are available upon request (Hitachi, Olympus, Beckman, Konelab, etc.).

CALIBRATION

Use the calibrator 1 and 2 for establish the calibration curve. Use saline or distilled water as calibration point 0 (0 µmol/l) where applicable or requested..

Assigned values are printed on the vials' label and are lot dependent.

Using only one calibrator (linear calibration) is an approximation of the actual calibration curve (non-linear calibration). Even if it is suggested to use both calibrators for establishing the calibration curve in order to achieve the best accuracy, using a linear calibration curve is possible, but it leads to a sensible reduction of the assay accuracy.

Calibration is stable up to 5 days. Anyway, it is suggested to recalibrate in case of the assay Quality Control gives non acceptable results.

PRECAUTIONS IN USE

The reagents contain inactive components such as preservatives (sodium azide or others), surfactants etc. The total concentration of these components is lower than the limits reported by 67/548/EEC and 88/379/EEC directives about classification, packaging and labelling of dangerous substances. However, the reagents should be handled with caution, avoiding swallowing and contact with skin, eyes and mucous membranes. The use of laboratory reagents according to good laboratory practice is recommended.

WASTE MANAGEMENT

Please refer to the local legal requirements.

ANALYTICAL PERFORMANCES

Sensitivity

Test sensitivity, in terms of limit of detection, is 0.4 µmol/l.

Linearity

The assay is linear up to 50 µmol/l.

Precision

Precision has been evaluated by testing replicates of four samples at different Homocysteine concentration.

The obtained results are reported in the following tables.

"Intra-assay" Precision (within-Run)

Sample	n	Mean (µmol/l)	SD	%CV
Sample 1	40	7.0	0.32	4.57
Sample 2	40	12.0	0.22	1.83
Sample 3	40	15.6	0.47	3.01
Sample 4	40	29.0	0.70	2.41

"Inter-assay" Precision (between-Run)

Campione	n	Mean (µmol/l)	SD	%CV
Sample 1	40	7.0	0.42	6.00
Sample 2	40	12.0	0.59	4.92
Sample 3	40	15.6	0.80	5.12
Sample 4	40	29.0	0.75	2.59

Correlation

A correlation study comparing the present method and a commercial one gave the following results:

$$y = 0.94x + 1.05 \mu\text{mol/l} \quad r = 0.99$$

Interferences

Interferences were tested with a serum sample spiked with various concentrations of substances normally present in the serum. In the following table, concentrations of different substances with interferences less than 10% are reported:

Bilirubin	40 mg/dl
Bilirubin conjugated	40 mg/dl
Triglycerides	1000 mg/dl
Ascorbic acid	10 mmol/l
hemoglobin	500 mg/dl
Cystathionine	0.1 mmol/l

Patients taking methotrexate, carbamazepine, phenytoin, nitrous oxide, anticonvulsants, or 6-azuridine triacetate may have high levels of HCY, due to the interference of this drugs with the Homocysteine metabolism.

Addition of 3-deazadenosine to inhibit HCY production in red cells has been suggested. However, the HCY assay can not use samples containing 3-deazadenosine since it inhibits one of the key enzyme used in the assay.

BIBLIOGRAPHY

1. Eikelboom JW et al. Ann Intern Med 131:363-75 (1999)
2. Scott J, Weir D. Q J Med 89: 561-3 (1996)
3. Nygard O, N. Engl J Med 337(4) : 230-6 (1997)
4. Seshadri S. et al. N. Engl. J. Med. 346:477-483 (2002)
5. Mc Lean R. et al. N. Engl. J. Med. 350: 2042-2049 (2004)
6. Refsum H. Clinical Laboratory News May 2002, pp 2-14
7. Guttormsen AB et al. Am. J. Nutr. 124(10) :1934-41 (1994)
8. Vilaseca et al. Clin. Chem. 43 : 690-692 (1997)
9. Faure –Delanef et al. Am. J. Hum. Genet. 60 : 999-1001 (1997)
10. NCCLS Document, « Procedures for the collection of arterial blood specimens », Appr. Std., 3rd Ed. (1999).
11. EU-Dir 1999/11 Commission Directive of 8 March 1999 adapting to technical progress the principles of good laboratory practice as specified in Council Directive 87/18/EEC.

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storage temperature

In vitro medical device

Producer