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MULTICOMPONENT CALBRATOR SERUM FOR CLINICAL CHEMISTRY ASSAYS. IVD

Product characteristics

BIOCAL H is a human lyophilized calibrator serum. The constituents concentrations and activities have been selected to ensure an optimum calibration for use in manual and automatic analyzers.

Reagents

Human serum joined to biological additives and bacteriostatics agents. The components concentrations/activities are lot-specific. The exact value is indicated in the value sheet.

Calibration

Information available on request.

Preparation

- Reconstitute, by pipetting, with 5.0 ml of bidistilled water at +15-25°C, into vial.
- Replace the stopper and let it stand for 30 minutes out of bright light before use.
- Swirl gently several times to ensure the complete dissolution. Prior to use, mix by inverting the vial. Do not shake the vial as the formation of foam should be avoided.
- It's advisable to leave stand the reconstituted serum for 1 hour at +15-25°C before use.
- Improper handling or storage can compromise the results. Inaccurate reconstitution can cause wrong results.

Storage and stability

The control serum is stable until the expiration date indicated on the label if stored tightly closed at 2-8°C, protected from light and contamination.

- Reconstituted serum stability:

8 hours at	+ 15-25°C
7 days at	+ 2-8°C
30 days at	-25-15°C
- The bicarbonate stability is 8 hours in the closed bottle and 1 hour in the open vial.
- The bilirubin stability, if stored protected from light, is 4 days at +2-8°C. Do not store at +15-25°C and do not freeze.
- Alkaline phosphatase values in the reconstituted serum tend to rise during the stability period.
- For total and prostatic acid phosphatase the material should be stabilised by adding 1 drop (25-30 µl) of 0.7M acetic acid solution to 1 ml of the serum. After stabilisation total and prostatic acid phosphatase are stable for at least 2 hours at +25°C, 2 days at +4°C and 1 month when frozen once at -20°C.

Note

Bacterial contamination of the reconstituted serum can cause reductions in the stability of many components.

Precaution in use

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of components is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances. However the reagent should be handled with caution, according to good laboratory practice. Avoid swallowing and contacting with skin, eyes and mucous membranes.

 Human control is obtained using only blood of donors tested by an FDA method and found non-reactive for HbsAg and negative for antibodies to HIV-1/2 and HCV. However, no known test method can offer complete assurance that products derived from human blood will not transmit infectious diseases, this product should be handled as potentially infectious biological material.

Waste Management

Please refer to the local legal requirements.

SYMBOLS

-  Read instruction for use
-  Biological risk
-  CE mark (requirement of 98/79 CE regulation)
-  Storing temperature limits
-  In vitro medical device
-  Producer