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Fabbricante

Kit for quantitative determination of Mucoprotein (AGP/ α 1-Acid Glycoprotein) in human serum - Immunturbidimetric end-point method

PRINCIPLE

As an early acute phase reactant, AGP is especially useful in monitoring tumor recurrence. Levels are also helpful in differentiating acute phase responses (elevated levels) from estrogen effects (normal or depressed levels). In addition it is an excellent protein to assay along with Haptoglobin in assessing in vivo hemolysis. An elevated AGP level but normal Haptoglobin suggests an acute phase response with mild to moderate in vivo hemolysis.

SAMPLE

Fresh serum.

Note

If the test can not be carried out on the same day, the serum may be stored at 2 - 8°C for 48 hours. If stored for a longer period, the sample should be frozen.

REAGENTS

R1 Buffer: Phosphate buffered saline pH 7.43; Polyethylene glycol 60 g/l; sodium azide 0.95 g/l

R2 Antiserum: Phosphate buffered saline pH 7.43; Polyclonal goat anti-human AGP (variable); sodium azide 0.95 g/l.

Note

According to the application the R1 reagent could be in excess.

PRECAUTIONS AND WARNINGS

- Polyethylene glycol is non biohazardous.

 Each donor unit used in the preparation of the standards and controls was found to be negative for the presence of HIV1 and HIV2 antibodies, as well as for the hepatitis B surface antigen and anti-hepatitis C antibodies, using a method approved by the FDA. However this product should be handled as potentially infectious biological material.

Attention: Reagents contain Sodium azide. Avoid swallowing and contacting with skin, eyes and mucous membranes. Sodium azide has been reported to form lead or copper azide in laboratory plumbing which may explode on percussion. Flush drains with water thoroughly after disposing of fluids containing sodium azide.

REFERENCE VALUES

Men	50 - 130 mg/dl (IFCC)
Women	40 - 120 mg/dl

This range is given for orientation only. Each Laboratory should establish its own reference values.

PREPARATION OF REAGENTS

Reagents are liquid and ready for use.

STORAGE AND STABILITY

- The reagents are stable until expiry date when kept at 2-8°C.
- Stability in the instrument is at least 4 weeks if contamination is avoided.
- Do not freeze.

AUXILIARY EQUIPMENT

Reagents required but not supplied:

- 0.9 g % sodium chloride
- Calibrators and controls: standard low and high and protein control.

QUALITY CONTROL

Run controls every kit use and verify that the value are within the range of reference value.

Suggested control sera:

Protein Control Serum Ref. 23792
 Protein Control Serum Ref. 23795

WASTE MANAGEMENT

Please refer to the local legal requirements.

MANUAL PROCEDURE

Sample/Control: dilute 1:10 in saline 9 g/L.

Reference curve: generate a reference curve by diluting the standard high level 1:10, 1:20, 1:40, 1:80, 1:160 in saline 9 g/L. Use saline 9 g/L as zero point.

Test: Mix 30 μ L diluted samples, standards and control(s) with 900 μ L buffer. Read optical density (OD1) of samples, standards and control(s) at 340 nm. Add 60 μ L of AGP Antiserum. Mix and incubate for 5 minutes at room temperature. Read optical density (OD2) of samples, standards and control(s) at 340 nm.

Calculate Δ OD's, plot a standardcurve and read the concentration of controls and samples.

Note

Automatic procedure are available upon request.

ANALYTICAL PERFORMANCES

The performance characteristics were measured on Cobas Mira analyzer.

Measuring Range : 0 - 300 mg/dL

Detection limit: 4 mg/dl

Hook effect: > 600 mg/dl

Sensitivity: 0.0023 ABS units/concentration unit

Precision %CV

	Low	Medium	High
Intra-Run	4.66	1.14	2.45
Inter-Run		2.55	

Accuracy mg/dl

Control	Assigned value	Measured value
Biorad 1	46 (37-55)	44.4
Biorad 2	110 (88-132)	102.3

Specificity: monospecific.

Interferences:

Haemoglobin does not interfere up to 1000 mg/dl.
 Sodium citrate does not interfere up to 1000 mg/dl.
 Heparin does not interfere up to 50 mg/dl.
 Bilirubin does not interfere up to 20 mg/dl.
 Triglycerides do not interfere up to 2500 mg/dl.

Limitations: none.

Comparison with Nephelometry: $y = 1.1554x - 9.1048$
 $r = 0.9958$

Stability at 4° C: at least 3 years after production.

BIBLIOGRAPHY

Schmid, K. in FW putman, Editor, The plasma Proteins, Vol 1, second edition, Academic Press, new York, 2975, pp184-228
 Jonhson, A.M. et al., J. Clin. Invest., 48 (1969) 2293
 Dati, F. et al., Lab. Med. 13 (1989) 87

SYMBOLS

Read instruction for use

Biological risk

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer