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Fabbricante

Kit for quantitative detection of Transferrin (TRF) in serum Immunoturbidimetric end-point Method

PRINCIPLE

Iron binding protein join the iron preventing iron diseases and and the renal elimination

Higher levels can be detected with iron deficiency ,pregnancy, absumtion of estrogens and lipoidal nefritis.

Lower levels can be detected in ereditary deficence, testosterone absumption, infections, acute flammatory process, tumors, nefritis, bad nutrition.

SAMPLE

Serum

Note: the test must be performed during the day , otherwise the sample must be kept at 2-8° C up to 48 hours or frozen for longer periods.

REAGENTS

R1: Phosphate buffer pH 7,43

R2: Antiserum in saline phosphate solution pH 7,43; goat policlonal antibody anti human transferrin.

Note

According to the application the R1 reagent could be in excess.

PRECAUTION IN USE

 Each donor used for standard and controls was tested and found negative for HIV (1/2) and Epatitis B and C . Moreover, treat as potentially infected material.

Attention: reagents contain Sodium Azide (0.095%) as preservant. Avoid ingestion and contact with skin eyes and mucosa.

Sodium azide can generate explosion in contact with plumb.

REFERENCE VALUES

SERUM:	170 – 340 mg/dl (IFCC)
--------	------------------------

This range is a media. Each laboratory must check its own reference.

STORAGE AND STABILITY

- Reagents are stable up to the expiry date if stored at 2-8°C.
- On board stability is 28 days, if protected from direct light and contamination. Do not freeze.

PREPARATION OF REAGENTS

Reagents are liquid and ready to use.

AUXILIARY EQUIPMENT

- Sodium Chloride 0.9 g %
- Calibrators and controls: standard, low, elevated and proteic control

QUALITY CONTROL

controls must be used each kits verifying that the value are included in the reference range.

Suggested sera :

Protein Control Serum Ref. 23795

Protein Control Serum Ref. 23792

MANUAL PROCEDURE

Samples/controls: dilute 1:10 in saline 9 g/l.

Reference curve: create a reference curve with a dilution of standard 1:10, 1:20, 1:40, 1:80, 1:160 in saline 9 g/l. Use saline solution 9 g/l as zero point.

Test: mix 25 µl of standard diluted, controls and samples with 900 µl of buffer. Read optical density (OD1) of samples, controls standard at 340 nm. Add 60 µl of TRF antibody. Mix and incubate 5 minutes at 15-25°C. Read optical density (OD2) of samples, controls standard at 340 nm. Calculate delta of density ΔOD and trace on the reference curve to read concentration of samples and controls.

ANALYTICAL PERFORMANCES

Performances was evaluated on Cobas Mira.

Range of measure: 0 - 500 mg/dl

Sensitivty: Limit 40 mg/dl

Prozone Effect: > 1400 mg/dl

Sensitivity: 0.00093 ABS

Precision %CV

	Low	Normal	High
Intra-Run	4.78	1.30	0.78
Inter-Run	4.64	2.46	

Accuracy mg/dl

Control	teoric Value	obtained Value
Biorad 1	142 (124-168)	153
Biorad 2	328 (262-394)	354

Specificity: monospecific.

Interferences haemoglobin does not interfere up to 1000 mg/dl, sodium citrate does not interfere up to 1000 mg/dl, heparin does not interfere up to 50 mg/dl, bilirubin does not interfere up to 20 mg/dl, triglycerides do not interfere up to 2500 mg/dl.

Limits: none

Stability at 4° C: 3 years after production.

Correlation with other turbidimetric test: $y = 0.8866x + 40.481$
 $r = 0.9756$

BIBLIOGRAPHY

Dati, F. et al., lab. Med. 13,87 (1989).

SYMBOLS

Read instruction for use

Biological risk

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer