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BIBLIOGRAFIA

Lars-Olof Hanson et al. Current Opinion in Infectious diseases 1997; 10:196 - 201.
M.M. Pepys. The lancet 1981; March 21:653 - 656.
Chetana Vaishnavi, Immunology and Infectious Diseases 196;6 : 139 - 144.
Yoshitsugy Hokama et al. Journal of Clinical Laboratory Status 1987 ; 1 : 15 - 27.
Yamamoto S et al. Veterinary Immunology and Immunopathology 1993; 36:257 - 264.
Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACC Press. 1995.

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Agglutination latex test without predilution

PRINCIPLE

The CRP latex is a slide agglutination test for the qualitative and semiquantitative detection of C-reactive protein (CRP) in human serum.

Latex particles coated with goat IgG anti-human CRP are agglutinated when mixed with samples containing CRP.

CLINICAL SIGNIFICANCE

CRP is an acute phase protein present in normal serum, which increases significantly after most forms of tissue injuries, bacterial and virus infections, inflammation and malignant neoplasia.

During tissue necrosis and inflammation resulting from microbial infections, the CRP concentration can rise up to 300 mg/l in 12-24 hours.

REAGENTS

R1 Latex particles coated with goat IgG anti-human CRP, pH 8.2

R2 Positive control red cap: human serum with a CRP concentration > 20 mg/l (in the kit with controls)

R3 Negative control blue cap: animal serum (in the kit with controls)

- Slide with 6 divisions

- Stirrer

SAMPLES

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C.

Samples with presence of fibrin should be centrifuged before testing.

Do not use highly hemolyzed or lipemic samples.

PRECAUTION IN USE

 Components from human origin have been tested and found to be negative for the presence of HbsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.

Caution: the reagents contain Sodium Azide (0.095%) as preservative. Avoid swallowing and contacting with skin, eyes and mucous membranes.

CALIBRATION

The CRP latex sensitivity is calibrated to the reference Material CRM 470/RPPHS.

STORAGE AND STABILITY

All the kit components are ready to use, and will remain stable until the expiration date printed on the label, when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not freeze: frozen reagents could change the functionality of the test

Reagents deterioration: presence of particles and turbidity.

ADDITIONAL EQUIPMENT

Mechanical rotator with adjustable speed at 80/100 r.p.m.

QUALITATIVE PROCEDURE

1. Allow the reagents and samples to reach room temperature. The sensitivity of the test may be reduced at low temperatures.
2. Place 50µl of the sample (note 1) and one drop of each positive and negative controls into separate circles on the slide test.
3. Swirl the CRP latex reagent gently before using and add one drop (50 µl) next to the sample to be tested.
4. Mix the drop with a stirrer, spreading them over the entire surface of the circle. Use different stirrer for each sample.
5. Place the slide on a mechanical rotator at 80-100 r.p.m. for 2 minutes. False positive results could appear if the test is read later than two minutes.

PROCEDIMENTO SEMIQUANTITATIVO

1. Make serial two fold dilutions of the sample in 9 g/l saline solution.
2. Proceed for each dilution as in the quantitative method.

READING AND INTERPRETATION

Examine macroscopically the presence or absence of visible agglutination immediately after removing the slide from the rotator.

The presence of agglutination indicates a CRP concentration equal or greater than 6 mg/l (note 2 and 3).

The titer in the semi-quantitative method, is defined as the highest dilution showing a positive result.

CALCULATIONS

The approximate CRP concentration in the patient sample is calculated as follows:

$$6 \times \text{CRP titer} = \text{mg/l}$$

QUALITY CONTROL

Positive and negative controls are recommended to monitor the performance of the procedure, as well as a comparative pattern for a better result interpretation.

REFERENCE VALUES

Adults	< 6 mg/l
--------	----------

Each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

1. Analytical sensitivity: 6 (5-10) mg/l under the described assay conditions.
2. Prozone effect: no prozone effect was detected up to 1600 mg/l.
3. Diagnostic sensitivity: 95,6%
4. Diagnostic specificity: 96,2%

INTERFERENCES

Hemoglobin does not interfere up to concentration of 10 g/l.

Bilirubin does not interfere up to concentration of 20 mg/dl.

Lipemia does not interfere up to concentration of 10 g/l.

Rheumatoid factors interfere up to concentration of 100 IU/ml.

Other substances may interfere.

NOTES

1. High CRP concentration samples may give negative results (prozone effect). Re-test the sample again using a drop of 20 µl.
2. The strength of agglutination is not indicative of the CRP concentration in the samples tested.
3. Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

BIBLIOGRAPHY

Lars-Olof Hanson et al. Current Opinion in Infectious diseases 1997; 10:196 - 201.

M.M. Pepys. The lancet 1981; March 21:653 - 656.

Chetana Vaishnavi, Immunology and Infectious Diseases 196;6 : 139 - 144.

Yoshitsugy Hokama et al. Journal of Clinical Laboratory Status 1987 ; 1 : 15 - 27.

Yamamoto S et al. Veterinary Immunology and Immunopathology 1993; 36:257 - 264.

Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACC Press. 1995.

SYMBOLS

Red instruction for use

Biological risk

CE mark (requirement of 98/78 CE regulation)

Storing temperature limits

In vitro medical device

Producer