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Fabbricante

---

**Rapid slide test for detection of anti-n-DNA**

---

**PRINCIPLE**

Lupus Erythematosus (LE) is an autoimmune disease in which antibodies against deoxyribonucleic acid (DNA), as well as other nuclear components, are produced.

In 75-80% of those patients suffering from LE, it can be detected what has been called the LE cell, the formation of which is believed to be due to those antibodies.

When the latex reagents, which has been coated with native deoxyribonucleic acid (n-DNA), is brought into contact with the serum sample, a clear agglutination will be observed if anti-n-DNA antibodies are present (positive result).

**REAGENTS**

**R1** Latex reagent: suspension of polystyrene latex particles coated with deoxyribonucleoproteins in a stabilization buffer.

**R2** Positive control: pool of human sera known to have a positive reaction with anti n-DNA latex reagent.

**R3** Negative control: pool of human sera known to have a negative reaction.

- Slide
- Stirrer

**SAMPLES**

Fresh serum or stored at 2-8°C for no longer than 48 h. It is necessary to freeze the sample when the assay is to be carried out after that period of time. Discard contaminated or hemolyzed sera.

**PRECAUTION IN USE**

 *Components from human origin have been tested and found to be negative for the presence of HbsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.*

**Caution:** the reagents contain Sodium Azide (0.095%) as preservative. Avoid swallowing and contacting with skin, eyes and mucous membranes.

**STORAGE AND STABILITY**

When stored at 2-8°C the components of this kit will remain stable until the expiration date stated on the label. Do not freeze. Discard those controls in which, in spite of containing sodium azide, a microbial growth is apparent.

**Reagents deterioration:** presence of particles and turbidity.

**PROCEDURE**

1. Bring reagents and samples to room temperature prior to assay.
2. Place one drop (40 µl) on the sample onto a slide circle.
3. Shake well the latex reagent and add one drop over the serum sample. Mix with the stirrer.
4. Tilt the slide and observe the presence or absence of agglutination after 3 minutes.

**READING AND INTERPRETATION**

**Positive:** agglutination after 3 minutes

**Negative:** no agglutination after 3 minutes

**SEMIQUANTITATIVE PROCEDURE**

Using physiological saline, dilute the specimens 1:2, 1:4, 1:8, 1:16... and proceed as stated in procedure.

**Titer:** last dilution that produces latex agglutination.

**QUALITY CONTROL**

Positive and negative control sera should be included in each test series.

**PERFORMANCE CHARACTERISTICS**

The reagent is designed to agglutinate in the presence of levels of anti-nuclear antibodies (ANA) commonly found in systemic lupus erythematosus.

Test gave positive reactions with 92% of the samples having titers of ANA  $\geq 1/20$  by immunofluorescence method.

No false positive reactions were observed with latex reagent in testing performance of samples from healthy individuals.

**INTERFERENCES**

Highly lipemic or haemolysed samples may interfere the agglutination reaction.

**LIMITATIONS OF THE PROCEDURE**

Those patients suffering from a variety of connective tissue diseases, may agglutinate the reagent.

Latex anti-n-DNA reagent does not give a diagnosis of the disease. Results have to be elevated within a total clinical study.

**BIBLIOGRAPHY**

Christian, C.L., Menchez-Bryan, R. and Larson, D.L., (1958). Proc. Soc. Exptl. Biol. Med. 98, 820-823.

Friou, G.J., Finch, S.C., Detree, K.D. (1958). J. Immunol., 80, 324-329.

Chaomao, J.C., (1976). Ara. J. Med. Tech., 42, 154-157.

---

**SYMBOLS**

Red instruction for use

Biological risk

CE mark (requirement of 98/78 CE regulation)

Storing temperature limits

In vitro medical device

Producer