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Fabbricante

Agglutination latex test without predilution for determination of heterophile antibodies

PRINCIPLE

The IM latex is a slide agglutination test for the qualitative and semiquantitative detection of heterophile antibodies (HE) specific for infectious mononucleosis (IM).

Latex particles coated with antigenic extract of beef erythrocytes membranes are agglutinated when mixed with sample containing Im heterophile antibodies.

CLINICAL SIGNIFICANCE

Infectious mononucleosis is a viral disease caused by Epstein-Barr virus that affects the reticuloendothelial system and has a broad spectrum of clinical presentations, ranging from asymptomatic to severe. The patients usually develop transient IgM heterophile antibodies, have an abnormal white cell picture, and abnormal liver function.

Disease diagnosis is obtained through the detection of HE antibodies or Paul-Burnell antibodies, or antibodies anti-viral structural antigens. The former generally decrease along the disease course, while the later remain along the patient life.

REAGENTS

R1 Latex particles coated with antigenic extract of beef erythrocytes membranes, phosphate buffer, pH 7.2. Preservative.

R2 Positive control red cap: human serum with an anti-IM antibodies titer $\geq 1/4$. Preservative.

R3 Negative control blue cap: animal serum. Preservative.
- Slide with 6 divisions
- Stirrer

SAMPLES

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C.
Samples with presence of fibrin should be centrifuged.
Do not use highly hemolyzed or lipemic samples.

PRECAUTION IN USE

 *Components from human origin have been tested and found to be negative for the presence of HbsAg, HCV and antibody to HIV (1/2). However handle cautiously as potentially infectious.*

Caution: the reagents contain Sodium Azide (0.095%) as preservative. Avoid swallowing and contacting with skin, eyes and mucous membranes.

CALIBRATION

The reagent sensitivity has been standardized an Internal Control by comparing methods with the Davidsohn method.

STORAGE AND STABILITY

All the kit components are ready to use, and will remain stable until the expiration date printed on the label, when stored tightly closet at 2-8°C and contaminations are prevented during their use. Do not freeze: frozen reagents could change the functionality of the test

Reagents deterioration: presence of particles and turbidity.

ADDITIONAL EQUIPMENT

Mechanical rotator with adjustable speed at 80/100 r.p.m.
Vortex mixer and pipettes (50ul)

QUALITATIVE PROCEDURE

1. Allow the reagents and samples to reach room temperature. The sensitivity of the test may be reduced at low temperatures.
2. Place 50µl of the sample and one drop of each positive and negative controls into separate circles on the slide test.
3. Mix the IM latex reagent vigorously or on a vortex mixer before using and add one drop (50 µl) next to the sample to be tested.
4. Mix the drop with a stirrer, spreading them over the entire surface of the circle. Use different stirrer for each sample.
5. Place the slide on a mechanical rotator at 80-100 r.p.m. for 2 minutes. False positive results could appear if the test is read later than two minutes.

SEMIQUANTITATIVE PROCEDURE

1. Make serial two fold dilutions of the sample in 9 g/l saline solution.
2. Proceed for each dilution as in the quantitative method.

READING AND INTERPRETATION

Examine macroscopically the presence or absence of visible agglutination immediately after removing the slide from the rotator.

The presence of agglutination indicates a titer $\geq 1/28$ of the specific anti-IM antibodies by the Davidsohn method.

The titer in the semi-quantitative method, is defined as the highest dilution showing a positive result.

QUALITY CONTROL

Positive and negative controls are recommended to monitor the performance of the procedure, as well as a comparative pattern for a better result interpretation. All results different from the negative control results will be considered as a positive.

PERFORMANCE CHARACTERISTICS

1. Analytical sensitivity: titer equal to 1/28 by the Davidsohn method under the described assay conditions.
2. Prozone effect: no prozone effect was detected up to 1/256 titer.
3. Diagnostic sensitivity: 100%
4. Diagnostic specificity: 100%

INTERFERENCES

Hemoglobin does not interfere up to concentration of 10 g/l.

Bilirubin does not interfere up to concentration of 20 mg/dl.

Lipemia does not interfere up to concentration of 10 g/l.

Rheumatoid factors do not interfere up to concentration of 300 IU/ml.

Other substances may interfere.

LIMITATIONS OF THE PROCEDURE

- False positive results may be obtained in some geographical areas where the "horse serum" is used as a prophylactic measure (vaccination).

- Patients suffering from leukaemia, Burkitt's lymphoma, pancreatic carcinoma, viral hepatitis, CMV infections and others, can result false positive reactions.

- False negative result have been reported in case of IM which persistently remain seronegative for IM heterophile antibodies or as a consequence of a delay IM heterophile antibodies response. In this case, repeat testing samples obtained at intervals of several days.

- Clinical diagnosis should not be made on findings of a single test result, but should integrate both clinical and laboratory data.

BIBLIOGRAPHY

Summaya CV et al. Manual of Clinical laboratory Immunology 4th ed p.568 Wasington Dc ASM, 1992.

Merlin JR et al. Human Pathol. 1986 ; 17 :2.

Paul JR et al. AM J Med Sci 1932 ; 183 :90.

Andiman WA. Manual of Clinical Laboratory Immunology, 3rd ed. P. 509 Wasigton, DC AMS 1986.

Henle W et al. Human Path 1974 ; 5 :551.

Barbara A Levey et al. Journal of Clinical Microbiology 1980;; 11:256-262.

Young DS. Effects of drugs on Clinical Laboratory test, 4th ed AACC press. 1995.

SYMBOLS

Red instruction for use

Biological risk

CE mark (requirement of 98/78 CE regulation)

Storing temperature limits

In vitro medical device

Producer