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Fabbricante

Rapid test on card for HCG detection in urine

PRINCIPLE

Human Chorionic gonadotropin (HCG) is a glycoprotein hormone secreted by the developing placenta shortly after fertilization. In normal pregnancy HCG may become detectable in serum since 7 days after conception and double every 1, 2-3 days. At the first missed menses the HCG concentration is about 100 mIU/ml, but the maximum value of 100-200.000 mIU/ml is detectable at the end of third month. After seven days until 100-200.000 mIU/ml at the end of the third month. HCG is a rapid test for qualitative detection in serum and urine. The method, using a combination of monoclonal antibodies conjugated with chromogen and polyclonal antibodies on membrane, allows to detect HCG levels greater than 25 mIU/ml within 5 minutes.

REAGENTS

Cards
Disposable plastic dropper

STORAGE AND STABILITY

- Cards are stable up to expiry date indicated on the label, stored at room temperature (15-25°C) or refrigerated (4-8°C).
- Do not freeze the test kit.

PREPARATION OF REAGENTS

Reagents are ready to use.

SAMPLE

Urine.

Urine:

- A first morning urine specimen is preferred since it contains the highest concentration of HCG.
- Collect the urine specimen in a clean container without preservatives.
- If testing is not immediate, the specimen should be stored at 4-8°C (or kept below 25°C) for up to 24 hours. Bring it to room temperature prior testing.
- The specimen can be stored at -20°C up to 72 hours. It must be thawed and brought to room temperature prior to testing. Avoid repeated freezing and thawing.

PROCEDURE


- Bring cards and specimens at room temperature before testing.
- Using the dropper provided, dispense 3-4 drops (~200 µl) of serum or urine in the sample well (⇒).

Read result when the control line in the (C) well is clearly marked. It's advisable read at most within 1 hour for urine and 10 minutes for serum.


IMPORTANT: While positive results can be detected within 2-3 minutes, for samples with low HCG values it is necessary wait 5-10 minutes to complete colorimetric reaction. Therefore, if negative result can be detected only after 10 minutes, positive test can be noticed shortly after control line in (T) well is clearly visible.

READING AND INTERPRETRATION

POSITIVE: in presence of a coloured band in T well confirmed by the control band in C zone.

NEGATIVE: only in presence of the control band in C zone and no visible band in T well.

ERROR: no distinct colour bands visible or one band in T well. It's advisable repeat the test for probable procedure imprecision.

NOTE

HCG level of m25 mIU/ml can be reached shortly after 7-10 days since conception. The maximum HCG level is about after 8 weeks since last menses and decrease during pregnancy. An high HCG level (> 200.000 mIU/ml) can give a wrong negative result. In case of suspicious pregnancy it's advisable repeat the test diluting the sample 1:1.

PERFORMANCES

Correlation-accuracy

A study was performed using a total of 150 positive and negative urine specimens, comparing with Hybritech Tandem II.

150 urine samples	Tests	
	SGM	Hybritech Tandem II
+	90	90
-	60	60

The results show a good correlation between SGM and Hybritech Tandem II. The accuracy study was also performed by determining the qualitative recovery of known amounts of human chorionic gonadotropin added to the negative urine pool.

RESULTS		
Final HCG concentration in the pool (en mIU/ml)	Expected	Observed
0	Negative	Negative
25	Positive	Positive
50	Positive	Positive
100	Positive	Positive

Specificity

Cross reactivity

The following concentration of homologous hormones are found to have no interferences with HCG test:

hTsh	1000 uIU/ml	WHO 68/38
hLH	500 mIU/ml	WHO 2nd IS 80/552
hFSH	100 mIU/ml	WHO 1st IS 83/575

Menopausal urine

A study was performed using specimens from 20 post-menopausal women because these urines frequently can interfere with pregnancy test due to cross reactivity with other gonadotropin hormones. All 20 urine specimens were negative when tested with HCG test.

Interferences

Potentially interfering substances were added to urine which have HCG levels of 0 and 50 mIU/ml. In each case, no interference with the expected results of test were observed:

Acetaminophen	20 mg/dl
Acetylsalicylic acid	20 mg/dl
Ascorbic acid	20 mg/dl
Atropine	20 mg/dl
Caffeine	20 mg/dl
Centesic acid	20 mg/dl
Glucose	2 g/dl
Haemoglobin	1 mg/dl
Ampicillin	20 mg/dl
Tetracycline	20 mg/dl

“Intra-Assay” precision (within-Run)

Within run precision was determined by using 10 replicates of three specimens containing 0, 25, 100 mIU/ml of HCG. The negative and positive values were correctly identified 100% of the time.

“inter-Assay” precision (between-Run)

Between run precision was determined by using the same three specimens of 0, 25, 100 mIU/ml of HCG in 10 independent assays and with three different lots of reaction device over a 6 months period. The negative and positive values were correct 100% of the time.

Sensitivity

- Test sensitivity is 25 mIU/ml according to WHO 3rd International Standard. Specimens containing less than 25 mIU/ml may also yield the positive results.
- Specimens containing high levels of HCG (1.000.000 mIU/ml) when tested gave positive results.
- Urine from healthy and non-pregnant women contain HCG levels not detectable by this test.
- This test is capable to detect pregnancy by the first day of the missed menstrual period, although the HCG levels are variable in the early pregnant women.

LIMITATIONS OF THE TEST

- In addition to pregnancy, HCG has been found in patients with trophoblastic disease. Since the HCG of trophoblastic neoplasms is similar to that found in pregnancy, these conditions, which include choriocarcinoma and hydatichiform mole, should be ruled out before a diagnosis of pregnancy is reached.
- A normal pregnancy can be distinguished from an ectopic pregnancy based on HCG levels alone. Also, spontaneous miscarriage may cause confusion in interpreting test results.
- A very early pregnancy containing an extremely low concentration of HCG can give a negative result. In this case, repeat the test after at least 48 hours.
- HCG levels remain detectable for several weeks after normal delivery, delivery by caesarean section, spontaneous or therapeutic abortion.
- Some serum specimens with high rheumatoid factor concentration (RF) may yield non specific positive results during testing. Such cases should be discriminated before testing.
- The physician should evaluate data obtained by the use of this kit in light of other clinical information.

WASTE MANAGEMENT

Please refer to the local legal requirements.

Note

The number indicated on the kit components identify the production lot. The number indicated on the external label identify the kit package. Any complaint must be furnished with one of these numbers.

BIBLIOGRAPHY

Braunstein, G.D., Vaitukaitis, J.L., Carbone, P.P., and Ross, G.T., Ann. Inter. Med., 78, 39-45 (1973).
 Morgan, F.J., Canfield, R.E., Vaitukaitis, J.L., and Ross, G.T., Endocrinology, 94, 1601-1606 (1974).
 Kohler, G. and Milstein, C., Nature, 256, 495-497 (1975).
 Braunstein, G.D., Rasor, J., Adler, D., Danzer, H., and Wade, M.E. am. J. Obstet. Gynecol., 126, 678-681 (1976).
 Thompson, R.J., Jackson, A. P., and Langlois, N., clin. Chem., 32, 476-481 (1986).
 Engvall, E., Methods in Enzymology, 70, 419-439 (1980).
 Rasor, J.L., and Braunstein, G.D., Obstet. Gyneco., 50, 553-558 (1977).
 Lenton, E.A., Neal, L.M., and Sulaiman, R., Fertility and Sterility, 37, 773-778 (1982).

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer