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visual card

For in vitro medical device  
31028 - 10 test

## Rapid test on card for methamphetamine detection in urine

### PRINCIPLE

Methamphetamine is a potent sympathomimetic agent with therapeutic applications, too. Acute higher doses lead to enhanced stimulation of the central nervous system and induce euphoria, alertness, and sense of increased energy and power. More acute responses produce anxiety, paranoia, psychotic behaviour, and cardiac dysrhythmias. The pattern of psychotic at high doses may be indistinguishable from schizophrenia.

Methamphetamine is excreted as amphetamine and oxidized and deaminated derivatives. However 10-20% of methamphetamine is excreted unchanged. Thus the presence of the parent compound in the urine indicates methamphetamine use.

The easy card methamphetamine test is a "competitive" immunoassay for the quantitative detection of methamphetamine through the competition between the methamphetamine placed on the membrane and the free one, eventually in the tested sample, along the coloured anti-methamphetamine conjugated. Enough methamphetamine quantity, over cut-off level, will fill all antibody binding sites, will prevent attachment of the coloured anti-methamphetamine conjugated and the formation of coloured band on the membrane. Therefore, the appearance of one coloured band indicates a positive result. Instead, the formation of two coloured bands means a negative result.

This method represents a valid, simple and fast alternative to classic methods as gas chromatography (GC) or mass spectrometry (MS) which require complicated equipment and longer execution times. However, this assay provides only a preliminary analytical test result and, in case of positivity, it must be confirmed by other chemical researches.

### REAGENTS

Cards (MET)	10
Specimens pipettes	10

### STORAGE AND STABILITY

Cards are stable up to expiry date showed on the package, stored at room temperature (15-25°C) or refrigerated (4-8°C) .

### PREPARATION OF REAGENTS

Reagents are ready to use.

### SAMPLE

Fresh urine.

### QUALITY CONTROLS

Good laboratory practice recommends the use of control materials to ensure proper kit performances. When testing the positive and negative controls, use the same assay procedure as with a urine specimen.

### PROCEDURE


- Allow specimens and cards to come to room temperature prior to testing.
- Fill the urine dropper with specimens and, by holding it vertically, dispense 6 drops (200 µl) of urine in the sample well A of the card.
- Wait 3-5 minutes to read the results. Do not read after 8 minutes because the intensity of the coloured bands may change, or a new line may appear which was no present at the end of the incubation time (5 minutes). Therefore, it's most important to record the results of assay immediately after the five minutes incubation period.

### Prudences


- Urine specimens and all materials coming in contact with them should be handled and disposed as if it is infectious and capable of transmitting infection. Avoid contact with skin.
- Avoid cross-contamination of urine samples by using a new specimen cup and specimen pipette for each urine sample.
- Test the sample as soon as possible, preferably the same day of specimen collection.
- The urines may be refrigerated at 2-8°C for 2 days or frozen at -20°C for longer periods.
- Specimen that were refrigerated must reach room temperature (15-25°C) prior to testing.
- Specimen previously frozen must be thawed and mixed thoroughly prior to testing.

### INTERPRETATION OF RESULTS

**POSITIVE:** One colour band appears on the control window (C). The methamphetamine concentration is over cut-off levels.

**NEGATIVE:** Two colour bands appear, even if with weak intensity. The methamphetamine concentration is under cut-off levels.

**ERROR:** If there is no distinct colour band visible both in the test window and test control window the test is inconclusive. It is recommended to repeat the test.

### Note

- A very faint line on the test region indicates that the methamphetamine in the sample is near or below the cut-off level for the test. However, if any line is observed at the test region, the sample is negative.
- Confirm eventually positive results by opportune confirmation assay.

## ANALYTICAL PERFORMANCES

### Specificity

- The test is specific for methamphetamine and metylendioxiamphetamine.
- The test specificity was tested by adding various drugs, metabolites and other compounds that are likely to be present in urine. All compounds were prepared in drug-free normal urine.

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Compound	Concentration (ng/ml)
D-amphetamine	50,000
Chloroquine	50,000
(+/-) Ephedrine	50,000
(+)-Methamphetamine	1,000
(-)-Methamphetamine	25,000
(+/-)-3,4-methyendedioxymethamphetamine	2,000
Procaine	10,000
β-Phenylethylamine	50,000
Ranitidine	50,000

- The following compounds were found not to cross-react when tested at concentrations up to 100 µg/ml:

Acetaminophen	Dexbrompheniramine
Acetone	Dextromethorphan
Albumin	Diazepam
L-Amphetamine	4-Dimethylaminoantipyrin
Ampicillin	Dopamine
Aspartame	Doxylamine
Aspirin	Ecgonine
Atropine	Ecgonine methyl ester
Benzocaine	(+/-)-Ephedrine
Benzoyllecgonine	(+)-Epinephryne
Bilirubin	Erythromycin
(+)-Brompheniramine	Ethanol
Caffeine	Furosemide
(+)-Chlorpheniramine	Glucose
(+/-)-Chlorpheniramine	Guaiacol Glyceyl Ether
Chlorpromazine	Hemoglobin
Hydrocodone	Creatine
Hydromaphone	Imipramine
Pentermine	(+ -)-Isoproterenol
Phencyclidine	Lidocaine
Pheniramine	Meperidine
Phenobarbital	Methaqualone
Phenothiazine	
(+/-)-3,4-L-Phenylephrinemethylenedioxyamphetamine	
(1R,2S)-(-)-N-Methyl-Ephedrine	β-Phenylethylamine

Methylphenidate	Promethazine
Morphine	Quinidine
Naltrexone	Riboflavin
(+)-Naproxen	Secobarbital
(+/-)-Norephedrine	Sodium Chloride
Nortriptyline	Sulindac
Oxalic Acid	Thioridazine
Oxazepam	Trifluoperazine
Oxycodone	Tyramine
Acetylsalicylic Acid	Penicillin G
Vitamin C	

### Correlation-accuracy

The comparison between the results obtained by the kit methamphetamine easy-card and a commercial immunoenzymatic kit on positive samples and negative samples, gives a superposition of results at 100%.

### Sensitivity

The test sensitivity is  $\geq 1000$  ng/ml.

### LIMITATIONS OF PROCEDURE

- The assay is designed for use with human urine only.
- A positive result with any of the tests indicates only the presence of a drug/metabolite and does not indicate or measure intoxication.
- There is a possibility that other substances and/or factors not listed may interfere with the test and cause false results, such as technical or procedural errors. See "Specificity" for lists of substances that will produce positive results or that do not interfere with tests performance.

### WASTE MANAGEMENT

Please refer to the local legal requirements.

### Note

The number showed on single product refer to the production lot. The number showed on the exterior label refer to the packaging. In case of claims, please refer to one of these two number.

### BIBLIOGRAPHY

Baselt, R.C. Disposition of toxic Drugs and Chemical in Man, Davis, Biomedical Communication, 1982.  
Urine teting for Drugs of Abuse, National Institute on Drugs Abuse (NIDA), Research Monograph 73, 1986.  
Fed. Register, Department of Health and Human Services, mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69, 11970-11979, 1988.  
McBay, A.J. Clin. Chem. 33, 33B-40B, 1987.  
The Pharmacological Basis of Therapeutics, Gilman, A.G. Goodman, L.S. and Gilman A. eds. MacMillan Publishing, New York, NY, 1980.

### SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer