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BIBLIOGRAFIA

Clauss A.: Acta Haemat 17, 237-246 (1957)

NCCLS: Collection, Transport and Processing of Blood Specimens for Coagulation Testing and Performance of Coagulation Assays; 2nd edition Approved Guideline. NCCLS document H21-A3, Wayne PA. (1998)

NCCLS: Procedure for determining fibrinogen in plasma. Approved guideline. NCCLS Document H30-A2. Wayne, PA, (2001)

Musgrave K. A., Bick R.L.: Qualità assurance in the Haemostasis Laboratory. In Bick R.L. et al. Editors: Hematology: Clinical and Laboratory Practice Vol. 2, 1309 - 1315, Mosby St. Louis MO (1993)

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Fabbricante

Kit for fibrinogen determination – Clauss method

**Reagent R1 is classified as Xn = noxious
CONTAINS NATRIUM AZOTURE**

R22: noxious for ingestion

R52/53: noxious for aquatic organism, may be produce long-term negative effects for per aquatic environment

S36/37: use protective clothes and suitable gloves

S60: this material and its container must be removed as dangerous waste

S61: do not waste in environment. Refer to special instructions/informative sheet for security

PRINCIPLE

FIBRINOGEN kit measure clotting time, exclusive expression of fibrinogen quantity in the sample, according to the original method by Clauss. In presence of high concentrations of Thrombin, the time required for clot formation in diluted plasma is inversely proportional to the fibrinogen concentration.

Note: the dilution of considered plasma serves to minimized the interferences of any anticoagulant as heparin.

REAGENTS

- R1** Imidazole buffer pH 7.4
- R2** Bovin thrombin ~ 100 NIH units/ml
- R3** Control plasma normal (N)
- R4** Control plasma pathological (P)
- R5** Fibrinogen Calibrator

SAMPLE

- Fresh blood in 3.2 % (0.109M) trisodium citrate.
- Collect the sample and, avoiding haemolysis and contaminations, mix immediately blood with anticoagulant (9 parts blood and 1 of 3.2 % trisodium citrate). Sample that do not reach the volume of 90% must be discarded.
- Centrifuge at 2500 x g for 15 minutes and separate the supernatant plasma.
- Run the test immediately (within 2 h) if stored at 22-24 °C.

Note

 Calibrator (R5) and human Controls (R3 and R4) are obtained using only blood of donors tested by an FDA method and found non-reactive for HbsAg and negative for antibodies to HIV-1/2 and HCV. However, no known test method can offer complete assurance that products derived from human blood will not transmit infectious diseases, this product should be handled as potentially infectious biological material.

- Every time infectious agents, chemical reagents, human or animal reagents, blood or other biological fluids are handled it's advisable follow common advices and all hygienic precautions as disposable gloves.

PREPARATION OF REAGENTS

- Reagent **R1** is liquid and ready to use.
- Reagent **R2:** reconstitute bovin thrombin vial with 2 ml of distilled water. Agitate gently until solution is complete.
With a photomeccanic device is suggested to reconstitute Bovin thrombin with 2 ml of Caolin instead of distilled water.
- Reagent **R3-R4-R5:** reconstitute with 1 ml of distilled water. Agitate gently until solution is complete and let it stand, before using, 15 minutes at room temperature.

Note

- It is essential that the reagent is at 15-25 °C when it is used.

STORAGE AND STABILITY

- Store the kit at 2-8 °C.
- Reagents, if used and stored according to good laboratory practice, are stable up to the expiry date showed on the label.
- Reconstituted **R2** is stable 7 days at 2-8 °C and 30 days if frozen within 4 hours from the reconstitution; in this case, when it is used, thaw rapidly at 37 °C. Do not refreeze.
- Reconstituted **R3, R4 and R5** are stable 8 hours stored at 2-8 °C.

PRECAUTION IN USE

Look at the caution advices indicated nearby the symbol.

Caution: the reagents contain Sodium Azide (0.095%) as preservative. Avoid swallowing and contacting with skin, eyes and mucous membranes.

WASTE MANAGEMENT

Please refer to the local legal requirements.

PROCEDURE

- Prepare, with imidazole buffer, 4 different scaled dilutions (1:5, 1:10, 1:15, 1:35) of calibrator for fibrinogen and test it as described below.
- Note:** the diluted plasma 1:10 represent 100% of the assigned value. The dilution factor indicates the relationship between the 1:10 dilution and the other dilutions.
- Dilute the supernatant (tested sample) and the possible quality controls (ricostituted with imidazole buffer in 1:10).
- Distribute in plastic test tubes and run the tests of different dilutions of calibrator for fibrinogen (for reference curve), of controls and sample, both diluted, following this plan:

Diluted sample	200 µl
Incubate for 4-6 minutes at 37 °C then add	
Bovin thrombin	100 µl
At the Thromboplastin addition measure the time of clot formation. (seconds).	

PREPARATION OF CALIBRATION CURVE

- Build, on bilogarithmic scale, the calibration curve placing on the abscissa the concentrations (mg/dl) of scaled dilutions of calibrator for fibrinogen and on ordinate the corresponding times of clot formation.

Example: calibrator for fibrinogen = 304 mg/dl

DILUTION	DILUTION FACTOR	FIBRINOGEN CONCENTR. (mg/dl)	COAGULATION MEAN TIME (sec)
1:5	10/5 = 2	304 x 2 = 608	6.6
1:10	10/10 = 1	304 x 1 = 304	12.3
1:15	10/15 = 0.67	304 x 0.67 = 204	18.6
1:35	10/35 = 0.29	304 x 0.29 = 88	24.2

- Obtain the fibrinogen concentrations of sample and control placing on the calibration curve their times of clot formation.

REFERENCE VALUES

Generally the normal reference range is 150- 350 mg/dl. However, It's advisable that every laboratory determine its own reference values.

ANALYTICAL PERFORMANCES

Accuracy

A low, a normal, and a high fibrinogen plasma were tested in multiple laboratories using FIBRINOGEN and the results were compared to results obtained using other manufacturer's reagents.

The results are reassumed as follows:

N	SAMPLE	FIBRINOGEN	N	OTHER REAGENTS
10	Low	144 mg/dl	196	163 mg/dl
10	Normal	294 mg/dl	196	297 mg/dl
10	High	488 mg/dl	390	474 mg/dl

Precision

A low, a normal and a high fibrinogen plasma were tested on multiple days using FIBRINOGEN on a photo-optical instrument.

Ten standard curves were determined on each test day.

The percent CV of each plasma were:

5.9% for low plasma

3.4% for normal plasma

2.9% for high plasma

QUALITY CONTROLS

- It's advisable, with sample or at least with every tests group, for lot changing or for important instrument adjustment run a normal and pathological control.

- The normal and pathological control included in the kit correspond to a lyophilised normal plasma and a correct plasma to simulate a plasma with moderate deficit.

LIMITATIONS

- Blood must be immediately added to trisodium citrate anticoagulant and gently mixed (DO NOT USE EDTA AND HEPARIN).

- Hemolysis can cause clotting factor activation and end point detection interference. Icteric and lipemic specimens may also be inappropriate for end point detection methods.

- The sample should only contact nonwetable surfaces.
- The ratio of blood to anticoagulant is usually 9:1 and results in a citrate concentration of 10.9 to 12.9 mmol/l
- Freezing and thawing plasma that contains residual cells will generate damaged cells membranes that can affect results.
- Acute inflammatory reactions can elevate circulating Factor I (fibrinogen).
- High Fibrinogen degradation products (FDP) may prolong clotting times especially when the fibrinogen level is below 150 mg/dl.
- In patients with qualitative fibrinogen abnormalities, the thrombin clotting time assay may indicate decreased fibrinogen. The quantitative fibrinogen results may be normal on these same samples if tested by other methods.
- Heparin does not interfere at therapeutic levels. However, very high heparin levels may cause low fibrinogen results. Batroxobin enzyme can be substituted for thrombin in this assay if heparin interference is suspected.
- High paraprotein levels, thrombin antibodies, and drugs that activate the fibrinolytic system can interfere with fibrinogen assays.

BIBLIOGRAPHY

Clauss A.: Acta Haemat 17, 237-246 (1957)

NCCLS: Collection, Transport and Processing of Blood Specimens for Coagulation Testing and Performance of Coagulation Assays; 2nd edition Approved Guideline. NCCLS document H21-A3, Wayne PA. (1998)

NCCLS : Procedure for determining fibrinogen in plasma . Approved guideline. NCCLS Document H30-A2. Wayne , PA, (2001)

Musgrave K. A., Bick R.L.: Qualità assurance in the Haemostasis Laboratori. In bick R.L. et al. Editors: Hematology: Clinical and Laboratori Practice Vol. 2, 1309 - 1315, Mosby St. Louis MO (1993)

SYMBOLS

Read instruction for use

Biological risk

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer