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
### Bibliografia

National Committee for Clinics Laboratory Standards. Collection, Transport and processing of blood specimens for coagulation testing and performances of coagulation assay. 2<sup>nd</sup> edition, Approved guideline NCCLS Document H21-A2 villanova PA (1991)

### SIMBOLOGIA

## MULTICOMPONENT CONTROL PLASMA FOR COAGULATION ASSAYS

### Principle

CONTROL PLASMA N is a lyophilised human control plasma and is a valuable tool for assuring quality in coagulation testing. It is better a normal and pathological control for every dosage group, for 40 samples of the same test and in case of operator change.

### Reagents

- Human plasma collected with sodium citrate as anticoagulant.  
- The single component concentration is reported on the table joined every package.

### Note

 Human control is obtained using only blood of donors tested by an FDA method and found non-reactive for HbsAg and negative for antibodies to HIV-1/2 and HCV. However, no known test method can offer complete assurance that products derived from human blood will not transmit infectious diseases, this product should be handled as potentially infectious biological material.

- Every time infectious agents, chemical reagents, human or animal reagents, blood or other biological fluids are handled it's advisable follow common advices and all hygienic precautions as disposable gloves.

### Storage and stability

- Lyophilized Control is stable up to expiry date showed on the label and stored at 2-8°C.  
- Reconstituted Control is stable 8 hours stored at 2-8°C. Before using shake slowly.

### Instructions for use

- Open the vial and add 1 ml of distilled water.  
- Close the vial, mix slowly to dissolve all the content and let it stand for 15 minutes at least at 15-25°C.  
- Treat the control in the analytical procedure as patient sample.  
- Before using, the control must be at 15-25°C.

### Note

- Control sometimes can be cloudy but this is a property of lyophilised matter. It's advisable discard the product only in case of evident contamination.

### Auxiliary materials not provided

- Precision pipettes, distilled water.

### Performances

- The Control is formulated to have PT, APTT and Fibrinogen values in the normal range.

### Precaution in use

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration of components is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances. However the reagent should be handled with caution, according to good laboratory practice.

### Waste management

Please refer to the local legal requirements.

### Bibliography

National Committee for Clinics Laboratory Standards. Collection, Transport and processing of blood specimens for coagulation testing and performances of coagulation assay. 2<sup>nd</sup> edition, Approved guideline NCCLS Document H21-A2 villanova PA (1991).

### SYMBOLS