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Fabbricante

Kit for prothrombin time determination

Quick method

PRINCIPLE

Prothrombin or Quick time is the time, in seconds, necessary for a plasma to clot at 37°C, once an external coagulation factor (thromboplastin) and Ca^{2+} has been added. Abnormally long times with regard to those obtained for healthy people indicate a deficiency in some of the components of the coagulation system.

The kit is a calcium thromboplastin, freeze-dried, that allows the determination of variations in those coagulation factors (especially factors VII and X) that have been depressed during the treatment with coumarinic derivatives, in vitamin K deficiency, hepatopathy or due to a hereditary abnormality in the coagulation system.

The reagent has been obtained from a rabbit brain extract, and it is standardized from lot to lot according to the technique of Hills-Ingram.

REAGENTS

R1 Calcium thromboplastin from rabbit brain extract

SAMPLE

- Fresh blood in 3.2 % (0.105M) trisodium citrate.
- Collect the sample and, avoiding haemolysis and contaminations, mix immediately blood with anticoagulant (9 parts blood and 1 of 3.2 % trisodium citrate).
- Centrifugate at 1500 x g for 15 minutes and separate the supernatant plasma.
- If testing is not completed immediately, the sample can be stored at +2-8°C for 24 hours, at -20°C for 15 days or at -70°C for 6 months.

Note

It's advisable follow the common instructions and take all necessary hygienic precautions (as using disposable gloves) in case of handling of infectious agents, chemical reagents, human or animal reagents, blood or other biological liquids.

PREPARATION OF REAGENTS

- Reconstitute the reagent vial content with 4 ml of distilled water.
- Mix slowly (do not invert or mix vigorously) and let it stand for 15 minutes at 15-25°C. Do not freeze.

Note

- It is essential that the reagent is at 15-25°C when it is used.

STORAGE AND STABILITY

- The kit is stable up to the expiry date showed on the label if stored at 2-8°C.
- The reconstituted reagent is stable for 7 days at 2-8°C and 8 hours at 37°C. Do not freeze.

AUXILIARY EQUIPMENT

Materials not included in the kit: precision pipette, stopwatch or timer, clotting instrument and controls.

PRECAUTION IN USE

The product is not classified as dangerous (DLg. N. 285 art. 28 l. n. 128/1998). The total concentration components is lower than the limits reported by 67/548 and 88/379 CE Regulations (and following modifications) about classification, packaging and labelling of dangerous substances. However the reagent should be handled with caution, according to good laboratory practice.

Caution: the reagents contain Sodium Azide (0.095%) as preservative. Avoid swallowing and contacting with skin, eyes and mucous membranes.

WASTE MANAGEMENT

Please refer to the local legal requirements.

PROCEDURE

- Dispense the supernatant (sample tested) into a plastic or siliconized glass test tube as follows:

| | |
|--|--------|
| Citrate plasma (supernatant) | 100 µl |
| Incubate for about 2 minutes at 37°C | |
| Reagent (prewarmed at 37°C) | 200 µl |
| At the reagent addition (Thromboplastin) measure the time of clot formation. Duplicate the determination (the difference between two values must not be above 5%; otherwise perform a third test). | |

INTERPRETATION OF RESULTS

The prothrombin time can be expressed in three ways:

- Percentage (Quick time)
- Ratio
- I.N.R. (International Normalized Ratio)

Percentage

The values are obtained by calibration curve built using graduated dilution of a pool of normal plasma as in the table:

| Dilution Value % | None 100% | 1+1 50% | 1+2 33% | 1+3 25% | 1+7 12.5% |
|---------------------------|-----------|---------|---------|---------|-----------|
| Plasma | 0.5 ml | 0.5 ml | 0.5 ml | 0.5 ml | 0.5 ml |
| Physiological sol. | -- | 0.5 ml | 1 ml | 1.5 ml | 3.5 ml |

N.B. For every dilution determine the value twice or three times.

Ratio

The values are obtained by this formula: $\frac{\text{Prothrombin time (sample)}}{\text{Prothrombin time (normal control)}}$

I.N.R. = RATIO^{ISI}

An undesirable consequence of oral anticoagulant therapy may be a tendency to bleed unnecessarily. In order to maximize the desired therapeutic effects and minimize bleeding, the World Health Organization (WHO) has recommended a procedure to standardized testing and treatment. This procedure is based on the International Normalized Ratio (INR).

The 'I.N.R.' is calculated using the RATIO according to the following mathematical relationship:

$$I.N.R. = RATIO^{ISI}$$

For example, with an ISI of 1.00 and a RATIO of 3,2, the INR will be:

$$INR = (3,2)^1 = 3,2$$

The International Sensitivity Index (ISI) is a measure of a thromboplastin/instrument sensitivity to coagulation factors. ISI values are assigned by comparison to a primary reference material. High sensitivity reagents have low ISI values.

According to OMS recommendations, INR values above 5.5 place the patient at unnecessary risk for bleeding complications. It is generally advised that patients on stabilized oral anticoagulant therapy should be maintained at an INR of 2.0-3.5, depending on the clinical indication.

The lot specific ISI value for Thromboplastin can be found on the kit box label.

EXPECTED VALUES

In multi-center studies, when the product was evaluated on a normal population, the following results were obtained:

| N. SAMPLE | INSTRUMENT | MEAN PT (sec) | RANGE (± 2 DS) |
|-----------|-----------------------|---------------|----------------|
| 40 | MILA™ Electra 1000 C™ | 13.2 | 11.4 – 15.0 |
| 20 | MILA™ Electra 900 C™ | 13.7 | 12.4 – 15.0 |
| 61 | IL ACL™ 300 / 3000 | 10.5 | 8.9 – 12.1 |
| 20 | Amelung KC 10™ | 12.7 | 9.3 – 14.2 |
| 38 | TromboScreen 400C | 13.5 | 12.2 – 14.8 |
| 60 | TromboScreen 200 | 13.5 | 12.0 – 15.1 |

- These values are intended as guideline only.
- Each laboratory should establish a normal reference range (NRR) using instrumentation, blood collection methods, and testing techniques used in that laboratory.
- The NRR should be reestablished or at least verified when changing lot numbers of the same reagent. A new NRR should be established with any change in reagents, instrumentation, blood collection techniques, or anticoagulant.
- The clotting time of abnormal plasmas will depend on the ISI of the reagent lot in use.

ANALYTICAL PERFORMANCES

Precision

The precision was verified by testing a normal and an abnormal plasma on several different instruments.

The %CV, using 20 samples, are reassumed in the following table:

| SAMPLE | MILA ELECTRA 1000 C | TROMBO SCREEN 400 C | TROMBO SCREEN 200 | AMELUNG |
|----------|---------------------|---------------------|-------------------|---------|
| NORMAL | 1.1% | 1.9% | 1.9% | 2.9% |
| ABNORMAL | 2.8% | 2.5% | 2.3% | 1.1% |

Sensitivity

Factor sensitivity testing was performed by diluting a normal plasma with factor deficient plasma such that the final factor concentration ranged from 10-100% and using MLA 1000 C instrument.

| % FACTOR | PROTHROMBIN TIME | | | |
|----------|------------------|----------|------------|----------|
| | FACTOR II | FACTOR V | FACTOR VII | FACTOR X |
| 100 | 11.6 | 11.6 | 11.8 | 11.7 |
| 50 | 11.6 | 13.2 | 12.6 | 12.8 |
| 40 | 11.7 | 13.9 | 12.8 | 13.3 |
| 30 | 12.3 | 14.9 | 13.5 | 14.1 |
| 20 | 12.8 | 15.9 | 13.9 | 14.8 |
| 10 | 14.1 | 18.3 | 15.2 | 17.0 |

Correlation

The comparison between this product (y) and two other marketed thromboplastin reagents (x) performed on the Stago STA instrument, gives the following results:

With the first product

N = 49
Correlation PT
 $r = 0.98$
 $y = 1.16x + 1.3$

Correlation INR
 $r = 0.98$
 $y = 0.89x + 0.05$

With the second product

N = 49
Correlation PT
 $r = 0.95$
 $y = 1.01x + 2.20$

Correlation INR
 $r = 0.95$
 $y = 0.82x + 0.10$

Interferences

- Sodium oxalate, EDTA and heparin as anticoagulants.
- Substances as oral contraceptives, corticosteroids, EDTA, asparaginase, clofibrate, erythromycin, ethanol, tetracycline, and anticoagulants as heparine and warfarin may prolong the PT.
- The PT may be shorted by substances including antihistamines, butabarbital, caffeine, oral contraceptives, phenobarbital and vitamin K.

QUALITY CONTROLS

- It's advisable, with sample or at least with every tests group, for lot changing or for important instrument adjustment run a normal and pathological control:

REF 61004 Control Plasma N

REF 61005 Control Plasma P

- The first product correspond to a lyophilised normal plasma, the second correspond to a correct plasma to simulate a plasma with moderate deficit. In laboratories with many work of PT and/or APTT, it's advisable a normal and abnormal control every 40 samples at least.
- Each laboratory should establish a control group range to represent the allowable variation in day to day performance for each control.

LIMITATIONS

The biochemistry of coagulation involves a series of reactions that are influenced by many pre-test conditions. To obtain reproducible results these variables must be controlled as for example:

- The pH increase if plasma is open to air. Store plasma tightly closed.
- The product was designed to work at $37^{\circ}\text{C} \pm 0.5^{\circ}\text{C}$. Ensure that all heating elements are constantly at this temperature.
- All labware must be clean and free of trace amounts of detergents.

NOTE

- To maintain the suspension in optimal conditions, during the use, prepare a mechanism as a magnetic agitator.
- Lack of vacuum in vials, erratic results, quality control values outside established rangel, or product color variations could indicate deterioration.
- Turbid, icteric, lipemic or hemolyzed specimens may generate erroneous results.
- Use only plastic or siliconized glass containers.
- Freezing and thawing plasma that contains residual cells will generate damaged cells membranes that can affect results.
- Acute inflammatory reactions can shorten PT results because of elevated fibrinogen.
- Plasma samples with hematocrits outside the range of 20-55% may be improperly anticoagulated and should be adjusted appropriately.
- The product is suitable for use with manual, mechanical, photoptical, nephelometric etc; to have a correct use with the instrument it's necessary follow the instructions from instrument producer.

BIBLIOGRAPHY

- Erricchetti A. M., golden A., Ansell J.: mangement of Oral Anticoagulant Therapy. Experience with an Anticoagulation Clinic. Arch Inter Med 144, 1966-1968 (1984).
- Hirsh J., Dalen J. E., Deykin D., Poller L.: Oral Anticoagulants: Mechanism of Action, Clinical Effectiveness and Optimal Therapeutic Ramge. Chest 102 (Suppl) 312s-315s (1992).
- NCCLS: Collection, Transport and Processing of Blood Spocimens for Coagulation Testing and General Performance of Coagulation Assays; Approved Guideline. NCCLS document H21-A3, NCCLS Wayne PA. (1998).
- Palmer R.N., Grainick H.R.: Inhibition of the Cold Activation of Factor VII and the Protrombin Time. Am. J. Clin. Path 81, 618-622 (1984).
- Young D.S., Thomas D.W., Friedman R.B. et al: Effect of drugs on Clinical Laboratory Tests Clin Chem 18, 1041 (1972).
- NCCLSA: One-Stage Protrombin Time (PT) Test and Activated Partial Thromboplastin Time (APTT) Test; Approved Guideline. NCCLS document H47-A, NCCLS Wayne, PA. (1996).
- Dalen J.E., Hirsh J. : American College of Chest Physicians and the National Hearth, Lung, and Blood Institute National Conference on Antithrombotic Therapy. Arch. Inter. Med. 146, 462-472 (1986).
- Palaereti G., Cocchieri S., Poggi M. et al: Oral Anticoagulant Therapy Control; Evidence that the INR Expression Improves the interlaboratory Comparabilità of results. The Bologna Oral Anticoagulant Control Exercise. Tromb Haemostasis 58: 905-910 (1987).

SYMBOLS

Read instruction for use

CE mark (requirement of 98/79 regulation)

Storing temperature limits

In vitro medical device

Producer